



Chiral Shift Reagent for Determination of Enantiomeric Excess

Chirabite-AR 100ma

Advantages

- Available for both high- and low-field NMR spectrometers because C2184 contains no paramagnetic metals, which cause signal broadening.
- The enantiomeric purities of various compounds can be determined.
- NMR spectra showing chemical-shift nonequivalences can be obtained just by adding C2184 to the NMR tube containing a target sample in CDCl₃.

Applications

Chart 1.







Chart 2.





Chirabite-AR (7 mg) in CDCl₃ at 22°C

References T. Ema, D. Tanida, T. Sakai, Org. Lett. 2006, 8, 3773. https://doi.org/10.1021/ol0613665 T. Ema, D. Tanida, T. Sakai, J. Am. Chem. Soc. 2007, 129, 10591. https://doi.org/10.1021/ja073476s

For further information please refer to our website at www.TClchemicals.com. >>> TCI NMR shift





Ordering and **Customer Service**

TCI AMERICA

:800-423-8616 / 503-283-1681 Tel :888-520-1075 / 503-283-1987 E-mail: Sales-US@TCIchemicals.com

TCI EUROPE N.V.

: +32 (0)3 735 07 00 : +32 (0)3 735 07 01 E-mail: Sales-EU@TCIchemicals.com

TCI Deutschland GmbH

Tel : +49 (0)6196 64053-00 : +49 (0)6196 64053-01 E-mail: Sales-DE@TCIchemicals.com

Tokyo Chemical Industry UK Ltd. Tokyo Chemical Industry (India) Pvt. Ltd.

: +44 (0)1865 78 45 60 E-mail: Sales-UK@TCIchemicals.com

梯希爱(上海)化成工业发展有限公司

Tel :800-988-0390 / 021-67121386 : 021-6712-1385 E-mail: Sales-CN@TCIchemicals.com

: 1800 425 7889 / 044-2262 0909 E-mail: Sales-IN@TCIchemicals.com

TOKYO CHEMICAL INDUSTRY CO., LTD.

Tel: +81 (0)3-5640-8878

E-mail: globalbusiness@TCIchemicals.com

[•] Chemicals itemized in this brochure are for research and testing use only. Please avoid use other than by chemically knowledgeable professionals. • Information such as listed products and its specifications and so on are subject to change without prior notice. • The contents may not be reproduced or duplicated in whole or in part without permission of Tokyo Chemical Industry Co., Ltd.