



# MM3D-Fluorogenic Probe for Real-Time Measurement of ENGase Activity



MM3D is hydrolyzed by *endo*- $\beta$ -N-acetylglucosaminidase (ENGase) to generate fluorescence. MM3D is a useful tool for real-time measurement of the enzymatic activity of ENGase.

## Principle



*Endo*- $\beta$ -N-acetylglucosaminidase (ENGase) first hydrolyzes the chitobiose moiety of N-glycan. MM3D, a N-glycan which features both a fluorescent group (N-methyl anthraniloyl, MANT) and a quencher group (dinitrophenyl, DNP) in the structure. Cleavage of MM3D by ENGase removes the self-quenching effect, leading to an increase in fluorescence intensity. MM3D has shown to be useful for high throughput analysis with plate readers, and in the search for new types of ENGase or inhibitor of ENGase. These inhibitors can be served as therapeutic agents for NGLY1-deficiency, a very rare genetic disease.

\*Note: MM3D is not hydrolyzed by Endo-H which belongs to GH Family 18, and therefore is not suitable for measurement of Endo-H activity.


## References

- 1) N. Ishii, C. Sunaga, K. Sano, C. Huang, K. Iino, Y. Matsuzaki, T. Suzuki, I. Matsuo, *ChemBioChem* **2018**, 19, 660.
- 2) H. Fujihira, Y. Masahara-Negishi, M. Tamura, C. Huang, Y. Harada, S. Wakana, D. Takakura, N. Kawasaki, N. Taniguchi, G. Kondoh, T. Yamashita, Y. Funakoshi, T. Suzuki, *PLoS Genet.* **2017**, 13, e1006696.
- 3) C. Huang, Y. Harada, A. Hosomi, Y. Masahara-Negishi, J. Seino, H. Fujihira, Y. Funakoshi, T. Suzuki, N. Dohmae, T. Suzuki, *Proc. Natl. Acad. Sci. U.S.A.* **2015**, 112, 1398.

This product has been commercialized under the instruction of Prof. Ichiro Matsuo.

# MM3D-Fluorogenic Probe for Real-Time Measurement of ENGase Activity

## Measurement of Endo-M [A1651] activity using MM3D<sup>1)</sup>



(Data are provided by Prof. Ichiro Matsuo.)

Comparing the fluorescence intensity of MM3D with that of the hydrolyzed product, the fluorescence of MANT is decreased by 98% due to the quenching effect of the DNP group (a). However, since the cleaved product of MM3D is highly sensitive, the concentration of MM3D is sufficient in 0.5 μM to measure Endo-M [A1651] activity (b). Next, addition of DMSO to hinder the enzymatic reaction suppresses enzyme activity in a concentration-dependent manner (c). These results indicate that MM3D is suitable for screening for potential inhibitors of ENGase.

## Related Product

**endo-β-N-Acetylglucosaminidase (= Endo-M)**  
 Recombinant: from *Mucor hiemalis* expressed in *Candida boidinii*

1vial [A1651]

For further information please refer to our website at [www.TCIchemicals.com](http://www.TCIchemicals.com).

TCI glycoscience



## Ordering and Customer Service

### TCI AMERICA

Tel : 800-423-8616 / 503-283-1681  
 Fax : 888-520-1075 / 503-283-1987  
 E-mail : Sales-US@TCIchemicals.com

### TCI EUROPE N.V.

Tel : +32 (0)3 735 07 00  
 Fax : +32 (0)3 735 07 01  
 E-mail : Sales-EU@TCIchemicals.com

### TCI Deutschland GmbH

Tel : +49 (0)6196 64053-00  
 Fax : +49 (0)6196 64053-01  
 E-mail : Sales-DE@TCIchemicals.com

### Tokyo Chemical Industry UK Ltd.

Tel : +44 (0)1865 78 45 60  
 E-mail : Sales-UK@TCIchemicals.com

### 梯希爱(上海)化成工业发展有限公司

Tel : 800-988-0390 / 021-67121386  
 Fax : 021-6712-1385  
 E-mail : Sales-CN@TCIchemicals.com

### Tokyo Chemical Industry (India) Pvt. Ltd.

Tel : 1800 425 7889 / 044-2262 0909  
 E-mail : Sales-IN@TCIchemicals.com

### TOKYO CHEMICAL INDUSTRY CO., LTD.

Tel : +81 (0)3-5640-8878  
 E-mail : globalbusiness@TCIchemicals.com

• Chemicals itemized in this brochure are for research and testing use only. Please avoid use other than by chemically knowledgeable professionals. • Information such as listed products and its specifications and so on are subject to change without prior notice. • The contents may not be reproduced or duplicated in whole or in part without permission of Tokyo Chemical Industry Co., Ltd.