

Europium Fluorophore DTBTA-Eu³⁺-labeled Secondary Antibodies and Streptavidin

**Highly-sensitive Detection Probes
by Time-resolved Fluorometry**

Advantages

No cross talk of excitation light

- $\lambda_{\text{ex, max}} = 335 \text{ nm}$
- $\lambda_{\text{em, max}} = 616 \text{ nm}$

Sharpened emission spectrum



Large Stokes shift (the difference in wavelength between positions of the band maxima of the absorption and emission spectra)

Stable fluorescence in various aqueous buffers

Available in Tris, TE, PBS, etc. for wide use


Long fluorescent life time ($\tau = 1.02 \text{ ms}$)

Time-resolved fluorometric measurement can remove background fluorescence from the sample matrix and often gives detectability better than one order of magnitude compared to those of conventional fluorometric assays.




Comparison of secondary antibody conjugated to DTBTA-Eu³⁺ or FITC

**Time-resolved fluorometric measurement can remove background fluorescence!
To obtain a high SN ratio**



Goat Anti-Mouse IgG DTBTA-Eu³⁺ Conjugate [G0505]



Goat Anti-Mouse IgG FITC Conjugate [G0406]

<Assay condition>

Dilute the Mouse IgG to each concentration. Coat 96-well plates with diluted Mouse IgG. Block the plates with BSA/TBST. Incubate with Goat Anti-Mouse IgG Conjugates prepared from DTBTA or FITC at 2.5 µg/mL. After incubation, measure the fluorescence intensity on a plate reader. DTBTA-Eu³⁺; excitation=340 nm, emission=620 nm. Lag Time : 450 µsec. FITC; excitation=485 nm, emission=520 nm.

Europium Fluorophore DTBTA-Eu³⁺-labeled Secondary Antibodies and Streptavidin

DTBTA-Eu³⁺-labeled Secondary Antibodies and Streptavidin

Goat Anti-Mouse IgG DTBTA-Eu ³⁺ Conjugate	0.1mg/vial [G0505]
Goat Anti-Rabbit IgG DTBTA-Eu ³⁺ Conjugate	0.1mg/vial [G0506]
Streptavidin DTBTA-Eu ³⁺ Conjugate	0.1mg/vial [S0993]

Related Products

Europium Fluorescent Labeling Reagent and Related Products

ATBTA-Eu ³⁺	10mg [A2083]
Cyanuric Chloride	25g / 500g [C0460]
Anti-DTBTA-Eu ³⁺ Rabbit Antiserum	0.5mL [A2181]
Anti-DTBTA-Eu ³⁺ Rabbit Polyclonal Antibody	0.5mL [A2239]

Secondary Antibodies

Goat Anti-Mouse IgG	1mg/vial [G0386]
Goat Anti-Mouse IgG Biotin Conjugate	0.1mg/vial [G0387]
Goat Anti-Mouse IgG FITC Conjugate	0.1mg/vial [G0406]
Goat Anti-Mouse IgG R-PE Conjugate	0.1mg/vial [G0569]
Goat Anti-Mouse IgG HRP Conjugate	0.1mg/vial [G0407]
Goat Anti-Mouse IgM	1mg/vial [G0408]
Goat Anti-Mouse IgM Biotin Conjugate	0.1mg/vial [G0432]
Goat Anti-Mouse IgM FITC Conjugate	0.1mg/vial [G0453]
Goat Anti-Mouse IgM HRP Conjugate	0.1mg/vial [G0417]
Goat Anti-Rabbit IgG	1mg/vial [G0388]
Goat Anti-Rabbit IgG Biotin Conjugate	0.1mg/vial [G0389]
Goat Anti-Rabbit IgG FITC Conjugate	0.1mg/vial [G0452]
Goat Anti-Rabbit IgG R-PE Conjugate	0.1mg/vial [G0577]
Goat Anti-Rabbit IgG HRP Conjugate	0.1mg/vial [G0418]
Mouse Anti-Human IgG Fc	0.1mg/vial [M2977]
Sheep Anti-Chicken IgY	1mg/vial [S0998]
Sheep Anti-Chicken IgY Biotin Conjugate	0.1mg/vial [H1619]

Streptavidins

Streptavidin from <i>Streptomyces avidinii</i>	1mg/vial [S0951]
Streptavidin FITC Conjugate	0.1mg/vial [S0966]
Streptavidin R-PE Conjugate	0.1mg/vial [T3885]
Streptavidin HRP Conjugate	0.1mg/vial [S0972]
Streptavidin Maleimide Conjugate	0.5mg/vial [T3531]

For further information please refer to our website at www.TCIchemicals.com. >>>

TCI DTBTA or antibody



Ordering and Customer Service

TCI AMERICA

Tel : 800-423-8616 / 503-283-1681
Fax : 888-520-1075 / 503-283-1987
E-mail : Sales-US@TCIchemicals.com

TCI EUROPE N.V.

Tel : +32 (0)3 735 07 00
Fax : +32 (0)3 735 07 01
E-mail : Sales-EU@TCIchemicals.com

TCI Deutschland GmbH

Tel : +49 (0)6196 64053-00
Fax : +49 (0)6196 64053-01
E-mail : Sales-DE@TCIchemicals.com

Tokyo Chemical Industry UK Ltd.

Tel : +44 (0)1865 78 45 60
E-mail : Sales-UK@TCIchemicals.com

梯希爱(上海)化成工业发展有限公司

Tel : 800-988-0390 / 021-67121386
Fax : 021-6712-1385
E-mail : Sales-CN@TCIchemicals.com

Tokyo Chemical Industry (India) Pvt. Ltd.

Tel : 1800 425 7889 / 044-2262 0909
E-mail : Sales-IN@TCIchemicals.com

TOKYO CHEMICAL INDUSTRY CO., LTD.

Tel : +81 (0)3-5640-8878
E-mail : globalbusiness@TCIchemicals.com

• Chemicals itemized in this brochure are for research and testing use only. Please avoid use other than by chemically knowledgeable professionals. • Information such as listed products and its specifications and so on are subject to change without prior notice. • The contents may not be reproduced or duplicated in whole or in part without permission of Tokyo Chemical Industry Co., Ltd.