LIFE SCIENCE



Pre-Weighed Biotinylation Reagents

Biotin-LC-LC-NHS (2mg×5) Biotin-PEG₂-NHS (2mg×5) 1set [**B6096]** 1set [**B6097**]



Advantages

B6096 and **B6097** contain both a linker and an *N*-hydroxysuccinimidyl ester moiety, and easily react with amino group $(-NH_2)$ of proteins. Target samples can be biotinylated without weighing of the products during the preparation. **B6096** and **B6097** include 5 reagent vials, each containing 2 mg of respective reagent. The pre-aliquoted packaging prevents decline of the reagent reactivity over time by eliminating the need for repetitive opening of the vial.

Applications

Preparation :

Use of a 10 mM biotinylation solution is recommended. In order to efficiently biotinylate a sample, biotinylation solution should be used at a 15-fold molar excess over the amount of amine-containing protein. Make sure to calculate the 10 mM biotinylation solution amount (see example below).

Calculate : A μ L of 10 mM biotinylation solution for biotinylation 2 mg lgG (150,000 M.W.) 2 [mg lgG] x 10⁻³ [g/mg] x 1/150,000 [mol/g] x 15 [fold] = A [μ L of 10 mM biotinylation solution] x 10⁻⁶ [L/ μ L] x 10 [mmol/L] x 10⁻³ [mol/mmol]

A = 20 [μ L of 10 mM biotinylation solution]

Direction for Use :

- 1. Bring each product to room temperature.
- 2. Dissolve 2 mg of Biotin-LC-LC-NHS [B6096] in 350 μL of DMSO or DMF or 2 mg of Biotin-PEG₂-NHS [B6097] in 400 μL of PBS to prepare a 10 mM biotinylation solution.
- 3. Dissolve the sample (1-10 mg/mL) in an appropriate buffer such as PBS. Do not use buffers including amines (such as Tris).
- Add A μL of 10 mM biotinylation solution to the sample solution and incubate the mixed solution for 30 min at room temperature.
- 5. Remove unreacted and hydrolyzed reagent using desalting column or dialysis methods.

Example for use of **B6096**

The following is an example of the process used for the biotinylation of Mouse Anti-Human IgG Antibody [M2977] using Biotin-LC-LC-NHS [B6096].

Direction for Use :

- 1. Bring 2 mg of B6096 to room temperature and dissolve it in 350 µL of DMSO to prepare 10 mM Biotin-LC-LC-NHS solution.
- 2. Add 10 µL of 10 mM Biotin-LC-LC-NHS solution to 500 µL of M2977 (Conc.: 1 mg/mL Vol.: 0.5 mg).
- 3. Incubate the mixed solution at room temperature for 30 minutes .
- 4. Remove unreacted and hydrolyzed reagent by dialysis methods.

The biotinylation results of Mouse IgG Antibody [M2977] carried out through the above method using B6096 were compared to a similar antibody biotinylated using a product from another company.



The titer of the antibody biotinylated using B6096 was greater than that of the same antibody biotinylated using the product of another company.

Related Products

| DMSO | 25g / 500g [D0798] |
|--|------------------------------------|
| | 5 5 - |
| DMF | 25mL / 100mL / 500mL [D0722] |
| Sulfo-SMCC Sodium Salt | 20mg / 100mg <mark>[S0883]</mark> |
| Biotin-PEG ₂ -NHS | 25mg / 100mg <mark>[\$0955]</mark> |
| Biotin-LC-LC-NHS | 25mg / 100mg <mark>[S0956]</mark> |
| НАВА | 5g / 25g <mark>[H0586]</mark> |
| BSA Maleimide Conjugate (1mg×3) | 1set [B5944] |
| Horseradish Peroxidase Maleimide Conjugate (0.5mg×3) | 1set [H1621] |
| Streptavidin Maleimide Conjugate | 0.5mg/vial [T3531] |
| Streptavidin from Streptomyces avidinii | 1mg/vial [S0951] |
| Streptavidin FITC Conjugate | 0.1mg/vial [S0966] |
| Streptavidin HRP Conjugate | 0.1mg/vial [S0972] |
| Streptavidin DTBTA-Eu³⁺ Conjugate | 0.1mg/vial [S0993] |

For further information please refer to our website at www.TCIchemicals.com.

Tel

Ordering and Customer Service

TCI AMERICA

:800-423-8616 / 503-283-1681 Tel :888-520-1075 / 503-283-1987 Fax E-mail : Sales-US@TCIchemicals.com

TCI EUROPE N.V. :+32 (0)3 735 07 00 Tel Fax :+32 (0)3 735 07 01

E-mail : Sales-EU@TCIchemicals.com

TCI Deutschland GmbH Tel

Fax

:+49 (0)6196 64053-00 +49 (0)6196 64053-01 Fax E-mail : Sales-DE@TCIchemicals.com E-mail : Sales-CN@TCIchemicals.com

Tokyo Chemical Industry UK Ltd. Tokyo Chemical Industry (India) Pvt. Ltd. :+44 (0)1865 78 45 60 Tel E-mail : Sales-UK@TCIchemicals.com

梯希爱(上海)化成工业发展有限公司 Tel :800-988-0390 / 021-67121386

E-mail : Sales-IN@TCIchemicals.com TOKYO CHEMICAL INDUSTRY CO., LTD.

: 1800 425 7889 / 044-2262 0909

Tel : +81 (0)3-5640-8878 E-mail : globalbusiness@TCIchemicals.com

• Chemicals itemized in this brochure are for research and testing use only. Please avoid use other than by chemically knowledgeable professionals. • Information such as listed products and its specifications and so on are subject to change without prior notice. • The contents may not be reproduced or duplicated in whole or in part without permission of Tokyo Chemical Industry Co., Ltd.

:021-6712-1385