


Fluorous Chemistry



Fluorous Solvents


Fluorous Compounds

Fluorous Chemistry


Recently, fluorous chemistry has been studied intensively from the perspective of "Green Chemistry", as the products can be readily separated and the solvents used are reusable. The term "fluorous" was introduced as the analogue to the term aqueous, meaning dissolve in fluorocarbon solvents. Although highly fluorinated compounds (fluorous compounds) neither dissolve in common organic solvents nor in water, they dissolve well in fluorous solvents such as perfluoroalkane. Fluorous chemistry utilizes this property and a numerous application of this chemistry has been made.

1. Organic reaction using fluororous solvents


Although fluorous solvents are immiscible with water and common organic solvents, certain fluorous solvents have the properties to form a homogeneous solution with some organic solvents at elevated temperatures. They also have the properties that the boiling points are almost equal to those of the corresponding hydrocarbons regardless of their molecular weight, and the high solubility of many gases in these solvents. Taking the advantage of these properties, Horváth *et al.* accomplished the hydroformylation of olefins using a fluorous rhodium catalyst in perfluoromethylcyclohexane [P0846] and toluene in 1994.¹⁾ This was regarded as the origin of the fluorous chemistry. This reaction uses perfluoromethylcyclohexane and toluene as solvent, which form a biphasic system at room temperature. In this system the fluorous catalyst exists in the fluorous phase and the olefins in the organic phase. However, the two phases form a homogeneous solution when heated. The reaction then proceeds by introduction of carbon monoxide and hydrogen gases. When the reaction is complete and cooled, the two phase system reappears, where the resulting product is dissolved in the toluene phase and the fluorous catalyst in the fluorous phase, thus, making the catalyst and the product easily separable. The biphasic system using a fluorous solvent and an organic solvent is called Fluorous Biphasic System (FBS), and the multiple phase system is called Fluorous Multiphase System (FMS). The advantages of FBS and FMS are that the resulting product and the catalyst can be easily separated simply by separating the fluorous phase from the other phase after the reaction, and that the fluorous phase containing fluorous catalyst can be reusable after separation.





Zhu also reported the synthesis of carboxylic ester from methacrylic acid and decanol using p-toluenesulfonic acid in perfluoro(2-butyltetrahydrofuran) [P0867].²⁾ Although methacrylic acid, decanol and the acid catalyst dissolve in a fluorous solvent upon heating, the water generated in the reaction does not dissolve in the fluorous solvent. When the reaction mixture is cooled, the resulting carboxylic ester floats on the water phase, and the desired product can be separated.




A numerous oxidation reactions in biphasic system with fluorous solvents and organic solvents have also been studied with oxygen molecule. Knochel *et al.* have reported the oxidation of aldehydes, olefins and sulfides in the presence of a nickel complex catalyst with a fluorous β -diketone as ligand.³⁾ For the oxidation of aldehydes, perfluorodecalin [P0837] and toluene were used as solvents, and this system also was found to form a homogeneous solution upon heating. After the reaction was over and cooled to room temperature, the catalyst staying in the fluorous phase and the product in the organic phase were easily separated. Due to their strong solubility, fluorous solvents are suitable for many reactions that requires to gases reagents.




Various oxidation reaction of sulfides and olefins have also been studied similarly in the presence of isobutylaldehyde.³⁾ The solvents used in these reactions were perfluorooctyl bromide [H0946] and toluene, and this solvent system also formed a homogeneous solution upon heating.




Pozzi *et al.* have also reported an epoxidation of olefins using molecular oxygen catalyzed by fluorous porphyrin-cobalt complex in the presence of isobutylaldehydes.⁴⁾ This reaction was carried out in biphasic system of perfluorohexane [P0839] and acetonitrile by stirring the mixture at room temperature. When the reaction was complete, the catalyst and the product were separated as usual, and the fluorous phase containing the catalyst was reused.



The Wacker oxidation reaction using perfluoroctyl bromide [H0946] as fluorous solvent has also been reported.⁵⁾ Perfluoroctyl bromide and benzene form a homogeneous solution when heated. After the reaction is complete and cooled, the product is separated from the palladium catalyst complexed with fluorous β -diketone. The fluorous phase can be reused after separation.



Horváth and Gladysz *et al.* have reported a hydroboration in perfluoromethylcyclohexane [P0846] and toluene using a rhodium complex catalyst with fluorous ligands.⁶⁾ After the reaction was complete, the product was separated, and the fluorous phase containing the catalyst was reused.



2. Application to the synthesis of sugar chains and Combinatorial Chemistry

Curran *et al.* have introduced the use of fluorous substituents (fluorous tags) into non-fluorous substrates and the synthesis of isoxazoline by using this fluorous compound.⁷⁾ After the reaction,

the fluorous product was separated by extraction with dichloromethane, water, and perfluorohexane. Following this report, a numerous applications of this fluorous chemistry have been made in combinatorial chemistry.⁸⁾ And, Inazu *et al.* have applied this chemistry to the synthesis of oligosaccharide.⁹⁾ In this reaction, the fluorous tag was first introduced into the sugar molecule, and then glycoxylation followed. The desired oligosaccharide thus obtained was extracted with an organic solvent, water, and perfluorohexane.

As shown by the aforementioned examples, fluorous chemistry introduced by Horváth *et al.* has widely been applied in many areas of synthetic chemistry. Utilizing this chemistry, it is possible to isolate the desired product easily from the catalyst and the fluorous solvents. Furthermore, the separated fluorous solvents and the catalysts can be reused. A great deal of studies have been made on this subject, especially, because of its usefulness in term of Green Chemistry. It is also expected that this chemistry will be widely used in the application to the combinatorial chemistry where many compounds are handled at multiple steps.

References

- I. T. Horváth, J. Rábai, *Science* **1994**, 266, 72.
- D.-W. Zhu, *Synthesis* **1993**, 953.
- I. Klement, H. Lütjens, P. Knochel, *Angew. Chem. Int. Ed. Engl.* **1997**, 36, 1454.
- G. Pozzi, F. Montanari, S. Quici, *Chem. Commun.* **1997**, 69.
- B. Betzemeier, F. Lhermitte, P. Knochel, *Tetrahedron Lett.* **1998**, 39, 6667.
- J. J. J. Juliette, I. T. Horváth, J. A. Gladysz, *Angew. Chem. Int. Ed. Engl.* **1997**, 36, 1610; J. J. J. Juliette, D. Rutherford, I. T. Horváth, J. A. Gladysz, *J. Am. Chem. Soc.* **1999**, 121, 2696.
- A. Studer, S. Hadida, R. Ferritto, S. -Y. Kim, P. Jeger, P. Wipf, D. P. Curran, *Science* **1997**, 275, 823.
- D. P. Curran, S. Hadida, *J. Am. Chem. Soc.* **1996**, 118, 2531; D. P. Curran, M. Hoshino, *J. Org. Chem.* **1996**, 61, 6480; D. P. Curran, *Angew. Chem. Int. Ed.* **1998**, 37, 1174; D. P. Curran, Z. Luo, *J. Am. Chem. Soc.* **1999**, 121, 9069; Q. Zhang, Z. Luo, D. P. Curran, *J. Org. Chem.* **2000**, 65, 8866; Z. Luo, Q. Zhang, Y. Oderoatoshi, D. P. Curran, *Science* **2001**, 291, 1766; S. Darses, M. Puchéault, J.-P. Genét, *Eur. J. Org. Chem.* **2001**, 1121.
- T. Miura, Y. Hirose, M. Ohmae, T. Inazu, *Org. Lett.* **2001**, 3, 3947; T. Miura, T. Inazu, *Tetrahedron Lett.* **2003**, 44, 1819.
- Review
K. Ishihara, H. Yamamoto, *Kagaku To Kogyo* **2001**, 54, 1061; K. Ishihara, *Kagaku To Kogyo* **2002**, 55, 865; I. Ryu, H. Matsubara, *Kagaku* **2002**, 57(5), 20; S. Takeuchi, Y. Nakamura, *Kagaku* **2002**, 57(6), 16; K. Mikami, H. Matsuzawa, *Kagaku* **2002**, 57(7), 22; K. Ishihara, H. Yamamoto, *Kagaku* **2002**, 57(8), 30.

Keywords: fluorous chemistry, fluorous solvents, environmentally-friendly solvents

Fluorous Solvents





















E0485 5g 25g	P1755 1g	T1012 25mL	H1013 25g 500g	O0292 10g 50g
<chem>CF3(CF2)7CF3</chem> Perfluorononane CAS RN: 375-96-2	<chem>CF3(CF2)10CF3</chem> Perfluorododecane CAS RN: 307-59-5	<chem>CF3CF2CF2CF(CF3)2</chem> Perfluoroisohexane CAS RN: 355-04-4	 1,1,2,2,3,3,4-Heptafluorocyclopentane CAS RN: 15290-77-4	 Perfluorocyclopentene CAS RN: 559-40-0
P0846 25g 100g	P1420 25g	P0837 25g	H0085 5g 25g 250g	P0856 5g 25g
 Perfluoromethylcyclohexane CAS RN: 355-02-2	 Perfluoro(1,3-dimethylcyclohexane) CAS RN: 335-27-3	 Perfluorodecalin CAS RN: 306-94-5	 Perfluorobenzene CAS RN: 392-56-3	 Perfluorotoluene CAS RN: 434-64-0
H0946 5g 25g	D4484 25g 500g	P0867 25g	P1348 5g 25g	P0074 5g 25g 100g
<chem>CF3(CF2)7Br</chem> Perfluoro-n-octyl Bromide CAS RN: 423-55-2	<chem>CF3CF2CFCF(OCH3)CF3</chem> 1,1,2,2,3,4,5,5,5-Decafluoro-3-methoxy-4-(trifluoromethyl)pentane CAS RN: 132182-92-4	 Perfluoro-(2-butyltetrahydrofuran) CAS RN: 335-36-4	<chem>CF3CF2N(CF2CF3)CF2CF3</chem> Perfluorotriethylamine CAS RN: 359-70-6	<chem>CF3(CF2)3N(CF2)3CF3</chem> Perfluorotributylamine CAS RN: 311-89-7
P1051 25g				
 Perfluorotriamylamine CAS RN: 338-84-1				

Fluorous Compounds

P1102 5g 25g	T2496 5g	U0076 5g	H0846 5g 25g	N0601 5g 25g
<chem>CF3(CF2)5CH=CH2</chem> (Perfluoroethyl)ethylene CAS RN: 25291-17-2	<chem>CF3(CF2)4CHF2</chem> 1 <i>H</i> -Tridecafluorohexane CAS RN: 355-37-3	<chem>CF3(CF2)3CHF2</chem> 1 <i>H</i> -Undecafluoropentane CAS RN: 375-61-1	<chem>CF3(CF2)7CH=CH2</chem> (Perfluoro-n-octyl)ethylene CAS RN: 21652-58-4	<chem>CF3(CF2)3CH=CH2</chem> (Perfluorobutyl)ethylene CAS RN: 19430-93-4
D4128 5g 25g	D1101 25g 100g	D2891 5g	E0239 10g	N1040 5g 25g 100g
<chem>CHF2CH2OH</chem> 2,2-Difluoroethanol CAS RN: 359-13-7	<chem>H(CF2)6CH2OH</chem> 1,1,7-Trihydroperfluoroheptanol CAS RN: 335-99-9	<chem>HOCH2(CF2)6CH2OH</chem> 2,2,3,3,4,4,5,5,6,6,7,7-Dodecafluoro-1,8-octanediol CAS RN: 90177-96-1	<chem>H(CF2)10CH2OH</chem> 1 <i>H</i> ,1 <i>H</i> ,11 <i>H</i> -Eicosfluoro-1-undecanol CAS RN: 307-70-0	<chem>CF3(CF2)3(CH2)3OH</chem> 4,4,5,5,6,6,7,7,7-Nonafluoro-1-heptanol CAS RN: 83310-97-8

U0074 5g <chem>CF3(CF2)4CH2OH</chem> 1H,1H-Undecafluoro-1-hexanol CAS RN: 423-46-1	T3258 5g 25g <chem>CF3(CF2)5(CH2)3OH</chem> 1H,1H,2H,2H,3H,3H-Tridecafluoro-1-nonanol CAS RN: 80806-68-4	H0845 25g 250g <chem>CF3(CF2)7CH2CH2OH</chem> 2-(Perfluoro-n-octyl)ethanol CAS RN: 678-39-7	H1232 5g <chem>CF3(CF2)7CH2OH</chem> 1H,1H-Perfluoro-1-nonanol CAS RN: 423-56-3	H1349 5g 25g <chem>CF3 CH3-C(O)- CF3</chem> 1,1,1,3,3,3-Hexafluoro-2-methyl-2-propanol CAS RN: 1515-14-6
H0548 5g 25g <chem>CF3CF2CF2CH2OH</chem> 1H,1H-Heptafluoro-1-butanol CAS RN: 375-01-9	H1233 1g 5g <chem>HOCH2(CF2)8CH2OH</chem> 1H,1H,10H,10H-Hexadecafluoro-10-decanediol CAS RN: 754-96-1	H1035 25g <chem>H(CF2)8CH2OH</chem> 1H,1H,9H-Hexadecafluoro-1-nonanol CAS RN: 376-18-1	H1279 5g 25g <chem>F3C HO-C(C(O)- CF3)- CF3</chem> Perfluoropinacol CAS RN: 918-21-8	H0649 25g <chem>CF3CHFCF2CH2OH</chem> 2,2,3,4,4,4-Hexafluoro-1-butanol CAS RN: 382-31-0
H0746 1g 5g 25g <chem>HOCH2(CF2)3CH2OH</chem> 2,2,3,3,4,4-Hexafluoro-1,5-pentanediol CAS RN: 376-90-9	H0424 25g 100g 500g <chem>CF3 C(O)- CF3</chem> 1,1,1,3,3,3-Hexafluoro-2-propanol CAS RN: 920-66-1	N0814 5g <chem>CF3(CF2)8CH2OH</chem> 1H,1H-Perfluoro-1-decanol CAS RN: 307-37-9	N0692 1g 5g 25g <chem>CF3 F3C-C(O)- CF3</chem> Perfluoro-tert-butanol CAS RN: 2378-02-1	N0600 5g 25g <chem>CF3(CF2)3CH2CH2OH</chem> 2-(Perfluorobutyl)ethanol CAS RN: 2043-47-2
N0810 1g 5g 25g <chem>CF3(CF2)3CH2OH</chem> (Perfluorobutyl)methanol CAS RN: 355-28-2	O0294 5g 25g <chem>HOCH2(CF2)4CH2OH</chem> 2,2,3,3,4,4,5-Octafluoro-1,6-hexanediol CAS RN: 355-74-8	O0114 25g 100g 500g <chem>CHF2(CF2)3CH2OH</chem> 2,2,3,3,4,4,5,5-Octafluoro-1-pentanol CAS RN: 355-80-6	P0904 5g 25g <chem>CF3(CF2)6CH2OH</chem> 1H,1H-Perfluoro-1-octanol CAS RN: 307-30-2	P0845 25g <chem>CF3CF2CH2OH</chem> 1H,1H-Pentafluoro-1-propanol CAS RN: 422-05-9
T1701 5g 25g <chem>CF3(CF2)5CH2OH</chem> 1H,1H-Perfluoro-1-heptanol CAS RN: 375-82-6	T2528 5g 25g <chem>CF3(CF2)5CH2CH2OH</chem> 2-(Perfluorohexyl)ethanol CAS RN: 647-42-7	T0435 25g 100g 500g <chem>CF3CH2OH</chem> 2,2,2-Trifluoroethanol CAS RN: 75-89-8	T3381 1g 5g <chem>CF3(CF2)10CH2OH</chem> 1H,1H-Tricosafluoro-1-dodecanol CAS RN: 423-65-4	T0101 25g 100g 500g <chem>CHF2CF2CH2OH</chem> 2,2,3,3-Tetrafluoro-1-propanol CAS RN: 76-37-9
Fluorous Ethers		B1293 1g 5g <chem>CF3CH2OCH2CF3</chem> 2,2,2-Trifluoroethyl Ether CAS RN: 333-36-8	H1773 100g <chem>CF3(CF2)2OCH3</chem> 1,1,1,2,2,3,3-Heptafluoro-3-methoxypropane CAS RN: 375-03-1	B4169 5g <chem>BrCF2CF2O-C=CF2</chem> 2-Bromotetrafluoroethyl Trifluorovinyl Ether CAS RN: 85737-06-0
C2862 5g <chem>Cl F(F)C(F)OC(F)CH3</chem> 2-Chloro-1,1,2-trifluoroethyl Ethyl Ether CAS RN: 310-71-4	C0853 5g <chem>Cl F(F)C(F)OC(F)CH3</chem> 2-Chloro-1,1,2-trifluoroethyl Methyl Ether CAS RN: 425-87-6	D4484 25g 500g <chem>OCH3 CF3CF2CFCF(CF3)CF3</chem> 1,1,1,2,2,3,4,5,5-Decafluoro-3-methoxy-4-(trifluoromethyl)pentane CAS RN: 132182-92-4	D4472 1g 5g <chem>FC(F)C(F)OC(F)F</chem> Difluoromethyl 2,2,3-Tetrafluoropropyl Ether CAS RN: 35042-99-0	H1507 5g 25g <chem>CF3 F(F)C(F)OC(F)CH3</chem> 1,1,2,3,3,3-Hexafluoropropyl Methyl Ether CAS RN: 382-34-3
E1020 5g <chem>CF3 F(F)C(F)OC(F)CH3</chem> Ethyl 1,1,2,3,3,3-Hexafluoropropyl Ether CAS RN: 380-34-7	E0528 25g 500g <chem>CF3(CF2)3OCH2CH3 CF3>CFCF2OCH2CH3</chem> (mixture of isomers) Ethyl Nonafluorobutyl Ether (mixture of isomers) CAS RN: 813458-04-7	E1019 5g 25g <chem>CF3 F(F)C(F)OC(F)CH3</chem> Ethyl 1,1,2,2-Tetrafluoroethyl Ether CAS RN: 512-51-6	F0691 5g <chem>CF3 CF3>CHOCH2F</chem> Fluoromethyl 1,1,3,3,3-Hexafluoroisopropyl Ether CAS RN: 28523-86-6	P1226 5g <chem>CF3 CF3CF2CF2OCFCF2O-C=CF2</chem> 2-(Heptafluoropropoxy)-hexafluoropropyl Trifluorovinyl Ether CAS RN: 1644-11-7

Fluorous Chemistry





































H1610	5g 25g	H1611	5g	H1524	5g 25g	M1345	25g 500g	M2500	1g
 CF ₃ O(CF ₂) ₃ O—C=CF ₂ 1,1,2,2,3,3-Hexafluoro-1-(trifluoromethoxy)-3-[(1,2,2-trifluorovinyl)oxy]propane CAS RN: 40573-09-9		 CF ₂ =C—O(CF ₂) ₃ O—C=CF ₂ 1,1,2,2,3,3-Hexafluoro-1,3-bis[(1,2,2-trifluorovinyl)oxy]propane CAS RN: 13846-22-5		 CF ₃ CF ₃ —CHOCH ₃ Isoindoklon CAS RN: 13171-18-1		 CF ₃ (CF ₂) ₃ OCH ₃ F ₃ C—CFCF ₂ OCH ₃ (mixture of isomers) Methyl Nonafuorobutyl Ether CAS RN: 219484-64-7		 CF ₃ CF ₃ —C(F) ₂ —OCH ₃ Methyl 2,2,3,3,3-Pentafluoropropyl Ether CAS RN: 378-16-5	
M2514	25g	O0422	5g 25g	P1224	10g	H1624	25g		
 Methyl 1,2,2-Tetrafluoroethyl Ether CAS RN: 425-88-7		 1H,5H-Octafluoropentyl 1,1,2,2-Tetrafluoroethyl Ether CAS RN: 16627-71-7		 CF ₃ CF ₂ CF ₂ O—C=CF ₂ Perfluoropropoxyethylene CAS RN: 1623-05-8		 CF ₃ CF ₂ CF ₂ O—C—CF ₂ OCHFCF ₃ 1,1,2,2,3,3-Heptafluoro-3-[[1,1,2,3,3-hexafluoro-3-[(1,1,1,2,3,3-hexafluoro-3-[(1,2,2,2-tetrafluoroethoxy)propan-2-yl]oxy]propan-2-yl]oxy]propane CAS RN: 3330-14-1			
H1625	25g	T3069	5g 25g	I1044	25g				
 1,1,1,2,2,3,3-Heptafluoro-3-[[1,1,1,2,3,3-hexafluoro-3-[(1,1,1,2,3,3-hexafluoro-3-[(1,2,2,2-tetrafluoroethoxy)propan-2-yl]oxy]propan-2-yl]oxy]propane CAS RN: 3330-16-3		 1,1,2,2-Tetrafluoroethyl 2,2,3,3-Tetrafluoropropyl Ether CAS RN: 16627-68-2		 CF ₃ CF ₂ CF ₂ O—C—CF ₂ O—CHFCF ₃ 1,1,1,2,4,4,5,7,7,8,10,10,11,13,13,14,14,14,15,15-Icosafluoro-5,8,11-tris(trifluoromethyl)-3,6,9,12-tetraoxapentadecane CAS RN: 26738-51-2					
T3538	25g	T3057	5g 25g	D5223	5g 25g	T3493	5g 25g		
 1,1,1,2,4,4,5,7,7,8,10,10,11,13,13,14,16,16,17,17,18,18,18-Tricosfluoro-5,8,11,14-tetrakis(trifluoromethyl)-3,6,9,12,15-pentaoxaoctadecane CAS RN: 37486-69-4		 1,1,2-Tetrafluoroethyl 2,2,2-Trifluoroethyl Ether CAS RN: 406-78-0		 NC(CF ₂) ₅ O—C=CF ₂ 2,2,3,3,4,4,5,5,6,6-Decafluoro-6-[(1,2,2-trifluorovinyl)oxy]-hexanenitrile CAS RN: 120903-40-4		 NCCF ₂ CF ₂ O—C—CF ₂ O—C=CF ₂ 2,2,3-Tetrafluoro-3-[(1,2,2-trifluorovinyl)oxy]propan-2-yl]oxy]propionitrile CAS RN: 69804-19-9			
Fluorous Ketons		B1240	5g 25g	C0993	1g	D1729	5g	P1363	5g
		 1-Bromo-3,3,3-trifluoroacetone CAS RN: 431-35-6		 Chloropentafluoroacetone Monohydrate CAS RN: 6984-99-2		 CF ₃ (CF ₂) ₂ —C(=O)—CH ₂ —C(=O)—CH ₃ 2,2-Dimethyl-6,6,7,7,8,8-heptafluoro-3,5-octanedione CAS RN: 17587-22-3		 CF ₃ (CF ₂) ₄ —C(=O)—CH ₂ CH ₃ Ethyl Undecafluoroamyl Ketone CAS RN: 383177-55-7	
H0425	5g 25g	H0476	5g 25g	P1452	5g	U0071	5g	N1038	5g 25g
 Hexafluoroacetone Hydrate CAS RN: 10057-27-9		 Hexafluoroacetylacetone CAS RN: 1522-22-1		 Methyl Pentadecafluoroheptyl Ketone CAS RN: 754-85-8		 CF ₃ (CF ₂) ₄ —C(=O)—CH ₃ Methyl Undecafluoroamyl Ketone CAS RN: 2708-07-8		 CF ₃ —C(F) ₂ —C(=O)—CF ₂ CF ₃ Perfluoroethyl Perfluoroisopropyl Ketone CAS RN: 756-13-8	
D1631	1g 5g	Fluorous Carboxylic Acids		T0431	25g 100g 500g	P1125	25g 100g	H0024	25g 100g
		 Decafluorobenzophenone CAS RN: 853-39-4		 Trifluoroacetic Acid CAS RN: 76-05-1		 CF ₃ CF ₂ CF ₂ —C(=O)—OH Pentafluoropropionic Acid CAS RN: 422-64-0		 CF ₃ CF ₂ CF ₂ —C(=O)—OH Heptafluorobutyric Acid CAS RN: 375-22-4	
N0605	5g 25g	U0067	5g 25g	T1545	5g 25g	H0843	5g 25g	N0607	5g
 Nonafluorovaleric Acid CAS RN: 2706-90-3		 Undecafluorohexanoic Acid CAS RN: 307-24-4		 CF ₃ (CF ₂) ₅ —C(=O)—OH Tridecafluoroheptanoic Acid CAS RN: 375-85-9		 CF ₃ (CF ₂) ₇ —C(=O)—OH Heptadecafluorononanoic Acid CAS RN: 375-95-1		 CF ₃ (CF ₂) ₈ —C(=O)—OH Nonadecafluorodecanoic Acid CAS RN: 335-76-2	

H1234 1g Heneicosfluoroundecanoic Acid CAS RN: 2058-94-8	T2492 1g 5g Tricosfluorododecanoic Acid CAS RN: 307-55-1	H1502 1g 2H,2H,3H,3H-Heptadecafluoroundecanoic Acid CAS RN: 34598-33-9	T1621 5g 25g Tetrafluorosuccinic Acid CAS RN: 377-38-8	H0658 5g 25g Hexafluoroglutaric Acid CAS RN: 376-73-8
O0260 5g 25g Octafluoroadipic Acid CAS RN: 336-08-3	D2465 5g 25g Dodecafluorosuberic Acid CAS RN: 678-45-5	H0892 5g 25g Hexadecafluorosebatic Acid CAS RN: 307-78-8	P1004 1g Pentafluorophenylacetic Acid CAS RN: 653-21-4	Fluorous Carboxylic Anhydrides
T2478 1g 5g Tetrafluorosuccinic Anhydride CAS RN: 699-30-9	H0745 5g 25g Hexafluoroglutaric Anhydride CAS RN: 376-68-1	T0433 20mL 100mL 400mL Trifluoroacetic Anhydride CAS RN: 407-25-0	P0566 5g 25g Pentafluoropropionic Anhydride CAS RN: 356-42-3	H0337 10g Heptafluorobutyric Anhydride CAS RN: 336-59-4
D4164 5g Difluoroacetic Anhydride CAS RN: 401-67-2	Fluorous Carboxylic Halides		H0508 5g 25g Heptafluorobutryl Chloride CAS RN: 375-16-6	U0075 5g Undecafluorohexanoyl Fluoride CAS RN: 355-38-4
Fluorous Carboxylic Esters		D3589 1g 5g Dimethyl Hexafluoroglutamate CAS RN: 1513-62-8	D3590 1g 5g Dimethyl Octafluoroadipate CAS RN: 3107-98-0	D3588 1g 5g Dimethyl Tetrafluorosuccinate CAS RN: 356-36-5
D2498 5g 25g Ethyl Difluoroacetate CAS RN: 454-31-9	E1018 25g Ethyl 4,4-Difluoroacetoacetate CAS RN: 352-24-9	E0547 1g Ethyl 3-Ethoxy-2,2-difluoro-3-hydroxypropanoate CAS RN: 141546-97-6	H1038 5g Ethyl Heptadecafluorononanoate CAS RN: 30377-52-7	H0594 5g 25g Ethyl Heptafluorobutyrate CAS RN: 356-27-4
N0689 5g Ethyl Nonfluorovalerate CAS RN: 424-36-2	E1022 5g 25g Ethyl 5H-Octafluorovalerate CAS RN: 2795-50-8	P1062 5g Ethyl Pentafluoropropionylacetate CAS RN: 663-35-4	T0432 25g 100g 500g Ethyl Trifluoroacetate CAS RN: 383-63-1	T1697 25g 2,2,2-Trifluoroethyl Trifluoroacetate CAS RN: 407-38-5
E0830 1g 5g Ethyl 4,4,4-Trifluorobutyrate CAS RN: 371-26-6	E0772 5g 25g Ethyl 4,4,4-Trifluorocrotonate CAS RN: 25597-16-4	H1033 5g 25g Methyl Heptafluorobutyrate CAS RN: 356-24-1	M2022 5g 25g Methyl Heptafluoroisobutyrate CAS RN: 680-05-7	M1916 5g Methyl Nonadecafluorodecanoate CAS RN: 307-79-9

Fluorous Chemistry

M1912 5g <chem>CF3(CF2)3C(=O)OC</chem> Methyl Nonfluorovalerate CAS RN: 13038-26-1	P1453 5g <chem>CF3(CF2)6C(=O)OC</chem> Methyl Pentadecafluoro-octanoate CAS RN: 376-27-2	M1917 5g 25g <chem>CF3(CF2)10C(=O)OC</chem> Methyl Tricosafauro-dodecanoate CAS RN: 56554-52-0	M1914 5g 25g <chem>CF3(CF2)5C(=O)OC</chem> Methyl Tridecafluoro-heptanoate CAS RN: 14312-89-1	M2496 1g 5g <chem>CC(C(F)(F)C)C(=O)OC</chem> Methyl 2-(Trifluoromethyl)-3,3,3-trifluoropropionate CAS RN: 360-54-3
M1913 5g <chem>CF3(CF2)4C(=O)OC</chem> Methyl Undecafluoro-hexanoate CAS RN: 424-18-0	M2030 5g <chem>CC(C(F)(F)C)C(=O)OC</chem> Methyl 2,5-Bis(trifluoromethyl)-3,6-dioxaundecafluorononanoate (mixture of isomers) CAS RN: 26131-32-8	B5785 1g 5g <chem>CC=COC(=O)CC(C(F)(F)C)COC(=O)C=CH</chem> 1,6-Bis(acryloyloxy)-2,2,3,4,4,5,5-octafluorohexane CAS RN: 2264-01-9	A1330 10g <chem>CC=COC(=O)CC(C(F)(F)C)COC(=O)C=CH</chem> 1H,1H,2H,2H-Heptadecafluorodecyl Acrylate CAS RN: 27905-45-9	T3451 5g 25g <chem>CC=COC(=O)CC(C(F)(F)C)COC(=O)C=CH</chem> 1H,1H,2H,2H-Tridecafluoro-n-octyl Acrylate (stabilized with HQ + MEHQ) CAS RN: 17527-29-6
T3414 5g 25g 100g <chem>CC(C(F)(F)C)C(=O)OC</chem> 1H,1H,2H,2H-Tridecafluoro-n-octyl Methacrylate (stabilized with HQ + MEHQ) CAS RN: 2144-53-8	H1703 5g <chem>CC=COC(=O)OC</chem> 1H,1H,2H,2H-Heptadecafluorodecyl Methacrylate (stabilized with MEHQ) CAS RN: 1996-88-9	N0977 25g <chem>CC=COC(=O)OC</chem> 1H,1H,2H,2H-Nonafluorohexyl Acrylate CAS RN: 52591-27-2	N1014 5g 25g <chem>CC(C(F)(F)C)C(=O)OC</chem> 1H,1H,2H,2H-Nonafluorohexyl Methacrylate CAS RN: 1799-84-4	O0318 5g 25g <chem>CC=COC(=O)OC</chem> 1H,1H,5H-Octafluoropentyl Acrylate CAS RN: 376-84-1
00481 5g 25g <chem>CC(C(F)(F)C)C(=O)OC</chem> 1H,1H,5H-Octafluoropentyl Methacrylate CAS RN: 355-93-1	P1754 1g 25g <chem>CC=COC(=O)OC</chem> 1H,1H-Pentadecafluoro-n-octyl Acrylate CAS RN: 307-98-2	N1107 5g 25g <chem>CC=COC(=O)OC</chem> 1H,1H,2H,2H-Nonafluorohexyl Acrylate CAS RN: 2591-27-2	H1674 5g 25g <chem>CC(C(F)(F)C)C(=O)OC</chem> 2,2,3,3,4,4,4-Heptafluorobutyl Methacrylate CAS RN: 13695-31-3	N1143 5g 25g <chem>CC=COC(=O)OC</chem> 1H,1H,2H-Nonafluorohexyl Methacrylate (stabilized with TBC) CAS RN: 1799-84-4
Fluorous Alkyl Halides		B2333 1g 5g <chem>BrCF2CH=CH2</chem> 3-Bromo-3,3-difluoropropene CAS RN: 420-90-6	B3222 5g 25g <chem>BrCF2CF2CH=CH2</chem> 4-Bromo-3,3,4,4-tetrafluoro-1-butene CAS RN: 18599-22-9	D3572 1g 5g <chem>Br(CF2)6Br</chem> 1,6-Dibromododecafluorohexane CAS RN: 918-22-9
D3573 5g <chem>Br(CF2)4Br</chem> 1,4-Dibromooctafluorobutane CAS RN: 335-48-8	D2804 5g 25g <chem>Cl(CF2)8Cl</chem> 1,8-Dichlorohexadecafluoro-octane CAS RN: 647-25-6	D2333 10g <chem>I(CF2)6I</chem> Dodecafluoro-1,6-diiodohexane CAS RN: 375-80-4	H0844 5g 25g <chem>CF3(CF2)9I</chem> Heneicosafluorodecyl Iodide CAS RN: 423-62-1	D3587 1g 5g <chem>Br(CF2)8Br</chem> 1,8-Dibromohexadecafluoro-octane CAS RN: 812-58-8
H0946 5g 25g <chem>CF3(CF2)7Br</chem> Perfluoro-n-octyl Bromide CAS RN: 423-55-2	H0689 5g <chem>CF3CF2CF2Br</chem> Heptafluoropropyl Bromide CAS RN: 422-85-5	H0596 5g 25g <chem>CF3CF2CF2I</chem> Perfluoropropyl Iodide CAS RN: 754-34-7	N0499 25g 100g 500g <chem>CF3CF2CF2CF2I</chem> Nonafluorobutyl Iodide CAS RN: 423-39-2	H1084 5g 25g <chem>CF3(CF2)7CH2CH2I</chem> 1H,1H,2H,2H-Perfluorodecyl Iodide CAS RN: 2043-53-0
D2329 5g 25g <chem>I(CF2)4I</chem> Octafluoro-1,4-diiodobutane CAS RN: 375-50-8	P1753 5g <chem>CF3(CF2)6Br</chem> Pentadecafluoroheptyl Bromide CAS RN: 375-88-2	T2482 1g 5g <chem>CF3(CF2)5CH2I</chem> 1H,1H-Tridecafluoroheptyl Iodide CAS RN: 212563-43-4	T2479 5g 25g <chem>CF3(CF2)5Br</chem> Tridecafluorohexyl Bromide CAS RN: 335-56-8	T1098 5g 25g <chem>CF3(CF2)5I</chem> Tridecafluorohexyl Iodide CAS RN: 355-43-1

Fluorous Chemistry

T0859 0.1mL  2,4,6-Tris(perfluoropropyl)-1,3,5-triazine CAS RN: 915-76-4	T0858 0.1mL  2,4,6-Tris(pentafluoroethyl)-1,3,5-triazine CAS RN: 858-46-8	T3041 1g 5g  Tris(1,1,3,3,3-hexafluoro-2-propyl) Phosphate CAS RN: 66489-68-7	T3203 5g 25g  TTFFPa CAS RN: 358-63-4	P1134 10g  Tris(1H,1H,5H-octafluoropentyl) Phosphate CAS RN: 355-86-2
T3353 1g 5g  Tris(1,1,3,3,3-hexafluoro-2-propyl) Phosphite CAS RN: 66470-81-3	T3991 5g 25g  Tris(2,2-trifluoroethyl)-Phosphite CAS RN: 370-69-4	T2484 1g 5g  Tris(pentafluorophenyl)-phosphine CAS RN: 1259-35-4	B3428 1g  1,2-Bis[bis(pentafluorophenyl)-phosphino]ethane CAS RN: 76858-94-1	T2313 1g 5g  Tris(pentafluorophenyl)-borane CAS RN: 1109-15-5
T1730 100mg  5,10,15,20-Tetrakis(pentafluorophenyl)porphyrin CAS RN: 25440-14-6	O0319 100mg  Perfluoro-5,10,15,20-tetraphenyl-21H,23H-porphyrin CAS RN: 121399-88-0	P0935 5g 25g  Pentafluorobenzonitrile CAS RN: 773-82-0	P0918 5g 25g  Pentafluoroanisole CAS RN: 389-40-2	P1408 5g 25g  2,3,4,5,6-Pentafluorotoluene CAS RN: 771-56-2
P1188 5g 25g  Pentafluoriodobenzene CAS RN: 827-15-6	B1116 5g 25g  Bromopentafluorobenzene CAS RN: 344-04-7	P0850 25g  Chloropentafluorobenzene CAS RN: 344-07-0	P0861 5g 25g  Pentafluorobenzenethiol CAS RN: 771-62-0	P0922 5g 25g  Pentafluoroaniline CAS RN: 771-60-8
P1228 5g 25g  Pentafluoronitrobenzene CAS RN: 880-78-4	P1904 1g 5g  Pentafluorophenylboronic Acid CAS RN: 1582-24-7	P2231 1g 5g  Pentafluorophenyl 4-Nitrobenzenesulfonate CAS RN: 244633-31-6	P0862 5g 25g  2,3,4,5,6-Pentafluorostyrene CAS RN: 653-34-9	P0925 5g  1-(Pentafluorophenyl)-ethanol CAS RN: 830-50-2
P1242 1g  Pentafluorophenylethoxydimethylsilane CAS RN: 71338-73-3	T1542 1g 5g  2,3,5,6-Tetrafluoro-4-(trifluoromethyl)benzenethiol CAS RN: 651-84-3	T1541 10g  4-Trifluoromethyl-2,3,5,6-tetrafluorobromobenzene CAS RN: 17823-46-0	T1529 1g 5g  4-Aminoheptafluorotoluene CAS RN: 651-83-2	T1983 5g  Perfluoro-p-cresol CAS RN: 2787-79-3
P0926 5g 25g  Pentafluoropyridine CAS RN: 700-16-3	P3155 1g 5g  2,3,4,5,6-Pentafluorobenzylphosphonic Acid CAS RN: 137174-84-6	P0921 5g 25g  2,3',4',5',6'-Pentafluoro-acetophenone CAS RN: 652-29-9	P0849 5g 25g  Pentafluorobenzyl Alcohol CAS RN: 440-60-8	P0852 5g 25g  Pentafluorophenylhydrazine CAS RN: 828-73-9
P2244 1g 5g  2,3,4,5,6-Pentafluorobenzyl Cyanide CAS RN: 653-30-5	T3012 1g 5g  Trimethyl(pentafluorophenyl)ilsilane CAS RN: 1206-46-8	P0854 1mL 5mL  Flophemesyl Chloride CAS RN: 20082-71-7	T3352 1g 5g  Trimethoxy(pentafluorophenyl)silane CAS RN: 223668-64-2	P1343 5g  3-(Pentafluorophenyl)-pentafluoro-1-propene CAS RN: 67899-41-6

T3134	1g 5g	T3030	1g 5g	B3604	5g	C2700	5g	P2900	200mg 1g
Triethoxy(pentafluorophenyl)silane CAS RN: 20083-34-5		Trichloro[3-(pentafluorophenyl)propyl]silane CAS RN: 78900-02-4		Bis(pentafluorophenyl) Carbonate CAS RN: 59483-84-0		Chlorodimethyl[3-(2,3,4,5,6-pentafluorophenyl)propyl]silane CAS RN: 157499-19-9		Poly(pentafluorostyrene) CAS RN: 26838-55-1	

Ordering and Customer Service

TCI AMERICA

Tel : 800-423-8616 / 503-283-1681
Fax : 888-520-1075 / 503-283-1987
E-mail : Sales-US@TCIchemicals.com

TCI EUROPE N.V.

Tel : +32 (0)3 735 07 00
Fax : +32 (0)3 735 07 01
E-mail : Sales-EU@TCIchemicals.com

TCI Deutschland GmbH

Tel : +49 (0)6196 64053-00
Fax : +49 (0)6196 64053-01
E-mail : Sales-DE@TCIchemicals.com

Tokyo Chemical Industry UK Ltd.

Tel : +44 (0)1865 78 45 60
E-mail : Sales-UK@TCIchemicals.com

梯希爱(上海)化成工业发展有限公司

Tel : 800-988-0390 / 021-67121386
Fax : 021-6712-1385
E-mail : Sales-CN@TCIchemicals.com

Tokyo Chemical Industry (India) Pvt. Ltd.

Tel : 1800 425 7889 / 044-2262 0909
E-mail : Sales-IN@TCIchemicals.com

TOKYO CHEMICAL INDUSTRY CO., LTD.

Tel : +81 (0)3-5640-8878
E-mail : globalbusiness@TCIchemicals.com

Availability, price or specification of the listed products are subject to change without prior notice. Reproduction forbidden without the prior written consent of Tokyo Chemical Industry Co., Ltd.