

Bradford Solution (Ready-to-use solution) [for Protein determination]

Product Information

Product No.	:	B5702
Volume	:	500 mL

Description

This product is supplied as a ready-to-use solution for protein assay based on the method of Bradford. This product contains Coomassie Brilliant Blue G-250 (CBB G-250). When the dye containing CBB G-250 binds proteins, the absorption maximum of the dye shifts from 465 to 595 nm linearly with the quantity of the protein. This product requires protein standard solution (such as BSA).

 TCI AMERICA

 Tel:
 800-423-8616 • 503-283-1681

 Fax:
 888-520-1075 • 503-283-1987

 TCI EUROPE N.V.
 Tel:

 Tel:
 : 432 (0)3 735 07 00

 Fax:
 : 432 (0)3 735 07 01

様希愛(上海)化成工业发展有限公司 Tel:800-988-0390 * 021-67121386 Fax:021-6712-1385 **TCI Chemicals (India) Pvt. Ltd.** Tel:+91(0)44-2262 9909 Fax:+91(0)44-2262 8902



Direction for Use

Material

Standard protein solution: Protein solution whose concentration is known

Procedure

1. Bring the Bradford solution (Product No. B5702) to room temperature.

2. Gently mix the Bradford solution.

3. Prepare standard protein solutions (see Table 1).

4. Mix standard protein solutions or unknown protein samples with the Bradford solution (see Table 1).

5. Incubate for 5 minutes at room temperature.

6. Measure absorbance at 595 nm. The absorbance should be measured within 60 minutes of the start of the reaction.

7. Calculate the protein concentration.

Calculation

a. Calculate the corrected blank absorbance at 595 nm which is the difference between the average absorbance of the blank solution and that of standard solution.

b. Create a standard curve by plotting the corrected blank absorbance at 595 nm for each standard solution or protein sample against its concentration in μ g/mL.

c. Use the standard curve and determine the protein concentration of each unknown protein sample.

Assay	test tube	micro plate	micro assay	
Measurement	0.1 - 1.0	0.1 - 1.0	1.0 – 25	
range	mg/mL	mg/mL	µg/mL	
Protein standard or sample solution	20 µL	4 µL	500 μL	
Product No. P2575	1 mL	200 µL	500 µL	
Reaction	Incubate for 5 minutes at room temperature.			
Measurement	Within 1 hour, measure absorbance at 595 nm.			

Table 1 : Volume for test tube or micro plate assay

Table 2 : Compatible substance concentrations in protein sample

Substances at the following concentrations in the sample solutions do not affect the reaction results.

reaction results.			
Substance	concentration		
Buffers			
Glycine	100 mM		
Tris	2 M		
HCI	100 mM		
HEPES	100 mM		
MES	100 mM		
MOPS	100 mM		
PIPES	100 mM		
Glucose	1 M		
Sucrose	25 %		
Fructose	1 M		
Salts			
(NH ₄) ₂ SO ₄	1 M		
KCI	1 M		
MgCl ₂	50 mM		
CaCl ₂	10 mM		
NiCl ₂	10 mM		
ZnCl ₂	10 mM		
NaCl	2 M		
NaOH	100 mM		
NaH ₂ PO ₄	500 mM		
NaN ₃	0.50 %		
Chelating Ag	jents		
EDTA	100 mM		
EGTA	10 mM		
Sodium citrate	200 mM		
Solvents			
Acetone	10 %		
DMSO	10 %		
Ethanol	10 %		
Methanol	10 %		
Glycerol	10 %		
Detergents			
SDS	0.05 %		
Triton X-100	0.10 %		
Tween-20	0.10 %		
Denaturants			
DTT	100 mM		
Glutathione	1 mg / mL		
2-Mercaptoethanol	1 M		
Guanidine-HCI	e-HCI 1 M		
Urea 3 M			

東京仁並工業株式合社 Tel: 03-3668-0489 * 06-6228-1155 Fax: 03-3668-0520 * 06-6228-1158 TOKYO CHEMICAL INDUSTRY CO., LTD. Tel: +81 (0)3-5640-8878 Fax: ++81 (0)3-5640-8902 TCI AMERICA Tel: 800-423-861 • 503-283-1681 Fax: 888-520-1075 • 503-283-1987 TCI EUROPE N.V. Tel: +32 (0)3 735 07 00 Fax: +32 (0)3 735 07 01 梯希愛(上海)化成工业发展有限公司 Tel: 800-988-0390 ● 021-67121386 Fax: 021-6712-1385
TCI Chemicals (India) Pvt. Ltd. Tel: +91 (0)44-2262 0909 Fax: +91 (0)44-2262 8902