


Reusable Photosensitive Materials Showing Reversible Solid-Liquid Phase Transition



- R = C₆H₁₃ 4,4'-Bis(hexyloxy)-3-methylazobenzene
500mg / 5g [B4596]
- R = C₁₀H₂₁ 4,4'-Bis(decyloxy)-3-methylazobenzene
500mg / 5g [B4597]
- R = C₁₂H₂₅ 4,4'-Bis(dodecyloxy)-3-methylazobenzene
500mg / 5g [B4598]



Visible light (450nm)


(Figures are provided by Advanced Industrial Science and Technology (AIST), Japan.)

Advantage

Can be repeatedly liquefied from solid form by UV light and solidified from liquid form by visible light at room temperature


Application 1

Light Adhesion



Two glass slides containing this material adhere to each other by irradiation with 450nm light. After phase transition, irradiation with 365 nm light reliquifies the material and breaks adhesion.

Light Non-adhesion



These materials show a maximal adhesion value of 60 N cm⁻² (0.6 MPa). B4598 (R=C₁₂H₂₅) exhibits the highest adhesion value.

(provided from AIST)

Reference
Y. Norikane et al., *J. Photopolym. Sci. Technol.* **2016**, 29, 149.
<https://doi.org/10.2494/photopolymer.29.149>


*These materials are under invention licenses by the Advanced Industrial Science and Technology (AIST), Japan. If you intend to use these materials for industrial purposes, please contact us.

Reusable Photosensitive Materials Showing Reversible Solid-Liquid Phase Transition

Application 2

Photoresist for lithography

These materials can be easily blown and washed-off, making them applicable to simple lithography processes.




References

Y. Norikane *et al.*, *Org. Lett.* **2014**, *16*, 5012.
<https://doi.org/10.2494/photopolymer.29.149>

Y. Norikane, National Institute of Advanced Industrial Science and Technology, PCT Int. Appl. WO2013081155, **2013**.

(Figures are provided by AIST.)

Photoresponses of crystalline film in polarized light microscopy



*Each sample thickness is ca. 1 μm . Photoresponse time dependent on film thickness and light intensity.

(Figures are provided by AIST.)

Ordering and Customer Service

TCI AMERICA

Tel : 800-423-8616 / 503-283-1681
 Fax : 888-520-1075 / 503-283-1987
 E-mail : Sales-US@TCIchemicals.com

TCI EUROPE N.V.

Tel : +32 (0)3 735 07 00
 Fax : +32 (0)3 735 07 01
 E-mail : Sales-EU@TCIchemicals.com

TCI Deutschland GmbH

Tel : +49 (0)6196 64053-00
 Fax : +49 (0)6196 64053-01
 E-mail : Sales-DE@TCIchemicals.com

Tokyo Chemical Industry UK Ltd.

Tel : +44 (0)1865 78 45 60
 E-mail : Sales-UK@TCIchemicals.com

梯希爱(上海)化成工业发展有限公司

Tel : 800-988-0390 / 021-67121386
 Fax : 021-6712-1385
 E-mail : Sales-CN@TCIchemicals.com

Tokyo Chemical Industry (India) Pvt. Ltd.

Tel : 1800 425 7889 / 044-2262 0909
 E-mail : Sales-IN@TCIchemicals.com

TOKYO CHEMICAL INDUSTRY CO., LTD.

Tel : +81 (0)3-5640-8878
 E-mail : globalbusiness@TCIchemicals.com

* Chemicals itemized in this brochure are for research and testing use only. Please avoid use other than by chemically knowledgeable professionals. * Information such as listed products and its specifications and so on are subject to change without prior notice. * The contents may not be reproduced or duplicated in whole or in part without permission of Tokyo Chemical Industry Co., Ltd.