



Hole Transport Materials for Stable Perovskite Solar Cells

TOP-HTMs



TOP-HTM-α1
1g / 5g / 25g
[B5672]

TCI Original & Practical
Hole Transport Materials




TOP-HTM-α2
1g / 5g / 25g
[T3722]

Advantages

- Realizes a high PCE, with or without additives.¹⁾
- Produces a highly stable perovskite solar cell with low cost.
- It is possible to achieve even higher PCE and superior device stability by tuning the compositions of the perovskite layer.^{2,3)}

Device Fabrication Process



1. In a glove box filled with N₂ gas, hole transport layers are deposited on the perovskite layer by spin-coating (slope 5 s, 4000 rpm 30 s, slope 5 s).
2. The resulting film is dried on a hot plate at 70 °C for 30 minutes.
3. A metal electrode (Au, etc.) is thermally deposited on the hole transport layer.
4. The solar cell devices are stored in air with ~20% relative humidity to promote oxygen doping.


For more details, see the reference 1.

Application: Device performance of MAPbI₃-based solar cells ¹⁾

HTMs are dissolved in solvents at a concentration of 40 mg/mL.

- **With additives**
Chlorobenzene is used as a solvent. LiTFSI and TBP are added to the HTM solution.
- **Without additives**
1,1,2,2-Tetrachloroethane is used as a solvent.

Device Structure




HTM	PCE (%)	Stability (200 h)
TOP-HTM-α1 with additive	13.1	-
TOP-HTM-α1	15.0	0.50
TOP-HTM-α2 with additive	18.6	0.46
TOP-HTM-α2	16.6	0.69
Spiro-OMeTAD with additive	18.4	<0.10
Spiro-OMeTAD	12.1	-

Application: Improvement of perovskite layer ²⁾

TOP-HTM-α2 is dissolved in 1,1,2,2-tetrachloroethane without additives.

PCE and device stability of perovskite solar cells using TOP-HTM-α2 are gotten better by changing composition of perovskite layer from MAPbI₃ to FAPbI₃.

Device Structure



HTM / Perovskite	PCE (%)	Stability (400 h)
TOP-HTM-α2 / MAPbI ₃	14.31	0.75
TOP-HTM-α2 / FAPbI₃	16.86	0.90
Spiro-OMeTAD with additive / FAPbI ₃	15.11	0.69
Spiro-OMeTAD / FAPbI ₃	10.44	0.87

Application: Control of hole transport layer (HTL) morphology ³⁾


TOP-HTM- α 2 is dissolved in chlorobenzene or dichloromethane. LiTFSI and TBP as additives are added to the solution.

The composition of perovskite layer is



Perovskite solar cells using TOP-HT- α 2 with additives show superior PCE with improvement of HTL morphology by using dichloromethane instead of chlorobenzene.

Device Structure



HTM	Solvent	PCE (%)
TOP-HTM- α 2 with additive	Chlorobenzene	18.67
TOP-HTM-α2 with additive	Dichloromethane	20.18
Spiro-OMeTAD with additive	Chlorobenzene	18.85
Spiro-OMeTAD with additive	Dichloromethane	18.55

References

- 1) Additive-free, Cost-Effective Hole-Transporting Materials for Perovskite Solar Cells Based on Vinyl Triarylamines
H. Nishimura, I. Okada, T. Tanabe, T. Nakamura, R. Murdey, A. Wakamiya, *ACS Appl. Mater. Interfaces* **2020**, *12*, 32994.
<https://doi.org/10.1021/acsami.0c06055>
- 2) Experimental investigation of additive free-low-cost vinyl triarylamines based hole transport material for FAPbI₃-based perovskite solar cells to enhance efficiency and stability
A. Kumar, S. Singh, M. K A Mohammed, D. S Ahmed, *Mater. Res. Express* **2023**, *10*, 044003.
<https://doi.org/10.1088/2053-1591/accd41>
- 3) An Alternative to Chlorobenzene as a Hole Transport Materials Solvent for High-Performance Perovskite Solar Cells
S. H. Lee, S. B. Lim, J. Y. Kim, S. Lee, S. Y. O, M. Kim, *Crystals* **2023**, *13*, 1667.
<https://doi.org/10.3390/cryst13121667>

Hole Transport Materials for Stable Perovskite Solar Cells: TOP-HTMs

Related Products

Formamidinium Hydroiodide (= FAI) (99.99%, trace metals basis) [for Perovskite precursor]	1g / 5g / 25g [F1263]
Methylamine Hydroiodide (= MAI) (Low water content)	1g / 5g / 25g / 100g [M2556]
Spiro-OMeTAD	1g / 5g [T3672]
Lithium Bis(trifluoromethanesulfonyl)imide (= LiTFSI)	25g / 250g [B2542]
4-tert-Butylpyridine	5g / 25g [B0388]
Lead(II) Iodide (99.99%, trace metals basis) [for Perovskite precursor]	1g / 5g / 25g / 100g / 1kg [L0279]
Lead(II) Bromide [for Perovskite precursor]	1g / 5g / 25g [L0288]
Lead(II) Chloride (purified by sublimation) [for Perovskite precursor]	1g / 5g [L0291]
Lead(II) Chloride [for Perovskite precursor]	1g / 5g / 25g [L0292]

For further information please refer to our website at www.TCIchemicals.com.



TCI perovskite



Ordering and Customer Service

TCI AMERICA

Tel : 800-423-8616 / 503-283-1681
Fax : 888-520-1075 / 503-283-1987
E-mail : Sales-US@TCIchemicals.com

TCI EUROPE N.V.

Tel : +32 (0)3 735 07 00
Fax : +32 (0)3 735 07 01
E-mail : Sales-EU@TCIchemicals.com

TCI Deutschland GmbH

Tel : +49 (0)6196 64053-00
Fax : +49 (0)6196 64053-01
E-mail : Sales-DE@TCIchemicals.com

Tokyo Chemical Industry UK Ltd.

Tel : +44 (0)1865 78 45 60
E-mail : Sales-UK@TCIchemicals.com

梯希爱(上海)化成工业发展有限公司

Tel : 800-988-0390 / 021-67121386
Fax : 021-6712-1385
E-mail : Sales-CN@TCIchemicals.com

Tokyo Chemical Industry (India) Pvt. Ltd.

Tel : 1800 425 7889 / 044-2262 0909
E-mail : Sales-IN@TCIchemicals.com

TOKYO CHEMICAL INDUSTRY CO., LTD.

Tel : +81 (0)3-5640-8878
E-mail : globalbusiness@TCIchemicals.com

• Chemicals itemized in this brochure are for research and testing use only. Please avoid use other than by chemically knowledgeable professionals. • Information such as listed products and its specifications and so on are subject to change without prior notice. • The contents may not be reproduced or duplicated in whole or in part without permission of Tokyo Chemical Industry Co., Ltd.