


Controlled Radical Polymerization Chain Transfer Agents for TERP

Organotellurium-mediated radical polymerization (TERP), developed by Yamago *et al.*, is a type of controlled radical polymerizations (CRPs) also known as reversible-deactivation radical polymerizations. CRPs provide polymers with narrow molecular weight distributions, allows easy control of molecular weights, and enables the synthesis of block copolymers. In addition, TERP is advantageous for the development of functional polymers due to its wide monomer functional groups scope and the possibility of end group transformations. Furthermore, compared to reversible addition-fragmentation chain transfer (RAFT) polymerization, TERP offers a broader range of monomers that can be controlled with the same chain transfer agent.


Chain Transfer Agents (CTAs)



Ethyl 2-Methyl-2-(methyltellanyl)-propanoate

100mg / 1g


New [E1508]



2-Methyl-2-(methyltellanyl)-propanenitrile

100mg / 1g


New [M3520]



Methyl(1-phenylethyl)tellane

100mg / 1g

New [M3521]




2-Methyl-2-(methyltellanyl)-propanoic Acid

100mg

New [M3728]

Related Products


Reaction Additives



Dimethyl Ditelluride

100mg / 1g

New [D6090]



Diphenyl Ditelluride

1g / 5g

[D2718]


Controlled Radical Polymerization Chain Transfer Agents for TERP

Applications


Polymerization by TERP can be carried out by thermal condition (Scheme 1), thermal condition with radical initiators (Scheme 2), or photoirradiation condition. For the polymerization of methacrylates, the addition of ditellurides is effective to improve the control of the molecular weight distribution (Scheme 3).

M3728, which is a chain-transfer agent having a carboxy group, is suitable for emulsion polymerization in aqueous media.⁵⁾


Scheme 1. Synthesis of polystyrene³⁾




Scheme 2. Synthesis of poly(*N*-isopropylacrylamide)³⁾



Scheme 3. Synthesis of poly(methyl methacrylate)⁴⁾



Scheme 4. Synthesis of polystyrene via emulsion polymerization⁵⁾



- References**
- 1) S. Yamago, *Chem. Rev.* **2009**, *109*, 5051. <https://doi.org/10.1021/cr9001269>
 - 2) S. Yamago, *Bull. Chem. Soc. Jpn.* **2020**, *93*, 287. <https://doi.org/10.1246/bcsj.20190339>
 - 3) S. Yamago, *J. Polym. Sci. A Polym. Chem.* **2006**, *44*, 1. <https://doi.org/10.1002/pola.21154>
 - 4) S. Yamago et al., *Beilstein J. Org. Chem.* **2013**, *9*, 1607. <https://doi.org/10.3762/bjoc.9.183>
 - 5) S. Yamago et al., *ACS Macro Lett.* **2022**, *11*, 1331. <https://doi.org/10.1021/acsmacrolett.2c00594>

Ordering and Customer Service

TCI AMERICA

Tel : 800-423-8616 / 503-283-1681
 Fax : 888-520-1075 / 503-283-1987
 E-mail : Sales-US@TCIchemicals.com

TCI EUROPE N.V.

Tel : +32 (0)3 735 07 00
 Fax : +32 (0)3 735 07 01
 E-mail : Sales-EU@TCIchemicals.com

TCI Deutschland GmbH

Tel : +49 (0)6196 64053-00
 Fax : +49 (0)6196 64053-01
 E-mail : Sales-DE@TCIchemicals.com

Tokyo Chemical Industry UK Ltd.

Tel : +44 (0)1865 78 45 60
 E-mail : Sales-UK@TCIchemicals.com

梯希爱(上海)化成工业发展有限公司

Tel : 800-988-0390 / 021-67121386
 Fax : 021-6712-1385
 E-mail : Sales-CN@TCIchemicals.com

Tokyo Chemical Industry (India) Pvt. Ltd.

Tel : 1800 425 7889 / 044-2262 0909
 E-mail : Sales-IN@TCIchemicals.com

TOKYO CHEMICAL INDUSTRY CO., LTD.

Tel : +81 (0)3-5640-8878
 E-mail : globalbusiness@TCIchemicals.com

* Chemicals itemized in this brochure are for research and testing use only. Please avoid use other than by chemically knowledgeable professionals. * Information such as listed products and its specifications and so on are subject to change without prior notice. * The contents may not be reproduced or duplicated in whole or in part without permission of Tokyo Chemical Industry Co., Ltd.