


# Oligosaccharide Replacement of a Therapeutic Antibody by using Endo-M and Glycosynthase

As a result of chemoenzymatic glycoengineering, multi-form *N*-linked oligosaccharides of a therapeutic antibody were replaced by the structure fine-defined oligosaccharide

In recent years, the expectation of a therapeutic benefit for antibody drugs is growing and the development of industrial technology for antibody production is required. However, heterogeneity of glycosylation of antibody drugs has long been left unsolved. In this context, TCI achieved introduction of a MONO-form oligosaccharide to a therapeutic antibody by using our enzyme products, "Endo-M" and "Glycosynthase".




Please harness TCI's [Endo-M], [Glycosynthase] and [Chemically Synthesized Oligosaccharides] for the benefit of your Research & Development.

# Oligosaccharide Replacement of a Therapeutic Antibody by using Endo-M and Glycosynthase

Cleavage of multi-form oligosaccharides by Endo-M and coupling of MONO-form oligosaccharide by Glycosynthase were conducted under non-reducing condition. Verification of enzymatic reaction was performed with capillary electrophoresis (Fig. 1) and SDS-PAGE (Fig. 2A). The terminal sialic acid of the chemoenzymatically-transferred *N*-linked oligosaccharide to the therapeutic antibody was detected by lectin blotting using a sialic acid binding lectin (Fig. 2B).


**Fig. 1 Verification of enzymatic reaction**




After coupling of a MONO-form oligosaccharide to the Endo-M-treated antibody, a peak shift of the H chain of the therapeutic antibody was observed while a peak shift of L chain was not.

**Fig. 2 Detection of sialylated oligosaccharide**

**A) SDS\_PAGE (15% gel)**



**B) Lectin blotting**



**A) SDS-PAGE**

After coupling of a MONO-form oligosaccharide to the Endo-M-treated antibody, a band shift of the H chain of the therapeutic antibody was observed while a band shift of L chain was not.

**B) Lectin blotting**

In this study, SSA (*Sambucus sieboldiana* agglutinin) was used. A band of H chain of the therapeutic antibody coupled with MONO-form oligosaccharide was significantly stained.

## Related Products

<b>endo-<math>\beta</math>-N-Acetylglucosaminidase (Endo-M)</b>	100m units/vial [A1651]
from <i>Mucor hiemalis</i> expressed in <i>Candida boidinii</i>	
<b>Glycosynthase (Endo-M-N175Q)</b>	100m units/vial [G0365]
from <i>Mucor hiemalis</i> expressed in <i>Escherichia coli</i>	
<b>Disialyloctasaccharide</b>	20mg [D4065]
<b>Sialylglycopeptide (SGP)</b>	10mg [S0523]

References K. Yamamoto, S. Kadowaki, J. Watanabe, H. Kumagai, *Biochem. Biophys. Res. Commun.* **1994**, 203, 244.

M. Umekawa, T. Higashiyama, T. Tanaka, M. Noguchi, A. Kobayashi, S. Shoda, W. Huang, L-X. Wang, H. Ashida, K. Yamamoto, *Biochim. Biophys. Acta, Gen. Subj.* **2010**, 1800, 1203.

M. Umekawa, C. Li, T. Higashiyama, W. Huang, H. Ashida, K. Yamamoto, L-X. Wang, *J. Biol. Chem.* **2010**, 285, 511.

## Ordering and Customer Service

### TCI AMERICA

Tel : 800-423-8616 / 503-283-1681  
Fax : 888-520-1075 / 503-283-1987  
E-mail : Sales-US@TCIchemicals.com

### TCI EUROPE N.V.

Tel : +32 (0)3 735 07 00  
Fax : +32 (0)3 735 07 01  
E-mail : Sales-EU@TCIchemicals.com

### TCI Deutschland GmbH

Tel : +49 (0)6196 64053-00  
Fax : +49 (0)6196 64053-01  
E-mail : Sales-DE@TCIchemicals.com

### Tokyo Chemical Industry UK Ltd.

Tel : +44 (0)1865 784560  
Fax : +44 (0)1865 784561  
E-mail : Sales-UK@TCIchemicals.com

### TCI Chemicals (India) Pvt. Ltd.

Tel : 1800 425 7889 / 044-2262 0909  
Fax : 044-2262 8902  
E-mail : Sales-IN@TCIchemicals.com

### 梯希爱(上海)化成工业发展有限公司

Tel : 800-988-0390 / 021-67121386  
Fax : 021-6712-1385  
E-mail : Sales-CN@TCIchemicals.com

### TOKYO CHEMICAL INDUSTRY CO., LTD.

Tel : +81 (0)3-5640-8878  
Fax : +81 (0)3-5640-8902  
E-mail : globalbusiness@TCIchemicals.com

Availability, price or specification of the listed products are subject to change without prior notice. Reproduction forbidden without the prior written consent of Tokyo Chemical Industry Co., Ltd.