


# Sulfated Hyaluronic Acids

## —Cell Culture Applications—



**Hyaluronic Acid, High-Sulfated**  
**Hyaluronic Acid, Low-Sulfated**

10mg [**H1739**]

10mg [**H1740**]

### Maintaining the undifferentiated state of hiPS cells



**Differentiated  
(control)**



**Undifferentiated  
(Hyaluronic Acid, High-Sulfated)**

Sulfated hyaluronic acid maintains the undifferentiated state of hiPS cells under feeder-free and FGF2-free conditions.<sup>1)</sup>


*These data were kindly provided by Prof. Shoko Nishihara and Dr. Taichi Miura, Soka University.*

We offer Sulfated Hyaluronic Acids prepared by chemically modifying hyaluronic acid with sulfate groups. The degree of this sulfation is controllable by adjusting the reaction, allowing us to offer two separate products with differing sulfation degrees currently, both of which are animal-derived component free. Recent studies have revealed that our sulfated hyaluronic acids can protect growth factors from degradation and demonstrate functions similar to heparin. Given these properties, the development of novel applications of sulfated hyaluronic acids within the wide variety of fields such as cell culture, tissue engineering and regenerative medicine are to be expected in the coming years.


Glycosaminoglycans (GAGs) are important polysaccharide components of the extracellular matrix found both within tissues as well as on the surface of cells. Hyaluronic acid, one of the GAGs, is found in the body in an unbound form, and comprises repeating disaccharide units of glucuronic acid (GlcA) and *N*-acetylglucosamine (GlcNAc). GAGs such as heparin have sulfate groups and physically interact with growth factors and other various functional molecules. Heparin binds to FGF2 to both protect it from degradation and stimulate receptor binding, resulting in physical responses including cell proliferation.<sup>2)</sup> Hyaluronic acid itself does not contain sulfate groups, however our artificially sulfated hyaluronic acids have been found to possess various functions.

## Binding of Sulfated Hyaluronic Acids to FGF2

**Binding of Various GAGs to FGF2 (Competitive ELISA)**



**Evaluation of sulfated hyaluronic acid binding to FGF2**




The FGF2 binding ability of heparin and sulfated hyaluronic acids was evaluated competitively and normalized to that of heparin.


Low-sulfated hyaluronic acid showed weaker, and highly sulfated hyaluronic acid showed stronger binding to FGF2 than heparin.

## Protection of FGF from Degradation by Sulfated Hyaluronic Acids

**FGF1**



**FGF2**



**An Example of Protective Effects on Each Protein**

	FGF1	FGF2	KGF	BMP2	IL-6
HA	—	—	—	—	—
HA-LS	++	++	+	+	—
HA-HS	++	+	+	—	+

Left: The ability of each sulfated hyaluronic acid to stabilize FGF1 and FGF2 was examined by mixing each pair together and incubating.

This data shows that sulfated hyaluronic acid protects unstable FGF from degradation.

HA: Hyaluronic Acid

HA-LS: Hyaluronic Acid, Low-Sulfated

HA-HS: Hyaluronic Acid, High-Sulfated


## Examples of Applications of Sulfated Hyaluronic Acids

Prof. Shoko Nishihara, Dr. Taichi Miura, *et al.* of Soka University successfully cultured human iPS cells in an undifferentiated state under feeder-free and FGF2-free conditions through the addition of highly sulfated hyaluronic acid to their cell culture medium (see photographs on first page). Their study showed that highly sulfated hyaluronic acid promotes FGF2 signaling in human iPS cells.<sup>1)</sup>

In addition, a more recent study revealed that administration of FGF1 in combination with highly sulfated hyaluronic acid ameliorates radiation-induced intestinal damage, and highly sulfated hyaluronic acid, unlike heparin, only possesses a weak anti-coagulant effect. These results suggest that highly sulfated hyaluronic acid may be a promising radioprotective agent with a low risk of hemorrhage which will be a problem in the combination of FGF1 and heparin (reported by Dr. Fumiaki Nakayama, Group Leader, and Dr. Taichi Miura of the Regenerative Therapy Research Group, Department of Radiation Regulatory Science Research, National Institute of Radiological Sciences (NIRS), and the Institute of Glycotechnology of Tokyo Chemical Industry Co., Ltd. (TCI)).<sup>3)</sup>

Sulfated hyaluronic acids have also been found to inhibit hyaluronidase, an enzyme that degrades hyaluronic acid.


Furthermore, it has been reported that sulfated hyaluronic acids inhibit the activity of heparanase, an enzyme that cleaves heparan sulfate and macromolecular heparin and which is essential for cancer cell invasion. It was also shown that highly sulfated hyaluronic acid inhibited heparanase-mediated extension protrusion of cancer cells in a 3D-culture system.<sup>4)</sup>



**Hypothesized Model of Action of Sulfated Hyaluronic Acid**

## Future Outlook

Given these various functions, our sulfated hyaluronic acids show potential to be used as alternatives for sulfated polysaccharides in various fields such as large-scale cell culture, regenerative medicine, wound healing, and cancer therapy.




**As an animal-derived component free biomaterial which enhances FGF activity to regulate cell differentiation and proliferation**

- As an additive to cell culture medium
- For industry applications that require sustainable large-scale cell culture
- As a heparin alternative
- Applicable to developmental biology, bioengineering, animal cell engineering
- To aid physiological processes that require tissue repair such as wound healing
- In regenerative medicine, including regeneration of skin and teeth

# Sulfated Hyaluronic Acids —Cell Culture Applications—


## Related Products

### Hyaluronic Acids



<b>R=H</b>	<b>Hyaluronic Acid from Bacteria</b>	1g	<b>[H1807]</b>
<b>R=Na</b>	<b>Hyaluronic Acid Sodium Salt from Bacteria</b>	100mg / 1g	<b>[H1791]</b>
<b>R=K</b>	<b>Hyaluronic Acid Potassium Salt from Bacteria</b>	1g	<b>[H1808]</b>
<b>R=H</b>	<b>Hyaluronic Acid from Cockscomb</b>	1g	<b>[H0595]</b>
<b>R=Na</b>	<b>Hyaluronic Acid Sodium Salt from Cockscomb</b>	100mg / 1g	<b>[H0603]</b>
<b>R=K</b>	<b>Hyaluronic Acid Potassium Salt from Cockscomb</b>	1g	<b>[H0652]</b>

### Oligo Hyaluronic Acids: $[GlcA\beta(1-3)GlcNAc\beta(1-4)]_n-GlcA\beta(1-3)GlcNAc$



<b>n=1</b>	<b>Hyaluronate Tetrasaccharide</b>	1mg / 5mg	<b>[H1284]</b>
<b>n=2</b>	<b>Hyaluronate Hexasaccharide</b>	1mg / 5mg	<b>[H1285]</b>
<b>n=3</b>	<b>Hyaluronate Octasaccharide</b>	1mg	<b>[H1148]</b>
<b>n=4</b>	<b>Hyaluronate Decasaccharide</b>	1mg	<b>[H1149]</b>
<b>n=5</b>	<b>Hyaluronate Dodecasaccharide</b>	1mg	<b>[H1666]</b>

### Other Glycosaminoglycans

<b>Heparin Sodium Salt from Hog intestine</b>	100mg / 1g	<b>[H0393]</b>
<b>Chondroitin Sulfate Sodium Salt</b>	25g / 100g	<b>[C0335]</b>
<b>Dermatan Sulfate Sodium Salt</b>	20mg / 100mg	<b>[D3672]</b>

#### Abbreviations and Synonyms

**FGF:** Fibroblast Growth Factor  
**FGF1:** acidic FGF (aFGF)  
**FGF2:** basic FGF (bFGF)

**KGF:** Keratinocyte Growth Factor (FGF7)  
**BMP2:** Bone Morphogenetic Protein 2  
**IL-6:** Interleukin-6

**CS:** Chondroitin Sulfate  
**KS:** Keratan Sulfate  
**HS:** Heparan Sulfate

- References**
- 1) Highly sulfated hyaluronic acid maintains human induced pluripotent stem cells under feeder-free and bFGF-free conditions  
T. Miura *et al.*, *Biochem. Biophys. Res. Commun.* **2019**, 518, 506. <https://doi.org/10.1016/j.bbrc.2019.08.082>
  - 2) Chapter 38. Proteins That Bind Sulfated Glycosaminoglycans  
D. Xu *et al.*, in *Essentials of Glycobiology*, 4th ed., ed. by A. Varki *et al.*, Cold Spring Harbor Laboratory Press, New York, 2022.  
<https://www.ncbi.nlm.nih.gov/books/NBK579914/>
  - 3) High-Sulfated Hyaluronic Acid Ameliorates Radiation-Induced Intestinal Damage Without Blood Anticoagulation  
T. Miura *et al.*, *Adv. Radiat. Oncol.* **2022**, 7, 100900. <https://doi.org/10.1016/j.adro.2022.100900>
  - 4) Sulfated Hyaluronan Binds to Heparanase and Blocks Its Enzymatic and Cellular Actions in Carcinoma Cells  
J. Shi *et al.*, *Int. J. Mol. Sci.* **2022**, 23, 5055. <https://doi.org/10.3390/ijms23095055>

#### Ordering and Customer Service

##### TCI AMERICA

Tel : 800-423-8616 / 503-283-1681  
Fax : 888-520-1075 / 503-283-1987  
E-mail : Sales-US@TCIchemicals.com

##### TCI EUROPE N.V.

Tel : +32 (0)3 735 07 00  
Fax : +32 (0)3 735 07 01  
E-mail : Sales-EU@TCIchemicals.com

##### TCI Deutschland GmbH

Tel : +49 (0)6196 64053-00  
Fax : +49 (0)6196 64053-01  
E-mail : Sales-DE@TCIchemicals.com

##### Tokyo Chemical Industry UK Ltd.

Tel : +44 (0)1865 78 45 60  
E-mail : Sales-UK@TCIchemicals.com

##### 梯希爱(上海)化成工业发展有限公司

Tel : 800-988-0390 / 021-67121386  
Fax : 021-6712-1385  
E-mail : Sales-CN@TCIchemicals.com

##### Tokyo Chemical Industry (India) Pvt. Ltd.

Tel : 1800 425 7889 / 044-2262 0909  
E-mail : Sales-IN@TCIchemicals.com

##### TOKYO CHEMICAL INDUSTRY CO., LTD.

Tel : +81 (0)3-5640-8878  
E-mail : globalbusiness@TCIchemicals.com

• Chemicals itemized in this brochure are for research and testing use only. Please avoid use other than by chemically knowledgeable professionals. • Information such as listed products and its specifications and so on are subject to change without prior notice. • The contents may not be reproduced or duplicated in whole or in part without permission of Tokyo Chemical Industry Co., Ltd.