

Mitochondrial Isolation Kit

Mitochondrial Isolation Kit

1 kit [M3527]

M3527 can be used to easily isolate mitochondria from cultured mammalian cells and tissue.

Isolated mitochondria can be used for downstream analyses such as Western blotting.

Kit Components

- **Mitochondrial Isolation Reagent A**
for approx. 30 samples
- **Mitochondrial Isolation Reagent B**
for approx. 30 samples
- **Mitochondrial Isolation Reagent C**
for approx. 30 samples



Application: Isolation from Mouse-Derived Cells

1. Collect cells by centrifugation and remove the supernatant without disturbing the cell pellet. Use a cell scraper for detaching adherent cells.
2. Resuspend cells in Mitochondrial Isolation Reagent A and vortex for 5 seconds. Add 400 μ L per 1×10^7 cells.
3. Incubate on ice for 2 minutes.
4. Add Mitochondrial Isolation Reagent B to the cell suspension and vortex for 5 seconds. Use 5 μ L per 400 μ L of suspension from step 2.
5. Incubate on ice for 5 minutes, vortexing once per minute.
6. Add Mitochondrial Isolation Reagent C and invert several times (do not vortex). Use a volume equal to the volume of Reagent A used in step 2.
7. Centrifuge at 700 x g for 10 minutes at 4°C.
8. Transfer the supernatant to a new tube and centrifuge at 12,000 x g for 15 minutes at 4°C.
9. Collect the supernatant (cytoplasmic fraction) and add 300 μ L of Mitochondrial Isolation Reagent C to the pellet (mitochondria).
10. Centrifuge at 12,000 x g for 5 minutes at 4°C and remove the supernatant.
11. Use mitochondria (pellet) for downstream experiments.




Figure 1.

Western blotting of **M3527**-isolated cytoplasmic fraction (C) and mitochondria (M), detected with VDAC1/2 antibody.

Mitochondrial proteins were extracted with RIPA Buffer.

M3527 allows for highly efficient purification of mitochondria

Mitochondrial Isolation Kit

Application: Isolation from Mouse Liver

1. Weigh 50-100 mg of mouse liver and wash with PBS(-).
2. Add 400 μ L of Mitochondrial Isolation Reagent A and homogenize.
3. Incubate on ice for 2 minutes.
4. Add 5 μ L of Mitochondrial Isolation Reagent B and vortex for 5 seconds.
5. Incubate on ice for 5 minutes, vortexing once per minute.
6. Add 400 μ L of Mitochondrial Isolation Reagent C and invert several times (do not vortex).
7. Centrifuge at 700 x g for 10 minutes at 4°C.
8. Transfer the supernatant to a new tube and centrifuge at 12,000 x g for 15 minutes at 4°C.
9. Collect the supernatant (cytoplasmic fraction) and add 300 μ L of Mitochondrial Isolation Reagent C to the pellet (mitochondria).
10. Centrifuge at 12,000 x g for 5 minutes at 4°C and remove supernatant.
11. Use mitochondria (pellet) for downstream experiments.




Figure 2.

Western blotting of **M3527**-isolated cytoplasmic fraction (C) and mitochondria (M), detected using VDAC1/2 antibody.

Mitochondrial proteins were extracted using RIPA Buffer.

M3527 allows for highly efficient purification of mitochondria from mouse liver.

Ordering and Customer Service

TCI AMERICA

Tel : 800-423-8616 / 503-283-1681
Fax : 888-520-1075 / 503-283-1987
E-mail : Sales-US@TCIchemicals.com

TCI EUROPE N.V.

Tel : +32 (0)3 735 07 00
Fax : +32 (0)3 735 07 01
E-mail : Sales-EU@TCIchemicals.com

TCI Deutschland GmbH

Tel : +49 (0)6196 64053-00
Fax : +49 (0)6196 64053-01
E-mail : Sales-DE@TCIchemicals.com

Tokyo Chemical Industry UK Ltd.

Tel : +44 (0)1865 78 45 60
E-mail : Sales-UK@TCIchemicals.com

梯希爱(上海)化成工业发展有限公司

Tel : 800-988-0390 / 021-67121386
Fax : 021-6712-1385
E-mail : Sales-CN@TCIchemicals.com

Tokyo Chemical Industry (India) Pvt. Ltd.

Tel : 1800 425 7889 / 044-2262 0909
E-mail : Sales-IN@TCIchemicals.com

TOKYO CHEMICAL INDUSTRY CO., LTD.

Tel : +81 (0)3-5640-8878
E-mail : globalbusiness@TCIchemicals.com

Availability, price or specification of the listed products are subject to change without prior notice. Reproduction forbidden without the prior written consent of Tokyo Chemical Industry Co., Ltd.