

Xanthine / Hypoxanthine Assay Kit

Xanthine / Hypoxanthine Assay Kit

1kit [X0094]

The two purine bases xanthine and hypoxanthine can be found in most tissues and biological fluids. Xanthine is produced during the digestion of purine to uric acid by the action of xanthine oxidase. This kit measures the xanthine / hypoxanthine concentration of biological samples upon comparison with a xanthine standard solution.

Advantages

- Enables sensitive measurement of Xanthine / Hypoxanthine maximum to 400 μM
- Quick procedure with results within 30 minutes
- Suitable for various biological samples

Kit Components

- 10 x Assay Buffer 1vial
- Advanced Enzyme II Solution 1vial
- Advanced Substrate Solution 1vial
- Enzyme Solution 1vial
- 20mM Xanthine Standard Solution 1vial

Each vial is sufficient for 96 tests.

Example: Xanthine / Hypoxanthine quantification by colorimetric assay

‡ Xanthine / Hypoxanthine quantification from 20 μM to 400 μM can be tested via the colorimetric assay.

1. Prepare 1 x Assay Buffer

Mix 20 mL of 10 x XO Assay Buffer with 180 mL of distilled water.

2. Prepare standard curve

2.1. Mix 10 μL of Xanthine Standard Solution with 40 μL of 1 x Assay Buffer to make a 4 mM Xanthine.

2.2. Prepare the xanthine standard solutions shown in the right table.

3. Prepare standards and samples

Aliquot 50 μL of standards and samples to separate wells of a 96 well plate.

4. Prepare reaction mixture

Mix an appropriate amount of each reagent for the number of wells to be assayed. The volumes below represent the amount of each reagent needed per reaction per well.

- Enzyme Solution : 1 μL
- Advanced Substrate Solution : 1 μL
- Advanced Enzyme II Solution : 1 μL
- 1 x Assay Buffer : 47 μL

5. Start reaction

Add 50 μL of reaction mixture (prepared in step 4) to each well and mix well. Incubate the reaction at room temperature in the dark for 30 minutes.

6. Measure Optical Density (OD)

Measure OD at a wavelength between 550 - 580 nm.

7. Calculate activity

7.1. Plot standards to obtain a standard curve and calculate the slope.

7.2. Use the slope, along with the OD value of the sample to calculate xanthine / hypoxanthine concentration based on the following formula:

$$\text{Xanthine / Hypoxanthine Concentration } (\mu\text{M}) = \frac{\text{OD}_{\text{Sample}} - \text{OD}_{\text{Blank}}}{\text{Slope}} \times \text{Dilution ratio}$$


[Experimental Data]

Slope = 0.002194, $\text{OD}_{\text{Mouse Spleen Cell Supernatant}} - \text{OD}_{\text{Blank}} = 0.0788$, Dilution ratio = 1

Based on the above data, the Xanthine / Hypoxanthine contained in Mouse Spleen Cell Supernatant is 35.92 μM .

(Measurement condition : room temperature, pH 7.5)

No.	4 mM Xanthine	1 x Assay Buffer	Final Conc. (μM)
1	0 μL	200 μL	0
2	6 μL	194 μL	120
3	12 μL	188 μL	240
4	20 μL	180 μL	400



Xanthine / Hypoxanthine Assay Kit

Example: Xanthine / Hypoxanthine quantification by fluorometric assay

‡ Xanthine / Hypoxanthine quantification up to 20 μM can be tested via the fluorescence assay.

1. Prepare 1 x Assay Buffer

Mix 20 mL of 10 x XO Assay Buffer with 180 mL of distilled water.

2. Plot standard curve

2.1. Mix 2 μL of 20 mM Xanthine Standard Solution with 198 μL of 1 x XO Assay Buffer to make a 200 μM xanthine solution.

2.2. Prepare the xanthine standard solutions shown in the right table.

3. Prepare standards and samples

Aliquot 50 μL of standards and samples to separate wells of a 96 well plate.

4. Prepare reaction mixture

Mix an appropriate amount of each reagent for the number of wells to be assayed. The volumes below represent the amount of each reagent needed per reaction per well.

- Enzyme Solution : 1 μL
- Advanced Substrate Solution : 1 μL
- Advanced Enzyme II Solution : 1 μL
- 1 x Assay Buffer : 47 μL

5. Start reaction

Add 50 μL of reaction mixture (prepared in step 4) to each well and mix well. Incubate the reaction at room temperature in the dark for 30 minutes.

6. Measure Fluorescent intensity (F)

Measure F under the condition of $\lambda_{\text{ex}} = 520 - 550 \text{ nm}$ and $\lambda_{\text{em}} = 585 - 595 \text{ nm}$.

7. Calculate activity

7.1. Plot standards to obtain a standard curve and calculate the slope.

7.2. Use the slope, along with the fluorescence intensity value of the sample to calculate xanthine / hypoxanthine concentration based on the following formula:

$$\text{Xanthine / Hypoxanthine Concentration } (\mu\text{M}) = \frac{F_{\text{Sample}} - F_{\text{Blank}}}{\text{Slope}} \times \text{Dilution ratio}$$

[Experimental Data]


Slope = 773.7638, $F_{\text{Mouse Serum}} - F_{\text{Blank}} = 9091$, Dilution ratio = 1

Based on the above data, the Xanthine / Hypoxanthine contained in mouse serum is 11.75 μM .

(Measurement condition : room temperature, pH 7.5)

No.	200 μM Xanthine	1 x Assay Buffer	Final Conc. (μM)
1	0 μL	200 μL	0
2	6 μL	194 μL	6
3	12 μL	188 μL	12
4	20 μL	180 μL	20

Fluorometric Assay Standard Curve



Ordering and Customer Service

TCI AMERICA

Tel : 800-423-8616 / 503-283-1681
Fax : 888-520-1075 / 503-283-1987
E-mail : Sales-US@TCIchemicals.com

TCI EUROPE N.V.

Tel : +32 (0)3 735 07 00
Fax : +32 (0)3 735 07 01
E-mail : Sales-EU@TCIchemicals.com

TCI Deutschland GmbH

Tel : +49 (0)6196 64053-00
Fax : +49 (0)6196 64053-01
E-mail : Sales-DE@TCIchemicals.com

Tokyo Chemical Industry UK Ltd.

Tel : +44 (0)1865 78 45 60
E-mail : Sales-UK@TCIchemicals.com

梯希爱(上海)化成工业发展有限公司

Tel : 800-988-0390 / 021-67121386
Fax : 021-6712-1385
E-mail : Sales-CN@TCIchemicals.com

Tokyo Chemical Industry (India) Pvt. Ltd.

Tel : 1800 425 7889 / 044-2262 0909
E-mail : Sales-IN@TCIchemicals.com

TOKYO CHEMICAL INDUSTRY CO., LTD.

Tel : +81 (0)3-5640-8878
E-mail : globalbusiness@TCIchemicals.com

• Chemicals itemized in this brochure are for research and testing use only. Please avoid use other than by chemically knowledgeable professionals. • Information such as listed products and its specifications and so on are subject to change without prior notice. • The contents may not be reproduced or duplicated in whole or in part without permission of Tokyo Chemical Industry Co., Ltd.