



Chiral Amine Catalyzing Enantioselective Fluorination of α -Branched Aldehydes



Advantages

- Catalyze the asymmetric fluorination of α -branched aldehydes using *N*-fluorobenzenesulfonimide (NFSI) as a fluorine source.
- The generated chiral tertiary fluorides can be converted into various optically-active compounds.

Application



Ethyl (11*bR*)-4-Amino-2,6-bis(3,5-di-*tert*-butylphenyl)-4,5-dihydro-3*H*-cyclohepta-[1,2-*a*:7,6-*a'*]dinaphthalene-4-carboxylate 50mg **[E1267]**

This product was produced by collaboration with Assoc. Prof. Kazutaka Shibatomi, Toyohashi University of Technology.

Related Products

3,5-Dinitrobenzoic Acid

25g / 500g **[D0824]**

***N*-Fluorobenzenesulfonimide (= NFSI)**

5g / 25g **[F0335]**

Sodium Borohydride

25g / 100g / 500g **[S0480]**

Triethyl Phosphonoacetate

25g / 100g / 500g **[D1523]**



Chiral Amine Catalyzing Enantioselective Fluorination of α -Branched Aldehydes

Introduction of the Researcher

Shibatomi Laboratory

Department of Environmental and Life Sciences, Toyohashi University of Technology



The members of the Shibatomi laboratory with Assoc. Prof. Kazutaka Shibatomi on the far left

Research Description

The Shibatomi group aims to develop the new synthetic methods of organic molecules, especially focusing on design and synthesis of novel chiral catalysts and their application to the asymmetric reactions. The Shibatomi group is also developing the efficient synthetic method for chiral pharmaceutical and agricultural compounds with the above-mentioned chiral catalysis. Recently, the Shibatomi group found highly enantioselective halogenation of carbonyl compounds and applied this method for the synthesis of a GPR119 agonist which is a potential drug for type 2 diabetes.

Experimental Procedure (synthesis of 2-fluoroalcohol; Ar=Ph, R¹=CH₃)

To a solution of **E1267** (20 mg, 0.026 mmol, 10 mol%) in toluene (0.54 mL) is added 3,5-dinitrobenzoic acid (5.5 mg, 0.026 mmol, 10 mol%), 2-phenylpropionaldehyde (52 mg, 0.39 mmol, 1.5 eq.), and *N*-fluorobenzenesulfonimide (NFSI) (0.26 mmol, 82 mg, 1 eq.) at 0 °C. The reaction mixture is stirred at 0 °C for 48 h, then poured into CH₃OH/CH₂Cl₂ (1:4, 1.3 mL) at 0 °C. To this solution, NaBH₄ (1.6 mmol, 6 eq.) is added, and the mixture is stirred at room temperature for 1 h. The reaction is quenched with saturated aq. NH₄Cl, and the mixture is extracted with Et₂O. The organic layer is dried over Na₂SO₄, and concentrated. The residue is purified by silica gel chromatography (eluent: hexane/ethyl acetate = 3/1) to give (*S*)-2-fluoro-2-phenylpropan-1-ol (34.5 mg, 86% yield based on NFSI, 95% ee) as a white solid.

Ordering and Customer Service

TCI AMERICA

Tel : 800-423-8616 / 503-283-1681
Fax : 888-520-1075 / 503-283-1987
E-mail : Sales-US@TCIchemicals.com

TCI EUROPE N.V.

Tel : +32 (0)3 735 07 00
Fax : +32 (0)3 735 07 01
E-mail : Sales-EU@TCIchemicals.com

TCI Deutschland GmbH

Tel : +49 (0)6196 64053-00
Fax : +49 (0)6196 64053-01
E-mail : Sales-DE@TCIchemicals.com

Tokyo Chemical Industry UK Ltd.

Tel : +44 (0)1865 784560
Fax : +44 (0)1865 784561
E-mail : Sales-UK@TCIchemicals.com

TCI Chemicals (India) Pvt. Ltd.

Tel : 1800 425 7889 / 044-2262 0909
Fax : 044-2262 8902
E-mail : Sales-IN@TCIchemicals.com

梯希爱(上海)化成工业发展有限公司

Tel : 800-988-0390 / 021-67121386
Fax : 021-6712-1385
E-mail : Sales-CN@TCIchemicals.com

TOKYO CHEMICAL INDUSTRY CO., LTD.

Tel : +81 (0)3-5640-8878
Fax : +81 (0)3-5640-8902
E-mail : globalbusiness@TCIchemicals.com

Availability, price or specification of the listed products are subject to change without prior notice. Reproduction forbidden without the prior written consent of Tokyo Chemical Industry Co., Ltd.