




# Chemiluminescence Reagents



Chemiluminescence Compounds

Chemiluminescence Enhancers

# Chemiluminescence Reagents

Chemiluminescence is the phenomenon of energy released as light when the excited state of molecules by oxidation returns to the ground state. Organic substances which exhibit chemiluminescence are *Cypridina* luciferin, firefly luciferin, oxalate, luminol, lucigenin, etc.

*Cypridina* luciferin analogs achieve chemiluminescence through reaction with superoxide ( $O_2^-$ ) or singlet oxygen ( $^1O_2$ ). Using this characteristic, CLA [A5307] and MCLA [A5309] have been used in the research of the functions of leukocytes. The maximum emission wavelengths of CLA and MCLA are 380nm and 465nm, respectively. FCLA [A5310] developed by Goto *et al.* is characterized by emitting light at a longer wavelength (532nm). Furthermore, Red-CLA [A5311] developed by Teranishi emits at an even longer wavelength. Red-CLA exhibits high emission intensity by reaction with superoxide, and can be used for efficient analysis of superoxide at the longest wavelength (610nm).

Firefly luciferin reacts with ATP in the presence of luciferase and magnesium ion to provide oxyluciferin via luciferyl-adenylic acid. The light with a wavelength of 562nm is emitted when activated molecules return to the ground state by decomposition of oxyluciferin. Using this characteristic, firefly luciferin is used in the trace detection of ATP and the activity measurement of nucleotide phosphatase and others.

Luminol reacts with hydrogen peroxide in the presence of metals such as iron or its complexes to emit strong blue light with a wavelength of 460nm. This reaction is called the luminol test and is applied to the identification of blood stains in forensic science. Luminol is also used for trace detection of hydrogen peroxide and metals which catalyze this reaction.

On the other hand, oxalates are oxidized by hydrogen peroxide etc. to produce 1,2-dioxetanediones. When these substances are decomposed, they transfer energy to coexisting fluorescent substances which are elevated to the excited state. These excited fluorescent substances emit light during relaxation to the ground state. Thus, the emission wavelength can be changed by choosing coexisting fluorescent substances. This technique is applied to HPLC detection systems.


## References

- 1) CLA
  - a) T. Kobayashi, K. Saga, S. Shimizu, T. Goto, *Agric. Biol. Chem.* **1981**, 45, 1403.
  - b) M. Nakano, K. Sugioka, Y. Ushijima, T. Goto, *Anal. Biochem.* **1986**, 159, 363.
- 2) MCLA
  - a) A. Nishida, H. Kimura, M. Nakano, T. Goto, *Clin. Chim. Acta* **1989**, 179, 177.
  - b) A. Takahashi, H. Totsune-Nakano, M. Nakano, S. Mashiko, N. Suzuki, C. Ohma, H. Inaba, *FEBS Lett.* **1989**, 246, 117.
  - c) A. Takahashi, M. Nakano, S. Mashiko, H. Inaba, *FEBS Lett.* **1990**, 261, 369.
  - d) T. Sakurai, K. Sugioka, M. Nakano, *Biochim. Biophys. Acta - Lipids and Lipid Metabolism* **1990**, 1043, 27.
  - e) S. Koga, M. Nakano, K. Uehara, *Arch. Biochem. Biophys.* **1991**, 289, 223.
- 3) CLA, MCLA
  - a) M. Nakano, *Tanpakushitu Kakusan Koso* **1988**, 33, 2699.
  - b) H. Sawada, M. Nakayama, *Kagaku to Kogyo (Tokyo)* **1988**, 41, 1163.
  - c) H. Sawada, M. Nakayama, *Yukagaku* **1989**, 38, 103.
  - d) M. Nakano, H. Kimura, *Rinsho Kensa* **1989**, 33, 256.
  - e) K. Imai, in *Bioluminescence and Chemiluminescence*, Hirokawa Publishing, **1989**.
  - f) H. Sawada, K. Masuyama, M. Nakayama, *Yukagaku* **1990**, 39, 47.
  - g) M. Nakano, *Methods Enzymol.* **1990**, 186, 227.
  - h) M. Nakano, *Methods Enzymol.* **1990**, 186, 585.
- 4) CLA, MCLA, FCLA
  - a) N. Suzuki, K. Suetsuna, S. Mashiko, T. Nomoto, Y. Toya, B. Yoda, H. Inaba, *Chem. Express* **1990**, 5, 537.
  - b) N. Suzuki, K. Suetsuna, S. Mashiko, B. Yoda, T. Nomoto, Y. Toya, H. Inaba, T. Goto, *Agric. Biol. Chem.* **1991**, 55, 157.
- 5) Red-CLA
  - a) K. Teranishi, Tokyo Kasei Kogyo Co., Ltd., JP Patent 4453528.
  - b) K. Teranishi, *Luminescence* **2007**, 22, 147.
- 6) Firefly luciferase
  - a) A. Lundin, *Methods Enzymol.* **2000**, 305, 346.
  - b) J. J. Lemasters, C. R. Hackenbrock, *Methods Enzymol.* **1979**, 56, 530.


## Chemiluminescence Compounds

### Cypridina Luciferin Analogs


A5307 10mg

CLA  
CAS RN: 19953-58-3


A5309 10mg

MCLA  
CAS RN: 128322-44-1

A5310 10mg


FCLA Free Acid  
CAS RN: 133118-06-6

A5311 1mg

Red-CLA  
CAS RN: 886840-56-8


### Firefly Luciferin

A5030 10mg 50mg


D-(-)-Luciferin  
CAS RN: 2591-17-5

### Luminols


A5301 1g 25g

Luminol  
CAS RN: 521-31-3

A5300 1g 10g


Isoluminol  
CAS RN: 3682-14-2

A5304 100mg 1g


ABEI  
CAS RN: 66612-29-1

### Oxalates


A5302 1g 5g

TCPO  
CAS RN: 1165-91-9


O0236 5g 25g

Bis[3,4,6-trichloro-2-(pentylloxycarbonyl)-phenyl] Oxalate  
CAS RN: 30431-54-0

A5303 1g 5g


DNPO  
CAS RN: 16536-30-4

A5305 1g


PFPO  
CAS RN: 16536-48-4

### Acridinium Salts

B1203 1g 5g


Lucigenin  
CAS RN: 2315-97-1

B4339 100mg 1g


MMT  
CAS RN: 469865-01-8

## Chemiluminescence Enhancers


S0887 25g 500g

Sodium Perborate  
Tetrahydrate  
CAS RN: 10486-00-7


B3910 1g 5g

4-Bromophenol  
CAS RN: 106-41-2


I0840 1g 5g

4-Iodophenol  
CAS RN: 540-38-5

M2325 100mg 1g

4-(2-Methyl-4-thiazolyl)-phenol  
CAS RN: 30686-73-8

D4178 100mg

4-(4,5-Diphenyl-1H-imidazol-2-yl)phenol  
CAS RN: 1752-94-9

---

**Ordering and  
Customer Service**

**TCI AMERICA**

Tel : 800-423-8616 / 503-283-1681  
Fax : 888-520-1075 / 503-283-1987  
E-mail : Sales-US@TCIchemicals.com

**TCI EUROPE N.V.**

Tel : +32 (0)3 735 07 00  
Fax : +32 (0)3 735 07 01  
E-mail : Sales-EU@TCIchemicals.com

**TCI Deutschland GmbH**

Tel : +49 (0)6196 64053-00  
Fax : +49 (0)6196 64053-01  
E-mail : Sales-DE@TCIchemicals.com

**Tokyo Chemical Industry UK Ltd.**

Tel : +44 (0)1865 784560  
Fax : +44 (0)1865 784561  
E-mail : Sales-UK@TCIchemicals.com

**TCI Chemicals (India) Pvt. Ltd.**

Tel : 1800 425 7889 / 044-2262 0909  
Fax : 044-2262 8902  
E-mail : Sales-IN@TCIchemicals.com

**梯希爱(上海)化成工业发展有限公司**

Tel : 800-988-0390 / 021-67121386  
Fax : 021-6712-1385  
E-mail : Sales-CN@TCIchemicals.com

**TOKYO CHEMICAL INDUSTRY CO., LTD.**

Tel : +81 (0)3-5640-8878  
Fax : +81 (0)3-5640-8902  
E-mail : globalbusiness@TCIchemicals.com

Availability, price or specification of the listed products are subject to change without prior notice. Reproduction forbidden without the prior written consent of Tokyo Chemical Industry Co., Ltd.