



# Cell Imaging Reagents

## Fluorescent Stains

**DAPI-2HCl** [for Biochemical Research] (Blue Fluorescence) 5mg [A2412]  
**Goat Anti-Mouse IgG FITC Conjugate** (Green Fluorescence) 0.1mg/1vial [G0406]

### Application



- (A) HeLa cells nuclei stained with 1  $\mu\text{g}/\text{mL}$  **A2412** (Blue).
- (B) The HeLa cells were incubated with the Mouse Anti  $\alpha$ -Tubulin Antibody, followed by staining with the secondary antibody **G0406** (Green) at 10  $\mu\text{g}/\text{mL}$ . **A2412** was used to stain the cell nuclei (Blue) at a concentration of 1  $\mu\text{g}/\text{mL}$ .

Laser Scanning Microscope: Olympus FLUOVIEW FV 3000

### Related Products

<b>Goat Anti-Rabbit IgG FITC Conjugate</b> (Green Fluorescence)	0.1mg/1vial [G0452]
<b>Goat Anti-Mouse IgM FITC Conjugate</b> (Green Fluorescence)	0.1mg/1vial [G0453]
<b>Streptavidin FITC Conjugate</b> (Green Fluorescence)	0.1mg/1vial [S0966]
<b>Goat Anti-Mouse IgG DTBTA-Eu<sup>3+</sup> Conjugate</b> (Red Fluorescence)	0.1mg/1vial [G0505]
<b>Goat Anti-Rabbit IgG DTBTA-Eu<sup>3+</sup> Conjugate</b> (Red Fluorescence)	0.1mg/1vial [G0506]
<b>Streptavidin DTBTA-Eu<sup>3+</sup> Conjugate</b> (Red Fluorescence)	0.1mg/1vial [S0993]
<b>ATBTA-Eu<sup>3+</sup></b> [DTBTA-Eu <sup>3+</sup> Labeling Reagent] (Red Fluorescence)	10mg [A2083]
<b>Bisbenzimidide H 33258 Hydrate</b> [for Biochemical Research] (Blue Fluorescence)	25mg [H1343]

## Staining Dye Solution

### Methylene Blue Solution (Methanol Solution) [for Cell Staining] 100mL [M2392]

#### Application

- (1) Culture cells in a 6-well plate
- (2) Remove medium from the plate and wash it with PBS(-) twice
- (3) Remove PBS(-) from it, add 1mL of **M2392** and stain cells for 15 minutes
- (4) Remove **M2392** from it and wash it with deionized water twice






Figure. NIH/3T3 cells stained by the above method  
Please adjust staining time and volume according to cells. Because some cells need to be fixed separately, preliminary tests should be performed.

## Cell Proliferation Assay Reagents


### Resazurin (Ready-to-use solution) [for Cell proliferation assay] 25mL [R0195]

#### Mechanism



Resazurin can be used quantitatively determine cell proliferation, viability, and cytotoxicity. Resazurin, when added to viable cells, is reduced by the cellular enzymatic or chemical reactions converting blue/non-fluorescent resazurin to highly fluorescent resorufin. The assay is simple to perform since the indicator is water-soluble and has low toxicity, thus eliminating the washing/fixing and extraction steps required in other commonly used cell proliferation assays.

#### Cell viability assay



#### Application

1. Add **R0195** at a volume equal to 10% of the cell culture media volume.
2. Return cells to the incubator and continue the incubation for 2-24 hours\*.
3. Measure the fluorescent intensity using 540-570 nm excitation and 590 nm emission wavelengths. Absorbance can be measured using a spectrophotometer set at 570 nm.

Resazurin may be added at any time point during the culture period. For measurement of cell proliferation, it is best to add resazurin during the log phase of growth.

## Substrates for Reporter Assays

Chromogenic Substrates for  $\beta$ -Galactosidase

## Generating insoluble dye

<b>X-Gal</b> (5-Bromo-4-chloro-3-indolyl $\beta$ -D-Galactopyranoside) ■ blue	200mg / 1g [B3201]
<b>Magenta-Gal</b> (5-Bromo-6-chloro-3-indolyl $\beta$ -D-Galactopyranoside) ■ red-purple	20mg / 100mg [B3469]
<b>Bluo-Gal</b> (5-Bromo-3-indolyl $\beta$ -D-Galactopyranoside) ■ dark-blue	20mg / 100mg [B3470]
<b>Salmon-Gal</b> (6-Chloro-3-indolyl $\beta$ -D-Galactopyranoside) ■ bright red-purple	20mg / 100mg [C2371]

## Generating soluble dye

<b>ONPG</b> (2-Nitrophenyl $\beta$ -D-Galactopyranoside) ■ yellow	1g / 5g / 25g [N0418]
<b>PNPG</b> (4-Nitrophenyl $\beta$ -D-Galactopyranoside) ■ yellow	1g / 5g [N0616]

Chromogenic Substrates for  $\beta$ -Glucuronidase

## Generating insoluble dye

<b>X-Gluc CHA Salt</b> (5-Bromo-4-chloro-3-indolyl $\beta$ -D-Glucuronide Cyclohexylammonium Salt) ■ blue	10mg / 100mg [B3620]
<b>X-Gluc Sodium Salt</b> (5-Bromo-4-chloro-3-indolyl $\beta$ -D-Glucuronide Sodium Salt) ■ blue	10mg / 100mg [B3621]

## Chemiluminescent Substrates for Luciferase

<b>D-(-)-Luciferin</b>	10mg / 50mg [A5030]
<b>CLA</b>	10mg [A5307]
<b>MCLA</b>	10mg [A5309]
<b>FCLA Free Acid</b>	10mg [A5310]
<b>Red-CLA</b>	1mg [A5311]

## Chemiluminescence Reagent for the Detection of Superoxide

Lucigenin

1g / 5g [B1203]

MMT [= 10,10'-Dimethyl-9,9'-biacridinium Bis(monomethyl Terephthalate)]

100mg / 1g [B4339]

MMT (B4339) is a specific probe having lucigenin-like chemiluminescence to superoxide among reactive oxygen species. Since amphiphilic MMT which is less hydrophilic than lucigenin possesses cell-permeability, MMT has been applicable for the detection of intramitochondrial superoxide production.

## Application

Figure shows the localization of MMT in mitochondria.





Figure. Fluorescence in mitochondria in mouse peritoneal neutrophils. (Provided by Prof. Kobayashi)

S. Sasaki, S. Yamada, M. Iwamura, Y. Kobayashi, *Free Radic. Biol. Med.* **2013**, 65, 1005.

For further information please refer to our website at [www.TCIchemicals.com](http://www.TCIchemicals.com). 

## Ordering and Customer Service

## TCI AMERICA

Tel : 800-423-8616 / 503-283-1681  
Fax : 888-520-1075 / 503-283-1987  
E-mail : Sales-US@TCIchemicals.com

## TCI EUROPE N.V.

Tel : +32 (0)3 735 07 00  
Fax : +32 (0)3 735 07 01  
E-mail : Sales-EU@TCIchemicals.com

## TCI Deutschland GmbH

Tel : +49 (0)6196 64053-00  
Fax : +49 (0)6196 64053-01  
E-mail : Sales-DE@TCIchemicals.com

## Tokyo Chemical Industry UK Ltd.

Tel : +44 (0)1865 784560  
Fax : +44 (0)1865 784561  
E-mail : Sales-UK@TCIchemicals.com

## TCI Chemicals (India) Pvt. Ltd.

Tel : 1800 425 7889 / 044-2262 0909  
Fax : 044-2262 8902  
E-mail : Sales-IN@TCIchemicals.com

## 梯希爱(上海)化成工业发展有限公司

Tel : 800-988-0390 / 021-67121386  
Fax : 021-6712-1385  
E-mail : Sales-CN@TCIchemicals.com

## TOKYO CHEMICAL INDUSTRY CO., LTD.

Tel : +81 (0)3-5640-8878  
Fax : +81 (0)3-5640-8902  
E-mail : globalbusiness@TCIchemicals.com

Availability, price or specification of the listed products are subject to change without prior notice. Reproduction forbidden without the prior written consent of Tokyo Chemical Industry Co., Ltd.