


# Reducing Agents



Aluminum Hydrides

Boranes

Borohydrides


Metal Hydrides

Silanes

Other Reduction Reagents

# Reducing Agents

Reduction is a chemical reaction in which the target substances receive electrons, and is one of the most fundamental reactions in organic chemistry. Reduction reactions include the deoxygenation reaction and the hydrogenation reaction. Well-known reducing agents include metal hydrides<sup>1)</sup> such as lithium aluminum hydride (= LiAlH<sub>4</sub>) [L0203], boranes for hydride reduction, and hydrazine [H0172] used in the Wolff-Kishner reduction. A disconnection reaction of a disulfide moiety into two thiols is also considered a reduction.



This brochure introduces a variety of reducing agents and catalysts for reduction. We hope that this brochure will be useful for your research in organic synthesis. Catalysts for hydrogenation are introduced in another brochure, "Hydrogenation Catalysts".

Caution: Many reducing agents may spontaneously ignite on contact with air, or may react violently with water to produce flammable gases. Sufficient safety measures, such as using safety shields, wearing protective equipment, and using extreme caution should be taken when working with these reagents as well as in their disposal.


## ● Reduction of Carbonyl Groups and Imino Group

The table below shows the reactivities of each reducing agent toward carbonyl compounds and imines. Please make use of this table as a standard for reactions.

Reagents	Reactions					
	Imine $\text{R}-\text{C}(\text{NR})=\text{H}$ ↓ $\text{R}-\text{CH}_2-\text{NR}$	Aldehyde $\text{R}-\text{C}(=\text{O})-\text{H}$ ↓ $\text{R}-\text{CH}_2-\text{OH}$	Ketone $\text{R}-\text{C}(=\text{O})-\text{R}$ ↓ $\text{R}-\text{CH}(\text{OH})-\text{R}$	Ester $\text{R}-\text{C}(=\text{O})-\text{OR}$ ↓ $\text{R}-\text{CH}_2-\text{OH}$	Amide $\text{R}-\text{C}(=\text{O})-\text{NR}_2$ ↓ $\text{R}-\text{CH}_2-\text{NR}_2$	Carboxylic Acid $\text{R}-\text{C}(=\text{O})-\text{OH}$ ↓ $\text{R}-\text{CH}_2-\text{OH}$
<b>NaBH<sub>3</sub>CN</b> [S0396]	High	Middle	Middle	Low	Low	Low
<b>NaBH(OAc)<sub>3</sub></b> [S0394]	High	Middle	Middle	Low	Low	Low
<b>NaBH<sub>4</sub></b> [S0480]	High	High	High	Middle	Low	Low
<b>LiBH<sub>4</sub></b> [L0186]	High	High	High	High	Low	Low
<b>LiAlH<sub>4</sub></b> [L0170]	High	High	High	High	Middle	Low
<b>THF · BH<sub>3</sub></b> [T2346]	High	Middle	Middle	Low	High	High
<b>Me<sub>2</sub>S · BH<sub>3</sub></b> [D1843]	High	Middle	Middle	Low	High	High
<b>PhNEt<sub>2</sub> · BH<sub>3</sub></b> [D2581]	High	Middle	Middle	Low	High	High


Reactivity : High Middle Low

Sodium tris(1,1,1,3,3,3-hexafluoroisopropoxy)borohydride (= NaBH(HFIP)<sub>3</sub>) [S0810] is a selective reducing agent developed by Toshima *et al.* Aldehydes are selectively reduced in the presence of ketones and other reducible functions using S0810 to afford the corresponding primary alcohols in high yields.<sup>2)</sup>





## Reductive Amination

Reductive amination is a synthetic method that converts aldehydes and ketones into an amino group with an amine and a reducing agent. First, carbonyl groups and an amine form imines or iminium salts and subsequent nucleophilic attack by a reducing agent gives the amine moiety. The Leuckart reaction<sup>3)</sup> and Eschweiler-Clarke reaction<sup>4)</sup> are known as classical methods and formic acid [F0513] is used as a reducing agent in these reactions. Recently, sodium cyanoborohydride [S0396] has been frequently used for reductive amination, in what is called the Borch reaction.<sup>5)</sup> However, this method has a problem in that S0396 has strong toxicity due to the cyano group. Meanwhile, Kikukawa *et al.* have reported a new method using 2-picoline borane [B3018].<sup>6)</sup> B3018 is less toxic than S0396 and can be applied in both aqueous and neat conditions.



## Luche Reduction


When  $\alpha,\beta$ -unsaturated ketones are reduced with sodium borohydride [B0480], 1,4-reduction, not 1,2-reduction, preferentially occurs. However, Luche *et al.* have found that 1,2-reduction preferentially occurred by adding cerium(III) chloride [C2058].<sup>7)</sup> The difference in reactivity can be explained with HSAB theory. A “hard” nucleophile is required for 1,2-reduction, whereas a hydride of S0480 is originally regarded as a “soft” nucleophile, resulting in 1,4-reduction. However, it is considered that the nucleophile turns “hard” in the presence of C2058 and alcohol. In addition, since the cerium cation works as a Lewis acid, which promotes the electrophilicity of the carbonyl group, the 1,2-addition preferentially occurs. Incidentally, an aldehyde is not reduced under this method because it forms an acetal, which is inert under these conditions.



## Silane Reduction

Hydrosilanes are utilized in reduction as a hydride or a hydrogen radical source since the hydrogen atom has lower electronegativity than silicon. For instance, triethylsilane [T0662] is known as a reducing agent in the presence of a metal catalyst<sup>8)</sup> and Lewis acid<sup>9)</sup> and can reduce carbonyl groups and hydroxy groups into methylene moieties. T0662 is also used in the dehalogenation<sup>10)</sup> reaction and in the reduction of olefins.<sup>11)</sup>

Tris(trimethylsilyl)silane (= TTMS) [T1463] is mainly used as a hydrogen radical source in the dehalogenation reaction.<sup>12)</sup> In this manner, hydrosilanes and tin hydrides are currently used complementarily.




## ● Wolff-Kishner Reduction and Clemensen Reduction

Both Wolff-Kishner reduction<sup>13)</sup> and Clemensen reduction<sup>14)</sup> are known as reactions to convert carbonyl groups into methylene groups. The Wolff-Kishner reduction can reduce carbonyl groups to methylene moieties in the presence of hydrazine [H0172], strong base, and alcohol. This reaction conventionally requires strong base and intense heat conditions, but an improved method to use a Lewis acid catalyst and a silylhydrazine derivative (Myers modification) has been developed,<sup>15)</sup> which proceeds at room temperature.


The Clemensen reduction is operated by the treatment of zinc under acidic conditions. Zinc amalgam had been used at first, but a non-aqueous modification using zinc powder and hydrogen chloride in organic solvent has been developed<sup>16)</sup> and is utilized widely since mercury in amalgam is highly harmful.

The Wolff-Kishner reduction and Clemensen reduction are operated under basic and acidic conditions, respectively, so the two reactions can be used complementarily.


### Wolff-Kishner reduction



### Myers modification




### Clemensen reduction



## ● Meerwein-Ponndorf-Verley Reduction

In the Meerwein-Ponndorf-Verley (MPV) reduction, ketones and aldehydes can be reduced to alcohols by the treatment of aluminum isopropoxide [A0246] in isopropyl alcohol [I0163] solvent with heat.<sup>17)</sup> This reaction is an equilibrium reaction, so an excess amount of A0246 is needed to bias the reaction toward the desired alcohol. This reaction has an advantage in that it does not affect other functional groups. Following the first report, a modified method using a catalytic amount of samarium(II) iodide<sup>18)</sup> [S0494] instead of an excess amount of A0246 and asymmetric MPV reactions<sup>19)</sup> was reported. Incidentally, the Oppenauer oxidation<sup>20)</sup> is regarded as an opposite reaction of MPV reduction and it proceeds under acetone solvent to oxidize alcohols to carbonyl groups.




### Sml<sub>2</sub>-catalyzed MPV reduction




## ● Staudinger Reaction

The Staudinger reaction is utilized to convert an azide to an amine<sup>21)</sup> and can be regarded as a reduction in the sense of the addition of a hydrogen atom. When an azide is treated with triphenylphosphine [T0519], an iminophosphorane is formed with the elimination of a nitrogen molecule. The iminophosphorane is hydrolyzed to give the amine moiety. In contrast, the iminophosphorane gives an imine via the aza-Wittig reaction when treated with aldehydes or ketones.<sup>22)</sup>



## ● Reducing Agent to Disconnect Disulfide Bond

Tris(2-carboxyethyl)phosphine hydrochloride [T1656] can reduce a disulfide bond to give two thiols.<sup>23)</sup> In addition, T1656 can remove the oxygen atom on *N*-oxides and sulfoxides and can be applied to the Staudinger reaction.<sup>24)</sup> In this manner, T1656 shows interesting effects in the reductions of heteroatoms.



## References

- Review: J. Málek, M. Černý, *Synthesis* **1972**, 217.
- Y. Kuroiwa, S. Matsumura, K. Toshima, *Synlett* **2008**, 16, 2523.
- R. Leuckart, *Ber.* **1885**, 18, 2341.
- a) W. Eschweiler, *Ber.* **1905**, 38, 880.  
b) H. T. Clarke, H. B. Gillespie, S. Z. Weisshaus, *J. Am. Chem. Soc.* **1933**, 55, 4571.
- a) R. F. Borch, M. D. Bernstein, H. D. Durst, *J. Am. Chem. Soc.* **1971**, 93, 2897.  
b) R. F. Borch, A. I. Hassid, *J. Org. Chem.* **1972**, 37, 1673.
- S. Sato, T. Sakamoto, E. Miyazawa, Y. Kikugawa, *Tetrahedron* **2004**, 60, 7899.
- a) J. L. Luche, *J. Am. Chem. Soc.* **1978**, 100, 2226.  
b) J.-L. Luche, L. Rodriguez-Hahn, P. Crabbé, *J. Chem. Soc., Chem. Commun.* **1978**, 601.
- M. Mirza-Aghayan, M. Kalantari, R. Boukherroub, *Appl. Organomet. Chem.* **2019**, 33, e4837.
- I. Smonou, *Synth. Commun.* **1994**, 24, 1999.
- R. P. Allen, B. P. Roberts, C. R. Willis, *J. Chem. Soc., Chem. Commun.* **1989**, 1387.
- A. R. Jesus, A. P. Marques, A. P. Rauter, *Pure Appl. Chem.* **2016**, 88, 349.
- C. Chatgililoglu, D. Griller, M. Lesage, *J. Org. Chem.* **1988**, 53, 3641.
- a) Kishner, *N. J. Russ. Phys. Chem. Soc.* **1911**, 43, 582.  
b) L. Wolff, *Ann.* **1912**, 394, 86.
- a) E. Clemmensen, *Ber.* **1913**, 46, 1837.  
b) E. Clemmensen, *Ber.* **1914**, 47, 51.  
c) E. Clemmensen, *Ber.* **1914**, 47, 681.
- M. E. Furrow, A. G. Myers, *J. Am. Chem. Soc.* **2004**, 126, 5436.
- a) S. Yamamura, S. Ueda, Y. Hirata, *Chem. Commun.* **1967**, 1049.  
b) S. Yamamura, Y. Hirata, *J. Chem. Soc. C* **1968**, 2887.  
c) M. Toda, M. Hayashi, Y. Hirata, S. Yamamura, *Bull. Chem. Soc. Jpn.* **1972**, 45, 264.
- a) H. Meerwein, R. Schmidt, *Ann.* **1925**, 444, 221.  
b) A. Verley, *Bull. Soc. Chim. Fr.* **1925**, 37, 537.  
c) W. Ponnendorf, *Angew. Chem.* **1926**, 39, 138.
- J. L. Namy, J. Soupe, J. Collin, H. B. Kagan, *J. Org. Chem.* **1984**, 49, 2045.
- D. A. Evans, S. G. Nelson, M. R. Gagne, A. R. Muci, *J. Am. Chem. Soc.* **1993**, 115, 9800.
- R. V. Oppenauer, *Rec. Trav. Chim. Pays-Bas* **1937**, 56, 137.
- H. Staudinger, J. Meyer, *Helv. Chim. Acta* **1919**, 2, 635.
- T. Sasaki, S. Eguchi, T. Okano, *J. Am. Chem. Soc.* **1983**, 105, 5912.
- J. A. Burns, J. C. Butler, J. Moran, G. M. Whitesides, *J. Org. Chem.* **1991**, 56, 2648.
- A. M. Faucher, C. Grand-Maitre, *Synth. Commun.* **2003**, 33, 3503.

## Aluminum Hydrides


**L0203** 25g 100g


 Lithium Aluminum Hydride  
(Powder)  
CAS RN: 16853-85-3


**S0467** 25g 100g 500g


$$\text{NaAl}(\text{OCH}_2\text{CH}_2\text{OCH}_3)_2\text{H}_2$$
 Sodium Dihydridobis(2-methoxyethoxy)aluminate  
(70% in Toluene, ca. 3.6mol/L)  
CAS RN: 22722-98-1


## Boranes


**B5082** 5g 25g

 Borane - Ammonia  
Complex  
CAS RN: 13774-81-7


**B1264** 25g 100g

 Borane - *tert*-Butylamine  
Complex  
CAS RN: 7337-45-3


**D1842** 25g 100g 500g

 Dimethylamine Borane  
CAS RN: 74-94-2


**B5544** 25g

 Borane - Diisopropylamine  
Complex  
CAS RN: 55124-35-1


**B5545** 5g 25g

 Dicyclohexylamine  
Borane  
CAS RN: 131765-96-3


**T1181** 25g

 Trimethylamine Borane  
CAS RN: 75-22-9


**T1180** 25g

 Triethylamine Borane  
CAS RN: 1722-26-5

**D2581** 25g 100g

*N,N*-Diethylaniline Borane  
CAS RN: 13289-97-9

**B1569** 25mL

 Borane - Pyridine Complex  
CAS RN: 110-51-0


**B3018** 5g 25g

 Borane - 2-Methylpyridine  
Complex  
CAS RN: 3999-38-0


**M0898** 5g 25g

 Borane - Morpholine  
Complex  
CAS RN: 4856-95-5

**T2346** 100mL 500mL

 Borane - Tetrahydrofuran  
Complex (8.5% in  
Tetrahydrofuran, ca. 0.9mol/L)  
CAS RN: 14044-65-6

**D1843** 25mL 100mL

 Dimethyl Sulfide Borane  
CAS RN: 13292-87-0

**B1827** 1g

 1,2-Bis(*tert*-butylthio)-  
ethane Borane  
CAS RN: 71522-78-6

**T1789** 25g

 Triphenylphosphine  
Borane  
CAS RN: 2049-55-0

## Borohydrides

**S0480** 25g 100g 500g

 Sodium Borohydride  
CAS RN: 16940-66-2

**L0186** 100mL


 Lithium Borohydride  
(ca. 4mol/L in Tetrahydrofuran)  
CAS RN: 16949-15-8

**P1681** 25g 100g

 Potassium Borohydride  
CAS RN: 13762-51-1


**S0394** 25g 100g


 Sodium  
Triacetoxylborohydride  
CAS RN: 56553-60-7


**S0396** 5g 25g 250g

 Sodium  
Cyanoborohydride  
CAS RN: 25895-60-7


**S0810** 5g


$$\left[ \text{B}(\text{OCH}(\text{CF}_3)_2)_3 \right]^- \text{Na}^+$$
 Sodium Tris(1,1,1,3,3,3-  
hexafluoroisopropoxy)-  
borohydride  
CAS RN: 139494-68-1


**L0190** 100mL 500mL

 Lithium Triethylborohydride  
(ca. 12% in Tetrahydrofuran, ca.  
1.0mol/L)  
CAS RN: 22560-16-3


**P1439** 1g 5g

 Potassium Bis(1-pyrazolyl)-  
borohydride  
CAS RN: 18583-59-0



**P1440** 1g 5g

 Potassium Tris(1-pyrazolyl)-  
borohydride  
CAS RN: 18583-60-3

**P1719** 5g 25g

 Benzyltrimethylammonium  
Cyanoborohydride Resin  
cross-linked with 10% DVB  
(30-50mesh) (2.6-3.0mmol/g)


**T0852** 5g 25g

 Tetramethylammonium  
Borohydride  
CAS RN: 16883-45-7

**T0837** 5g 25g

 Tetraethylammonium  
Borohydride  
CAS RN: 17083-85-1




**B3128** 5g 25g

 Benzyltriethylammonium  
Borohydride  
CAS RN: 85874-45-9






<b>T0917</b> 5g 25g	<b>T1553</b> 5g 25g
 <p>Tetrabutylammonium Borohydride CAS RN: 33725-74-5</p>	 <p>Tetramethylammonium Triacetoxymethylborohydride CAS RN: 109704-53-2</p>





## Metal Hydrides

<b>S0481</b> 100g 500g	<b>Z0010</b> 1g 5g 25g
<p>NaH</p> <p>Sodium Hydride (60%, dispersion in Paraffin Liquid) CAS RN: 7646-69-7</p>	 <p>Zirconocene Chloride Hydride CAS RN: 37342-97-5</p>




## Silanes




<b>T0662</b> 25mL 250mL	<b>T0661</b> 5g 25g	<b>T1334</b> 10g
 <p>Triethylsilane CAS RN: 617-86-7</p>	 <p>Triphenylsilane CAS RN: 789-25-3</p>	 <p>Trihexylsilane CAS RN: 2929-52-4</p>




<b>T1533</b> 5mL 25mL 100mL	<b>T1896</b> 1g 5g	<b>D1825</b> 25mL	<b>D2100</b> 25mL 100mL	<b>D2196</b> 5mL 25mL
 <p>Triisopropylsilane CAS RN: 6485-79-6</p>	 <p>1,1,2,2-Tetraphenylsilane CAS RN: 16343-18-3</p>	 <p>Methylphenylsilane CAS RN: 776-76-1</p>	 <p>Dimethoxy(methyl)silane CAS RN: 16881-77-9</p>	 <p>Dimethylphenylsilane CAS RN: 766-77-8</p>

<b>D2403</b> 25mL	<b>D2406</b> 5g 25g	<b>D2820</b> 5g 25g	<b>P1291</b> 5mL 25mL
 <p>Diethoxymethylsilane CAS RN: 2031-62-1</p>	 <p>Diphenylsilane (&gt;97.0%) CAS RN: 775-12-2</p>	 <p>Diphenylsilane (&gt;98.0%) CAS RN: 775-12-2</p>	 <p>Phenylsilane CAS RN: 694-53-1</p>



## Other Reduction Reagents

<b>D1071</b> 1g 5g 25g	<b>D3775</b> 1g 5g 25g	<b>D4311</b> 1g 5g
 <p>DL-Dithiothreitol CAS RN: 3483-12-3</p>	 <p>Hantzsch Ester CAS RN: 1149-23-1</p>	 <p>Di-tert-butyl 1,4-Dihydro-2,6-dimethyl-3,5-pyridinedicarboxylate CAS RN: 55536-71-5</p>

<b>D5792</b> 5g 25g 100g	<b>F1202</b> 25mL 100mL	<b>F0513</b> 300mL	<b>S0807</b> 500g	<b>H0172</b> 25mL 500mL
<p>Na</p> <p>SD Super Fine™ (Sodium 25wt% dispersion in mineral oil) CAS RN: 7440-23-5</p>	 <p>TEAF CAS RN: 15077-13-1</p>	 <p>Formic Acid CAS RN: 64-18-6</p>	 <p>Sodium Formate CAS RN: 141-53-7</p>	<p>H<sub>2</sub>NNH<sub>2</sub> · H<sub>2</sub>O</p> <p>Hydrazine Monohydrate CAS RN: 7803-57-8</p>

<b>H1221</b> 300mL	<b>I0777</b> 1g 5g	<b>S0494</b> 25mL 100mL	<b>T0519</b> 25g 100g 500g	<b>T1656</b> 1g 5g 25g
<p>HI</p> <p>Hydriodic Acid (57%) CAS RN: 10034-85-2</p>	 <p>IPNBSH CAS RN: 6655-27-2</p>	<p>Sml<sub>2</sub></p> <p>Samarium(II) Iodide (ca. 0.1mol/L in Tetrahydrofuran) CAS RN: 32248-43-4</p>	 <p>Triphenylphosphine CAS RN: 603-35-0</p>	 <p>Tris(2-carboxyethyl)phosphine Hydrochloride CAS RN: 51805-45-9</p>

## Reducing Agents

T1819	5g	C2058	5g 25g	A0246	100g 500g
					
Triphenylphosphine Selenide CAS RN: 3878-44-2		CeCl <sub>3</sub> Cerium(III) Chloride Anhydrous CAS RN: 7790-86-5		Aluminum Isopropoxide CAS RN: 555-31-7	

### Ordering and Customer Service

#### TCI AMERICA

Tel : 800-423-8616 / 503-283-1681  
Fax : 888-520-1075 / 503-283-1987  
E-mail : Sales-US@TCIchemicals.com

#### TCI EUROPE N.V.

Tel : +32 (0)3 735 07 00  
Fax : +32 (0)3 735 07 01  
E-mail : Sales-EU@TCIchemicals.com

#### TCI Deutschland GmbH

Tel : +49 (0)6196 64053-00  
Fax : +49 (0)6196 64053-01  
E-mail : Sales-DE@TCIchemicals.com

#### Tokyo Chemical Industry UK Ltd.

Tel : +44 (0)1865 784560  
Fax : +44 (0)1865 784561  
E-mail : Sales-UK@TCIchemicals.com

#### TCI Chemicals (India) Pvt. Ltd.

Tel : 1800 425 7889 / 044-2262 0909  
Fax : 044-2262 8902  
E-mail : Sales-IN@TCIchemicals.com

#### 梯希爱(上海)化成工业发展有限公司

Tel : 800-988-0390 / 021-67121386  
Fax : 021-6712-1385  
E-mail : Sales-CN@TCIchemicals.com

#### TOKYO CHEMICAL INDUSTRY CO., LTD.

Tel : +81 (0)3-5640-8878  
Fax : +81 (0)3-5640-8902  
E-mail : globalbusiness@TCIchemicals.com

Availability, price or specification of the listed products are subject to change without prior notice. Reproduction forbidden without the prior written consent of Tokyo Chemical Industry Co., Ltd.