





ポリイミド合成材料

Raw Materials for Polyimide Synthesis



カルボン酸無水物

ジアミン





























ジイソシアン酸エステル



























ポリイミド合成材料

ポリイミドはイミド結合を含む高分子化合物の総称で、一般的に優れた耐熱性、耐薬品性、電気絶縁性、機械強度を有します。さらに熱に対する膨張率も有機化合物としては低く、高熱時にも膨張による誤差がわずかであるため、工業製品に多く利用されています。例えば、これらの性質を利用してポリイミドは電子









機器のフレキシブルプリント配線基板、電子回路の絶縁材料や半導体の保護膜などに広く利用されています。また、ポリイミドは液晶分子を一方方向に配列させるため、配向膜材料としても用いられています。以下に、ポリイミド合成に使用される製品を示します。









































カルボン酸無水物









































B1770 5g 25g		C2262 1g 5g		C2842 1g 5g	
 <p><i>meso</i>-Butane-1,2,3,4-tetracarboxylic Dianhydride CAS RN: 17309-39-6</p>		 <p>1,2,3,4-Cyclobutanetetracarboxylic Dianhydride CAS RN: 4415-87-6</p>		 <p>1,2,3,4-Cyclobutanetetracarboxylic Dianhydride (purified by sublimation) CAS RN: 4415-87-6</p>	
T2697 1g 5g		T2698 5g 25g		C0857 5g 25g	
 <p>1,2,3,4-Tetramethyl-1,2,3,4-cyclobutanetetracarboxylic Dianhydride CAS RN: 64198-16-9</p>		 <p>Tricyclo[6.4.0.0^{2,7}]dodecane-1,8:2,7-tetracarboxylic Dianhydride CAS RN: 738-90-9</p>		 <p>1,2,3,4-Cyclopentane-tetracarboxylic Dianhydride CAS RN: 6053-68-5</p>	
C2920 1g 5g		B6198 5g 25g		B1317 25g 500g	
 <p>1,2,3,4-Cyclopentane-tetracarboxylic Dianhydride (purified by sublimation) CAS RN: 6053-68-5</p>		 <p>BODA CAS RN: 2754-40-7</p>		 <p>Bicyclo[2.2.2]oct-7-ene-2,3,5,6-tetracarboxylic Dianhydride CAS RN: 1719-83-1</p>	
B4261 5g 25g		C3110 1g 5g		C2419 5g 25g	
 <p>Bicyclo[2.2.2]oct-7-ene-2,3,5,6-tetracarboxylic Dianhydride (purified by sublimation) CAS RN: 1719-83-1</p>		 <p>3-(Carboxymethyl)-1,2,4-cyclopentanetricarboxylic Acid 1,4:2,3-Dianhydride CAS RN: 6053-46-9</p>		 <p>1,2,4,5-Cyclohexanetetracarboxylic Dianhydride CAS RN: 2754-41-8</p>	
D5723 5g 25g		C3731 5g 25g		C2919 1g	
 <p>HBPDA CAS RN: 122640-83-9</p>		 <p>CpODA CAS RN: 1312366-93-0</p>		 <p>1,2,4,5-Cyclohexanetetracarboxylic Dianhydride (purified by sublimation) CAS RN: 2754-41-8</p>	
D5956 1g		D1901 25g 500g		D2192 5g 25g	
 <p>DDSQ Bis(bicyclo[2.2.1]heptane-2,3-dicarboxylic Dianhydride) CAS RN: 948849-07-8</p>		 <p>5-(2,5-Dioxotetrahydrofuryl)-3-methyl-3-cyclohexene-1,2-dicarboxylic Anhydride CAS RN: 73003-90-4</p>		 <p>4-(2,5-Dioxotetrahydrofuran-3-yl)-1,2,3,4-tetrahydronaphthalene-1,2-dicarboxylic Anhydride CAS RN: 13912-65-7</p>	
B0040 25g 100g 500g		P2103 5g 25g		B1326 25g 100g	
 <p>Pyromellitic Dianhydride CAS RN: 89-32-7</p>		 <p>Pyromellitic Dianhydride (purified by sublimation) CAS RN: 89-32-7</p>		 <p>4,4'-Biphtalic Anhydride CAS RN: 2420-87-3</p>	
B4262 5g		B4629 5g		B0948 25g 100g 500g	
 <p>4,4'-Biphtalic Anhydride (purified by sublimation) CAS RN: 2420-87-3</p>		 <p>3,4'-Biphtalic Anhydride CAS RN: 36978-41-3</p>		 <p>4,4'-Carbonyldiphthalic Anhydride CAS RN: 2421-28-5</p>	
B4260 5g 25g		E1164 5g 25g		E1412 1g 5g	
 <p>4,4'-Carbonyldiphthalic Anhydride (purified by sublimation) CAS RN: 2421-28-5</p>		 <p>4,4'-(Ethyne-1,2-diyl)-diphthalic Anhydride CAS RN: 129808-00-0</p>		 <p>4,4'-(Ethyne-1,2-diyl)diphthalic Anhydride (purified by sublimation) CAS RN: 129808-00-0</p>	
H0771 5g 25g 100g		 <p>4,4'-(Hexafluoroisopropylidene)-diphthalic Anhydride CAS RN: 1107-00-2</p>			

















<p>H1438 1g 5g</p>  <p>4,4'-(Hexafluoroisopropylidene)-diphthalic Anhydride (purified by sublimation) CAS RN: 1107-00-2</p>	<p>F1192 5g 25g</p>  <p>9,9-Bis(3,4-dicarboxyphenyl)-fluorene Dianhydride CAS RN: 135876-30-1</p>	<p>O0237 25g 100g</p>  <p>4,4'-Oxydiphthalic Anhydride CAS RN: 1823-59-2</p>	<p>O0423 5g</p>  <p>4,4'-Oxydiphthalic Anhydride (purified by sublimation) CAS RN: 1823-59-2</p>	<p>O0384 1g 25g</p>  <p>3,4'-Oxydiphthalic Anhydride CAS RN: 50662-95-8</p>
<p>N1128 1g 5g</p>  <p>2,3,6,7-NTCDA CAS RN: 3711-01-1</p>	<p>N1247 1g</p>  <p>1,2,5,6-NTCDA CAS RN: 3711-03-3</p>	<p>N0369 25g 250g</p>  <p>NTCDA CAS RN: 81-30-1</p>	<p>N0755 1g 5g</p>  <p>NTCDA (purified by sublimation) CAS RN: 81-30-1</p>	<p>P0972 25g 100g 500g</p>  <p>3,4,9,10-Perylene-tetracarboxylic Dianhydride CAS RN: 128-69-8</p>
<p>E1413 5g 25g</p>  <p>TMEG CAS RN: 1732-96-3</p>	<p>P2625 5g 25g</p>  <p>TAHQ CAS RN: 2770-49-2</p>	<p>I0856 25g 100g</p>  <p>4,4'-(4,4'-Isopropylidene-diphenoxy)diphthalic Anhydride CAS RN: 38103-06-9</p>	<p>T4007 1g 5g</p>  <p>6-FAP-ATA CAS RN: 223255-30-9</p>	<p>E0480 25g 100g</p>  <p>EDTA Dianhydride CAS RN: 23911-25-3</p>
<p>M3617 1g 5g</p>  <p>Mellitic Trianhydride CAS RN: 4253-24-1</p>	<p>P1614 500g</p>  <p>Phthalic Anhydride CAS RN: 85-44-9</p>	<p>T2381 1g 5g</p>  <p>Tetrafluorophthalic Anhydride CAS RN: 652-12-0</p>	<p>T0071 25g 500g</p>  <p>Tetrachlorophthalic Anhydride CAS RN: 117-08-8</p>	<p>T0046 25g 100g</p>  <p>Tetrabromophthalic Anhydride CAS RN: 632-79-1</p>
<p>T3871 5g</p>  <p>Tetrabromophthalic Anhydride (purified by sublimation) CAS RN: 632-79-1</p>	<p>E0579 1g 5g</p>  <p>4-Ethynylphthalic Anhydride CAS RN: 73819-76-8</p>	<p>M2757 5g 25g</p>  <p>4-(1-Propynyl)phthalic Anhydride CAS RN: 1240685-26-0</p>	<p>P1437 5g 25g</p>  <p>4-Phenylethynylphthalic Anhydride CAS RN: 119389-05-8</p>	<p>C0046 25g 500g</p>  <p>Trimellitic Anhydride CAS RN: 552-30-7</p>
<p>T1079 25g 100g 500g</p>  <p>Trimellitic Anhydride Chloride CAS RN: 1204-28-0</p>				

ジアミン
















<p>B1083 5g 25g</p>  <p>1,4-Bis(aminomethyl)-cyclohexane (<i>cis</i>- and <i>trans</i>-mixture) CAS RN: 2549-93-1</p>	<p>B3852 25g 500g</p>  <p>Bis(aminomethyl)norbornane (mixture of isomers) CAS RN: 56602-77-8</p>	<p>C0813 5mL 25mL</p>  <p>1,3-Cyclohexanediamine (<i>cis</i>- and <i>trans</i>-mixture) CAS RN: 3385-21-5</p>	<p>C0814 5mL 25mL</p>  <p>1,4-Cyclohexanediamine (<i>cis</i>- and <i>trans</i>-mixture) CAS RN: 3114-70-3</p>	<p>B1005 25mL 500mL</p>  <p>1,3-Bis(aminomethyl)-cyclohexane (<i>cis</i>- and <i>trans</i>-mixture) CAS RN: 2579-20-6</p>
<p>M0699 25g 500g</p>  <p>4,4'-Methylenebis(cyclohexylamine) (mixture of isomers) CAS RN: 1761-71-3</p>	<p>M1422 25mL 500mL</p>  <p>4,4'-Methylenebis(2-methylcyclohexylamine) (mixture of isomers) CAS RN: 6864-37-5</p>	<p>I0228 25mL 500mL</p>  <p>Isophoronediamine (<i>cis</i>- and <i>trans</i>-mixture) CAS RN: 2855-13-2</p>		

<p>P1892 25g 500g</p>  <p>1,3-Phenylenediamine CAS RN: 108-45-2</p>	<p>P1911 25g 250g</p>  <p>1,3-Phenylenediamine CAS RN: 108-45-2</p>	<p>P0170 25g 250g</p>  <p>1,4-Phenylenediamine CAS RN: 106-50-3</p>	<p>D2183 25g 100g</p>  <p>2,5-Dimethyl-1,4-phenylenediamine CAS RN: 6393-01-7</p>	<p>T1457 5g 25g</p>  <p>2,3,5,6-Tetramethyl-1,4-phenylenediamine CAS RN: 3102-87-2</p>
<p>T2019 1g 5g</p>  <p>2,4,5,6-Tetrafluoro-1,3-phenylenediamine CAS RN: 1198-63-6</p>	<p>T1110 1g</p>  <p>2,3,5,6-Tetrafluoro-1,4-phenylenediamine CAS RN: 1198-64-7</p>	<p>T3239 1g 5g</p>  <p>5-Trifluoromethyl-1,3-phenylenediamine CAS RN: 368-53-6</p>	<p>D1888 5g 25g</p>  <p>4,6-Diaminoresorcinol Dihydrochloride CAS RN: 16523-31-2</p>	<p>D2022 5g 25g</p>  <p>2,5-Diamino-1,4-benzenedithiol Dihydrochloride CAS RN: 75464-52-7</p>
<p>A1431 5g 25g</p>  <p>3-Aminobenzylamine CAS RN: 4403-70-7</p>	<p>A1436 25g 250g</p>  <p>4-Aminobenzylamine CAS RN: 4403-71-8</p>	<p>D0127 25g 100g 500g</p>  <p><i>m</i>-Xylylenediamine CAS RN: 1477-55-0</p>	<p>D1018 25g 100g 500g</p>  <p><i>p</i>-Xylylenediamine CAS RN: 539-48-0</p>	<p>D0101 25g 100g 500g</p>  <p>1,5-Diaminonaphthalene CAS RN: 2243-62-1</p>
<p>D5689 1g 5g</p>  <p>1,5-Diaminonaphthalene (purified by sublimation) CAS RN: 2243-62-1</p>	<p>D0092 1g 5g 25g</p>  <p>2,7-Diaminofluorene CAS RN: 525-64-4</p>	<p>D3180 25g</p>  <p>2,6-Diaminoanthraquinone CAS RN: 131-14-6</p>	<p>T3098 25g</p>  <p><i>m</i>-Tolidine CAS RN: 84-67-3</p>	<p>T0253 1g 5g</p>  <p><i>o</i>-Tolidine CAS RN: 119-93-7</p>
<p>T1023 1g 5g 25g</p>  <p>TMB CAS RN: 54827-17-7</p>	<p>D1344 25g 100g</p>  <p><i>o</i>-Dianisidine CAS RN: 119-90-4</p>	<p>M0609 25g 500g</p>  <p>4,4'-Methylenebis(2-chloroaniline) CAS RN: 101-14-4</p>	<p>B0395 25g 500g</p>  <p>2,2'-Benzidinedisulfonic Acid (contains 30% Water at maximum) CAS RN: 117-61-3</p>	<p>D0077 5g 25g</p>  <p>3,3'-Diaminobenzidine CAS RN: 91-95-2</p>
<p>B1711 1g 5g 25g</p>  <p>2,2'-Bis(trifluoromethyl)benzidine CAS RN: 341-58-2</p>	<p>D1632 1g 5g</p>  <p>4,4'-Diaminooctafluorobiphenyl CAS RN: 1038-66-0</p>	<p>D0754 1g 5g</p>  <p>4,4'-Oxybis[3-(trifluoromethyl)aniline] CAS RN: 13138-48-2</p>	<p>D3390 1g 5g</p>  <p>4,4''-Diamino-<i>p</i>-terphenyl CAS RN: 3365-85-3</p>	<p>D1683 1g 5g</p>  <p>3,3'-Diaminodiphenylmethane CAS RN: 19471-12-6</p>
<p>D1684 1g 5g 25g</p>  <p>3,4'-Diaminodiphenylmethane CAS RN: 19430-83-2</p>	<p>M0220 25g 500g</p>  <p>4,4'-Diaminodiphenylmethane CAS RN: 101-77-9</p>	<p>D2046 25g</p>  <p>4,4'-Diamino-3,3'-dimethyldiphenylmethane CAS RN: 838-88-0</p>	<p>M3403 5g 25g</p>  <p>4,4'-Methylenebis(2,6-dimethylaniline) CAS RN: 4073-98-7</p>	<p>M1238 25g</p>  <p>4,4'-Methylenebis(2-ethyl-6-methylaniline) CAS RN: 19900-72-2</p>
<p>M1897 25g 500g</p>  <p>4,4'-Methylenebis(2,6-diethylaniline) CAS RN: 13680-35-8</p>	<p>D1682 1g 5g</p>  <p>3,3'-Diaminobenzophenone CAS RN: 611-79-0</p>	<p>D1561 1g 5g</p>  <p>4,4'-Diaminobenzophenone CAS RN: 611-98-3</p>	<p>E0346 5g 25g</p>  <p>4,4'-Ethylenedianiline CAS RN: 621-95-4</p>	<p>D1658 5g</p>  <p>4,4'-Diamino-2,2'-dimethylbibenzyl CAS RN: 54628-21-6</p>

<p>D0121 25g 100g 500g</p>  <p>4,4'-Diaminostilbene-2,2'-disulfonic Acid CAS RN: 81-11-8</p>	<p>B4429 5g 25g</p>  <p>2,2-Bis(3-amino-4-hydroxyphenyl)propane CAS RN: 1220-78-6</p>	<p>B1400 1g 5g 25g</p>  <p>2,2-Bis(4-aminophenyl)hexafluoropropane CAS RN: 1095-78-9</p>	<p>B1415 5g 25g</p>  <p>2,2-Bis(3-amino-4-hydroxyphenyl)hexafluoropropane CAS RN: 83558-87-6</p>	<p>B1562 25g 500g</p>  <p>1,4-Bis[2-(4-aminophenyl)-2-propyl]benzene CAS RN: 2716-10-1</p>
<p>B3853 25g</p>  <p>1,3-Bis[2-(4-aminophenyl)-2-propyl]benzene CAS RN: 2687-27-6</p>	<p>B2076 5g 25g</p>  <p>1,1-Bis(4-aminophenyl)cyclohexane CAS RN: 3282-99-3</p>	<p>B1549 25g 100g</p>  <p>9,9-Bis(4-aminophenyl)fluorene CAS RN: 15499-84-0</p>	<p>B2654 1g</p>  <p>9,9-Bis(4-aminophenyl)fluorene (purified by sublimation) CAS RN: 15499-84-0</p>	<p>B2690 5g</p>  <p>9,9-Bis(4-amino-3-chlorophenyl)fluorene CAS RN: 107934-68-9</p>
<p>B2693 5g 25g</p>  <p>9,9-Bis(4-amino-3-methylphenyl)fluorene CAS RN: 107934-60-1</p>	<p>D1924 5g 25g</p>  <p>3,4'-Diaminodiphenyl Ether CAS RN: 2657-87-6</p>	<p>O0088 25g 500g</p>  <p>4,4'-Diaminodiphenyl Ether CAS RN: 101-80-4</p>	<p>O0568 1g</p>  <p>3,3'-Oxydianiline CAS RN: 15268-07-2</p>	<p>O0542 5g 25g</p>  <p>4,4'-Oxybis[3-(trifluoromethyl)aniline] CAS RN: 344-48-9</p>
<p>B1370 5g 25g</p>  <p>1,3-Bis(3-aminophenoxy)benzene CAS RN: 10526-07-5</p>	<p>B1680 25g</p>  <p>1,3-Bis(4-aminophenoxy)benzene CAS RN: 2479-46-1</p>	<p>B1243 5g 25g</p>  <p>1,4-Bis(4-aminophenoxy)benzene CAS RN: 3491-12-1</p>	<p>B4428 5g 25g</p>  <p>1,4-Bis(4-amino-2-trifluoromethylphenoxy)benzene CAS RN: 94525-05-0</p>	<p>B1671 5g 25g</p>  <p>4,4'-Bis(4-aminophenoxy)biphenyl CAS RN: 13080-85-8</p>
<p>B6310 1g 5g</p>  <p>ET-2,7-Na CAS RN: 132257-69-3</p>	<p>B6326 5g 25g</p>  <p>43BAPOBP CAS RN: 105112-76-3</p>	<p>B6320 1g 5g</p>  <p>6FAPBP CAS RN: 138321-99-0</p>	<p>B6311 1g 5g</p>  <p>ABHQ CAS RN: 22095-98-3</p>	<p>B1551 25g 250g</p>  <p>2,2-Bis[4-(4-aminophenoxy)phenyl]propane CAS RN: 13080-86-9</p>
<p>B1398 5g 25g</p>  <p>2,2-Bis[4-(4-aminophenoxy)phenyl]hexafluoropropane CAS RN: 69563-88-8</p>	<p>B2928 25g</p>  <p>Bis(2-aminophenyl) Sulfide CAS RN: 5873-51-8</p>	<p>T0632 25g</p>  <p>Bis(4-aminophenyl) Sulfide CAS RN: 139-65-1</p>	<p>B0394 25g 500g</p>  <p>Bis(3-aminophenyl) Sulfone CAS RN: 599-61-1</p>	<p>D0089 25g 100g 500g</p>  <p>Bis(4-aminophenyl) Sulfone CAS RN: 80-08-0</p>
<p>B2859 5g 25g</p>  <p>Bis(3-amino-4-hydroxyphenyl) Sulfone CAS RN: 7545-50-8</p>	<p>B1681 25g</p>  <p>Bis[4-(3-aminophenoxy)phenyl] Sulfone CAS RN: 30203-11-3</p>	<p>B1682 25g</p>  <p>Bis[4-(4-aminophenoxy)phenyl] Sulfone CAS RN: 13080-89-2</p>	<p>T1419 25g</p>  <p>o-Tolidine Sulfone (contains 2,6-Dimethyl isomer) CAS RN: 55011-44-4</p>	<p>D2116 1g 5g</p>  <p>3,6-Diaminocarbazole CAS RN: 86-71-5</p>
<p>T2728 5g 25g</p>  <p>1,3,5-Tris(4-aminophenyl)benzene CAS RN: 118727-34-7</p>	<p>B4541 5g 25g</p>  <p>Bis(4-aminophenyl) Terephthalate CAS RN: 16926-73-1</p>	<p>E0077 25mL 500mL</p>  <p>Ethylenediamine Anhydrous CAS RN: 107-15-3</p>	<p>D0113 25mL 500mL</p>  <p>1,2-Diaminopropane CAS RN: 78-90-0</p>	<p>M1657 1g 5g</p>  <p>2-Methyl-1,3-propanediamine CAS RN: 2400-78-4</p>

D0239 25g 400g  1,4-Diaminobutane CAS RN: 110-60-1	D2102 5mL 25mL  1,2-Diamino-2-methylpropane CAS RN: 811-93-8	D0108 5mL 25mL  1,5-Diaminopentane CAS RN: 462-94-2	D2252 25mL 500mL  1,3-Diaminopentane CAS RN: 589-37-7	D0095 25g 500g  1,6-Diaminohexane CAS RN: 124-09-4
D0094 5g 25g  1,7-Diaminoheptane CAS RN: 646-19-5	D0085 25g 100g  1,10-Diaminodecane CAS RN: 646-25-3	D2075 1g 5g  1,11-Diaminoundecane CAS RN: 822-08-2	D0091 25g 100g 500g  1,12-Diaminododecane CAS RN: 2783-17-7	B1201 25mL 500mL  1,4-Butanediol Bis(3-aminopropyl) Ether CAS RN: 7300-34-7
D1571 25mL 100mL 500mL  Diethylene Glycol Bis(3-aminopropyl) Ether CAS RN: 4246-51-9	M1085 25mL  2,2'-Diamino-N-methyldiethylamine CAS RN: 4097-88-5	B0821 25mL 500mL  3,3'-Diamino-N-methylpropylamine CAS RN: 105-83-9	T1243 5mL 25mL 100mL  Tris(2-aminoethyl)amine CAS RN: 4097-89-6	T1540 5mL 25mL  Tris(3-aminopropyl)amine CAS RN: 4963-47-7
B1552 5g 25g  1,3-Bis(3-aminopropyl)-tetramethyldisiloxane CAS RN: 2469-55-8				

ジイソシアン酸エステル

ジイソシアン酸エステル					B1538 25mL 500mL  1,3-Bis(isocyanatomethyl)-cyclohexane (<i>cis</i> - and <i>trans</i> -mixture) CAS RN: 38661-72-2	I0314 25mL 500mL  Isophorone Diisocyanate (mixture of isomers) CAS RN: 4098-71-9	D2380 25g 100g 500g  3,3'-Dimethylnaphthidine CAS RN: 5124-30-1
					P1092 1g 5g  1,3-Phenylene Diisocyanate CAS RN: 123-61-5	P1640 25g 250g  1,4-Phenylene Diisocyanate CAS RN: 104-49-4	T0263 25g 500g  Toluene-2,4-diisocyanate CAS RN: 584-84-9
X0022 25g 100g 500g  <i>m</i> -Xylylene Diisocyanate CAS RN: 3634-83-1	B1442 25g 500g  1,3-Bis(2-isocyanato-2-propyl)benzene CAS RN: 2778-42-9	N0168 25g 100g  1,5-Diisocyanatonaphthalene CAS RN: 3173-72-6	D2051 25g 250g  4,4'-Diisocyanato-3,3'-dimethylbiphenyl CAS RN: 91-97-4	D1541 5g 25g  3,3'-Dichloro-4,4'-diisocyanatobiphenyl CAS RN: 5331-87-3			
D0897 25g 500g  Methylenediphenyl 4,4'-Diisocyanate CAS RN: 101-68-8	M3310 100g  Methylenediphenyl 4,4'-Diisocyanate (Granulated) CAS RN: 101-68-8	B1509 1g  2,2-Bis(4-isocyanatophenyl)-hexafluoropropane CAS RN: 10224-18-7	H1842 5g 25g  Hexamethylene Diisocyanate Isocyanurate Trimer CAS RN: 3779-63-3				

東京化成工業株式会社

試薬製品について

■本社営業部 〒103-0001 東京都中央区日本橋小伝馬町 16-12 T-PLUS 日本橋小伝馬町8階
Tel: 03-3668-0489 Fax: 03-3668-0520 E-mail: Sales-JP@TCIchemicals.com

■大阪営業部 〒541-0041 大阪府大阪市中央区北浜1-1-21 第2中井ビル1階
Tel: 06-6228-1155 Fax: 06-6228-1158 E-mail: osaka-s@TCIchemicals.com

スケールアップ、受託サービス(合成・開発・製造)について

□化成品営業部 〒103-0001 東京都中央区日本橋小伝馬町 16-12 T-PLUS 日本橋小伝馬町8階
Tel: 03-5651-5171 Fax: 03-5640-8021 E-mail: finechemicals@TCIchemicals.com

弊社製品取扱店

本誌掲載の化学品は試験・研究用にのみ使用するものです。化学知識のある専門家以外の方のご使用はお避けください。品目や製品情報等、掲載内容の変更を予告なく行う場合があります。内容の一部または全部の無断転載・複製はご遠慮ください。