



乙烯基化试剂


乙烯基硅烷衍生物




Vinyltrimethylsilane
25mL / 100mL
[V0067]




Triethylvinylsilane
5mL / 25mL
[T1703]



Vinyltrimethoxysilane
25mL / 100mL / 500mL
[V0042]




Triethoxyvinylsilane
25mL / 100mL / 500mL
[V0044]



2,4,6,8-Tetramethyl-2,4,6,8-tetravinylcyclotetrasiloxane
5g / 25g
[T2523]




(1-Fluorovinyl)methyldiphenylsilane
1g / 5g
[F0518]



V0067



Y. 98%

TASF: Tris(diethylamino)sulfonium Difluorotrimethylsilicate

Reference: Y. Hatanaka, T. Hiyama, *J. Org. Chem.* **1988**, 53, 918.

V0044

Y. ~ 98%


R = COMe, CO₂H
OMe, NH₂Reference: A. Gordillo, M. A. Ortuño, C. L.-Mardomingo, A. Lledos, G. Ujaque, E. de Jesus, *J. Am. Chem. Soc.* **2013**, 135, 13749.Reference: S. E. Denmark, C. R. Butler, *Org. Lett.* **2006**, 8, 63.

F0518 (1.5 eq.)

Y. 86%

Reference: T. Hanamoto, T. Kobayashi, *J. Org. Chem.* **2003**, 68, 6354.



乙烯基硼烷衍生物




2,4,6-Trivinylboroxin - Pyridine Complex
1g / 5g
[T2498]





Potassium Vinyltrifluoroborate
1g / 5g
[P1479]

Reference: F. Kerins and D. F. O'Shea, *J. Org. Chem.* **2002**, 67, 4968.R = CO₂Et or NO₂R = CO₂Et, Y. 84%
NO₂, Y. 88%Reference: S. Darses, G. Michaud, J.-P. Genet, *Eur. J. Org. Chem.* **1999**, 1875.




4,4,5,5-Tetramethyl-2-vinyl-1,3,2-dioxaborolane
(stabilized with Phenothiazine)

1g / 5g / 25g
[T2297]




Dibutyl Vinylboronate
(stabilized with Phenothiazine)

1g / 5g
[D3596]




其它



Tributylvinyltin

1g / 5g / 25g
[T1794]



更多信息，请查看我们的主页：www.TCIchemicals.com

乙烯基化

