


Conducting Polymer Research Reagents



Polymers

Thiophenes

Pyrroles

Carbazoles

Fluorenes

Aryls

Dichlorosilanes

Others



Conducting Polymer Research Reagents

In 1977, Shirakawa and co-workers reported that thin films of the semi-conducting polymer polyacetylene show a dramatic increase in electrical conductivity when doped with controlled amounts of iodine.¹⁾ Their reports triggered intensive R&D into the electrical conductivity of plastic materials. Studies of conducting polymers began and, as a result, many π -conjugated polymers such as polypyrrole²⁾, polythiophene³⁾, polyaniline⁴⁾ and polyphenylenevinylene⁵⁾ have been developed. Among these polymers, many of them have been utilized commercially. One example is an electrolytic condenser using polypyrrole. This condenser has characteristic features such as being compact and lightweight, as well as having high-capacity and high-frequency compliance. These features have achieved a downsizing and weight-saving in electronic devices, and are utilized currently in some mobile phones. Thereby, conducting polymers are widely used for electronic devices vital to our everyday life. These and many other achievements were the reasons Shirakawa received the Nobel Prize Award in 2000.

Recent studies have found that polyphenylenevinylene shows electroluminescence. Also the synthesis of a spiral



structure helical polyacetylene has been achieved.

Further active studies are underway, looking to discover new features and properties of these conducting polymers.

References













- 1) H. Shirakawa, E. J. Louis, A. G. MacDiarmid, C. K. Chiang, A. J. Heeger, *J. Chem. Soc. Chem. Commun.* **1977**, 578.
- 2) K. K. Kanazawa, A. F. Diaz, R. H. Geiss, W. D. Gill, J. F. Kwak, J. A. Logan, J. F. Rabolt, G. B. Street, *J. Chem. Soc. Chem. Commun.* **1979**, 854; M. Salmon, A. F. Diaz, J. A. Logan, K. Krounbi, J. Bargon, *Mol. Cryst. Liq. Cryst.* **1982**, 83, 265; M. Ogasawara, K. Funahashi, T. Demura, T. Hagiwara, K. Iwata, *Synth. Met.* **1986**, 14, 62.
- 3) M. Sato, S. Tanaka, K. Kaeriyama, *J. Chem. Soc. Chem. Commun.* **1985**, 713; S. Hotta, *Synth. Met.* **1987**, 22, 103.
- 4) T. Kobayashi, H. Yoneyama, H. Tamura, *J. Electroanal. Chem.* **1984**, 177, 281; J.-C. Chiang, A. G. MacDiarmid, *Synth. Met.* **1986**, 13, 193.
- 5) G. E. Wnek, J. C. W. Chien, F. E. Karasz, C. P. Lillya, *Polymer* **1979**, 20, 1441; D. R. Ganan, J. D. Capistran, F. E. Karasz, R. W. Lenz, S. Antoun, *Polymer* **1987**, 28, 567.








































Polymers











































Thiophenes









































3-Alkyl Thiophenes







































T3104 5g 25g  3-Tetradecylthiophene [110851-66-6]	H1462 1g 5g 25g  3-Hexadecylthiophene [119269-24-8]	O0245 1g 5g  3-Octadecylthiophene [104934-54-5]	E0973 200mg 1g  3-(2-Ethylhexyl)-thiophene [121134-38-1]	C1801 5g 25g  2-Chloro-3-methylthiophene [14345-97-2]
D3680 1g 5g  2,5-Dichloro-3-methylthiophene [17249-90-0]	B1025 5g 25g  2-Bromo-3-methylthiophene [14282-76-9]	B4273 1g  2-Bromo-4-methylthiophene [53119-60-1]	D4304 5g 25g  2,5-Dibromo-3-methylthiophene [13191-36-1]	C2613 5g  2-Chloro-3-hexylthiophene [817181-75-2]
B3507 5g 25g  2-Bromo-3-hexylthiophene [69249-61-2]	B3985 1g  2-Bromo-4-hexylthiophene [210705-84-3]	D3896 5g 25g  2,5-Dibromo-3-hexylthiophene [116971-11-0]	B3865 1g 5g  2-Bromo-3-hexyl-5-iodothiophene (stabilized with Copper chip) [160096-76-4]	H1294 1g 5g  4-Hexyl-2-(4,4,5,5-tetramethyl-1,3,2-dioxaborolan-2-yl)-thiophene [883742-29-8]
H1298 1g 5g 25g  3-Hexyl-2-(4,4,5,5-tetramethyl-1,3,2-dioxaborolan-2-yl)-thiophene [850881-09-3]	B3733 5g  2-Bromo-3-n-octylthiophene [145543-83-5]	D4053 5g  2,5-Dibromo-3-n-octylthiophene [149703-84-4]	B4584 1g  2-Bromo-5-iodo-3-n-octylthiophene [1085181-82-3]	O0411 1g 5g  4-n-Octyl-2-(4,4,5,5-tetramethyl-1,3,2-dioxaborolan-2-yl)-thiophene [405165-12-0]
O0410 1g 5g  3-n-Octyl-2-(4,4,5,5-tetramethyl-1,3,2-dioxaborolan-2-yl)-thiophene [405165-14-2]	B4783 1g 5g  2-Bromo-3-(2-ethylhexyl)-thiophene [303734-52-3]	B4585 1g 5g  2-Bromo-3-(2-ethylhexyl)-5-iodothiophene (stabilized with Copper chip) [1034352-30-1]	D4572 5g 25g  2,5-Dibromo-3-decylthiophene [158956-23-1]	B3691 1g 5g  2-Bromo-3-dodecylthiophene [139100-06-4]
D4031 5g  2,5-Dibromo-3-dodecylthiophene [148256-63-7]	B4586 1g  2-Bromo-3-dodecyl-5-iodothiophene (stabilized with Copper chip) [306769-48-2]	D4291 1g 5g  4-Dodecyl-2-(4,4,5,5-tetramethyl-1,3,2-dioxaborolan-2-yl)-thiophene [1173788-58-3]	D4290 1g 5g  3-Dodecyl-2-(4,4,5,5-tetramethyl-1,3,2-dioxaborolan-2-yl)-thiophene [960524-18-9]	<p style="text-align: center; font-size: 2em; color: #e91e63; font-weight: bold;">Other Thiophenes</p>
E0915 1g 5g  2-Ethynylthiophene [4298-52-6]	E0892 1g 5g  3-Ethynylthiophene [67237-53-0]	C1344 5g 25g  3-Cyanothiophene [1641-09-4]	M1093 1g 5g 25g  3-Methoxythiophene [17573-92-1]	C1192 5g 25g  3-Chlorothiophene [17249-80-8]
B1067 25g 250g  3-Bromothiophene [872-31-1]	I0380 1g 5g  3-Iodothiophene [10486-61-0]	B1541 5g 25g  3-Bromo-4-methylthiophene [30318-99-1]	B3385 5g 25g  4-Bromo-2-methylthiophene [29421-92-9]	B1540 5g 25g  2-Bromo-5-methylthiophene (stabilized with Copper chip + NaHCO ₃) [765-58-2]
































Conducting Polymer Research Reagents



































<p>B2934 1g 5g 25g</p>  <p>3-Bromothiophene-2-carboxaldehyde [930-96-1]</p>	<p>B2862 1g 5g</p>  <p>5-Bromo-2-thiopheneboronic Acid [162607-17-2]</p>	<p>B2302 1g 5g</p>  <p>5-Bromo-2-thiophenecarboxylic Acid [7311-63-9]</p>	<p>B3386 5g 25g</p>  <p>4-Bromo-2-thiophenecarboxylic Acid [16694-18-1]</p>	<p>B3868 1g 5g</p>  <p>2-(5-Bromo-2-thienyl)-2,3-dihydro-1H-naphtho[1,8-de][1,3,2]diazaborine [927384-46-1]</p>
<p>B3812 5g</p>  <p>2-Bromo-5-hexylthiophene [211737-28-9]</p>	<p>B3738 5g 25g</p>  <p>2-Bromo-5-n-octylthiophene [172514-63-5]</p>	<p>D4146 1g 5g</p>  <p>2,4-Dibromothiophene [3140-92-9]</p>	<p>D4376 1g</p>  <p>2,5-Dibromothiophene 1,1-Dioxide [89088-95-9]</p>	<p>D3999 5g 25g</p>  <p>3,5-Dibromo-2-methylthiophene [29421-73-6]</p>
<p>T0956 25g 250g</p>  <p>2,3,5-Tribromothiophene [3141-24-0]</p>	<p>D2611 5g</p>  <p>2,5-Diiodothiophene [625-88-7]</p>	<p>T1660 5g 25g</p>  <p>3-Thiophenemethanol [71637-34-8]</p>	<p>T0864 5g 25g</p>  <p>3-Thiophenecarboxaldehyde [498-62-4]</p>	<p>T0725 25mL 100mL 500mL</p>  <p>2-Thiophenecarboxaldehyde (stabilized with HQ) [98-03-3]</p>
<p>T1084 5g 25g</p>  <p>3-Thenoic Acid [88-13-1]</p>	<p>T1475 5g 25g</p>  <p>3-Thiopheneethanol [13781-67-4]</p>	<p>T1152 5g 25g</p>  <p>Thiophene-3-acetic Acid [6964-21-2]</p>	<p>T1354 5g 25g</p>  <p>3-Thiophenemalonic Acid [21080-92-2]</p>	<p>T2388 1g 5g</p>  <p>2,3-Thiophenedicarboxaldehyde [932-41-2]</p>
<p>D2378 100mg 1g</p>  <p>3,4-Diaminothiophene Dihydrochloride [90069-81-1]</p>	<p>P1126 1g 5g</p>  <p>3-Phenylthiophene [2404-87-7]</p>	<p>D3701 1g</p>  <p>3,4-Difluoro-2,5-bis(trimethylsilyl)thiophene [347838-12-4]</p>	<p>D2074 1g 5g</p>  <p>trans-1,2-Di(2-thienyl)ethylene [13640-78-3]</p>	<p>D2077 1g 5g</p>  <p>1,4-Di(2-thienyl)-1,3-butadiene (mixture of isomers) [23354-93-0]</p>
<p>D3712 1g 5g</p>  <p>1,4-Di(2-thienyl)-1,4-butanedione [13669-05-1]</p>	<p>D0148 5g 25g</p>  <p>Dibenzothiophene [132-65-0]</p>	<p>B3525 5g 25g</p>  <p>2-Bromodibenzothiophene [22439-61-8]</p>	<p>D4274 200mg</p>  <p>3,7-Dibromodibenzo[b,d]thiophene [83834-10-0]</p>	<p>D2245 1g 5g 25g</p>  <p>2,8-Dibromodibenzothiophene [31574-87-5]</p>
<p>D4153 5g 25g</p>  <p>Dibenzothiophene 5,5-Dioxide [1016-05-3]</p>	<p>D5052 200mg 1g</p>  <p>3,7-Dibromodibenzosulfolane [83834-12-2]</p>	<p>B3464 1g 5g</p>  <p>Benzo[1,2-b:4,5-b']dithiophene [267-65-2]</p>	<p>B4158 1g 5g</p>  <p>2,6-Bis(4,4,5,5-tetramethyl-1,3,2-dioxaborolan-2-yl)-benzo[1,2-b:4,5-b']dithiophene [861398-06-3]</p>	<p>D3954 1g 5g</p>  <p>2,6-Dibromobenzo[1,2-b:4,5-b']dithiophene [909280-97-3]</p>
<p>B3793 1g 5g</p>  <p>Benzo[1,2-b:4,5-b']dithiophene-4,8-dione [32281-36-0]</p>	<p>B4396 200mg 1g</p>  <p>4,8-Bis(n-octyloxy)benzo[1,2-b:4,5-b']dithiophene [1098102-94-3]</p>	<p>D4621 200mg</p>  <p>2,6-Dibromo-4,8-bis(n-octyloxy)benzo[1,2-b:4,5-b']dithiophene [1294515-75-5]</p>	<p>B4378 200mg</p>  <p>4,8-Bis(n-octyloxy)-2,6-bis(trimethylstannyl)benzo[1,2-b:4,5-b']dithiophene [1098102-95-4]</p>	<p>B4395 200mg 1g 5g</p>  <p>4,8-Bis(2-ethylhexyloxy)benzo[1,2-b:4,5-b']dithiophene [1160823-77-7]</p>

<p>D4524 200mg</p>  <p>2,6-Dibromo-4,8-bis(2-ethylhexyloxy)benzo[1,2-b:4,5-b']dithiophene [1226782-13-3]</p>	<p>D4592 200mg</p>  <p>2,6-Dibromo-4,8-bis([2-butyl-n-octyl]oxy)benzo[1,2-b:4,5-b']dithiophene [1336893-15-2]</p>	<p>B4685 200mg</p>  <p>4,8-Bis(2-butyl-n-octyloxy)-2,6-bis(trimethylstannyl)benzo[1,2-b:4,5-b']dithiophene [1271439-08-7]</p>	<p>B4488 200mg</p>  <p>2,6-Bis(trimethylstannyl)-4,8-bis([2-n-octyldodecyl]oxy)benzo[1,2-b:4,5-b']dithiophene [1320201-22-6]</p>	<p>B4437 200mg</p>  <p>4,8-Bis(5-(2-ethylhexyl)thiophen-2-yl)-2,6-bis(trimethylstannyl)benzo[1,2-b:4,5-b']dithiophene [1352642-37-5]</p>
<p>T2729 1g 5g</p>  <p>Thieno[2,3-b]thiophene [250-84-0]</p>	<p>B3692 1g</p>  <p>2-Bromothieno[2,3-b]thiophene [25121-81-7]</p>	<p>D4032 1g</p>  <p>2,5-Dibromothieno[2,3-b]thiophene [25121-86-2]</p>	<p>D4033 1g</p>  <p>3,4-Dibromothieno[2,3-b]thiophene [53255-78-0]</p>	<p>T2513 1g 5g</p>  <p>Thieno[3,2-b]thiophene [251-41-2]</p>
<p>B3226 200mg 1g</p>  <p>3-Bromothieno[3,2-b]thiophene [25121-83-9]</p>	<p>D3678 200mg 1g 5g</p>  <p>2,5-Dibromothieno[3,2-b]thiophene [25121-87-3]</p>	<p>D3679 200mg 1g</p>  <p>3,6-Dibromothieno[3,2-b]thiophene [392662-65-6]</p>	<p>B4536 200mg 1g</p>  <p>2,5-Bis(trimethylstannyl)thieno[3,2-b]thiophene [469912-82-1]</p>	<p>T2621 1g</p>  <p>Thieno[3,2-b]thiophene-2-boronic acid [160032-40-6]</p>
<p>T2772 1g</p>  <p>2,3,5,6-Tetrabromothieno[3,2-b]thiophene [124638-53-5]</p>	<p>D4725 200mg 1g</p>  <p>2,5-Di(2-thienyl)thieno[3,2-b]thiophene [21210-90-2]</p>	<p>O0436 100mg</p>  <p>n-Octyl 4,6-Dibromo-3-fluorothieno[3,4-b]thiophene-2-carboxylate [1160823-76-6]</p>	<p>E1040 100mg</p>  <p>2-Ethylhexyl 4,6-Dibromo-3-fluorothieno[3,4-b]thiophene-2-carboxylate [1237479-38-7]</p>	<p>D3641 200mg 1g</p>  <p>Dithieno[3,2-b:2',3'-d]thiophene [3593-75-7]</p>
<p>D3823 200mg 1g</p>  <p>Dithieno[3,2-b:2',3'-d]thiophene-2-boronic acid [183960-95-4]</p>	<p>D3799 200mg 1g</p>  <p>2,6-Dibromodithieno[3,2-b:2',3'-d]thiophene [67061-69-2]</p>	<p>D3800 200mg 1g</p>  <p>3,5-Dibromodithieno[3,2-b:2',3'-d]thiophene [502764-54-7]</p>	<p>E0741 5g 25g</p>  <p>3,4-Ethylenedioxythiophene [126213-50-1]</p>	<p>D4056 1g 5g</p>  <p>2,5-Dibromo-3,4-ethylenedioxythiophene [174508-31-7]</p>
<p>E1150 1g 5g</p>  <p>3,4-Ethylenedioxythiophene-2-carboxaldehyde [204905-77-1]</p>	<p>E0743 1g 5g</p>  <p>3,4-Ethylenedioxythiophene-2,5-dicarboxylic acid [18361-03-0]</p>	<p>T2897 1g</p>  <p>3,4-Thiophenedicarboxylic anhydride [6007-85-8]</p>	<p>D4542 200mg</p>  <p>4,6-Dibromothieno[3,4-c]furan-1,3-dione [1015423-45-6]</p>	<p>O0394 200mg 1g</p>  <p>N-n-Octyl-3,4-thiophenedicarboximide [773881-43-9]</p>
<p>D4219 200mg</p>  <p>2,5-Dibromo-N-n-octyl-3,4-thiophenedicarboximide [566939-58-0]</p>	<p>D5021 200mg 1g</p>  <p>2,5-Dibromo-N-(2-ethylhexyl)-3,4-thiophenedicarboximide [1231160-83-0]</p>	<p>D4487 200mg 1g</p>  <p>4,7-Di(2-thienyl)-2,1,3-benzothiadiazole [165190-76-1]</p>	<p>B3886 100mg 1g</p>  <p>4,7-Bis(5-bromo-2-thienyl)-2,1,3-benzothiadiazole [288071-87-4]</p>	<p>B4719 200mg</p>  <p>4,7-Bis(5-bromo-4-dodecyl-2-thienyl)-2,1,3-benzothiadiazole [1179993-72-6]</p>
<p>B4542 100mg</p>  <p>4,7-Bis(5-bromo-2-thienyl)-2-n-octyl-2H-benzotriazole [1254062-41-3]</p>	<p>B4453 200mg</p>  <p>4,7-Bis(5-trimethylstannyl-2-thienyl)-2,1,3-benzothiadiazole [1025451-57-3]</p>	<p>D3969 1g 5g</p>  <p>3,6-Di(2-thienyl)-2,5-dihydropyrrolo[3,4-c]pyrrole-1,4-dione [850583-75-4]</p>	<p>O0431 200mg 1g</p>  <p>4-n-Octyl-4H-dithieno[3,2-b:2',3'-d]pyrrole [141029-75-6]</p>	<p>D4564 200mg</p>  <p>4,4-Di-n-octyldithieno[3,2-b:2',3'-d]silole [1160106-12-6]</p>





































<p>D4799 100mg</p>  <p>2,6-Dibromo-4,4-di-<i>n</i>-octyldithieno[3,2-<i>b</i>:2',3'-<i>d</i>]-silole [1160106-14-8]</p>	<p>B4102 1g 5g</p>  <p>4,4-Bis(2-ethylhexyl)-4<i>H</i>-cyclopenta[2,1-<i>b</i>:3,4-<i>b'</i>]-dithiophene [365547-20-2]</p>	<p>D4568 200mg</p>  <p>2,6-Dibromo-4,4-bis(2-ethylhexyl)-4<i>H</i>-cyclopenta[2,1-<i>b</i>:3,4-<i>b'</i>]-dithiophene [365547-21-3]</p>	<p>B3906 1g</p>  <p>3,3'-Bithiophene [3172-56-3]</p>	<p>B1276 1g 25g</p>  <p>2,2'-Bithiophene [492-97-7]</p>	
<p>O0383 1g 5g</p>  <p>5-<i>n</i>-Octyl-2,2'-bithiophene [93164-73-9]</p>	<p>D2755 1g 5g</p>  <p>5,5'-Dibromo-2,2'-bithiophene [4805-22-5]</p>	<p>D4184 200mg 1g</p>  <p>5,5'-Dibromo-4,4'-dihexyl-2,2'-bithiophene [214493-03-5]</p>	<p>D4181 1g 5g</p>  <p>4,4'-Dibromo-2,2'-bithiophene [51285-60-0]</p>	<p>D3798 1g 5g 25g</p>  <p>3,3'-Dibromo-2,2'-bithiophene [51751-44-1]</p>	
<p>D3941 1g 5g</p>  <p>3,3'-Dibromo-5,5'-bis(trimethylsilyl)-2,2'-bithiophene [207742-50-5]</p>	<p>D3928 1g 5g</p>  <p>3,4'-Dihexyl-2,2'-bithiophene [135926-93-1]</p>	<p>D4180 1g 5g</p>  <p>3,3'-Dihexyl-2,2'-bithiophene [125607-30-9]</p>	<p>D4183 1g 5g</p>  <p>5,5'-Dibromo-3,3'-dihexyl-2,2'-bithiophene [170702-05-3]</p>	<p>T2699 1g 5g 25g</p>  <p>3,3',5,5'-Tetrabromo-2,2'-bithiophene [125143-53-5]</p>	
<p>B1874 1g 5g</p>  <p>5-Bromo-2,2'-bithiophene-5'-carboxaldehyde [110046-60-1]</p>	<p>B1838 1g 5g 25g</p>  <p>2,2'-Bithiophene-5-carboxaldehyde [3779-27-9]</p>	<p>T2518 1g 5g 25g</p>  <p>5-(4,4,5,5-Tetramethyl-1,3,2-dioxaborolan-2-yl)-2,2'-bithiophene [479719-88-5]</p>	<p>F0548 1g</p>  <p>5'-Formyl-2,2'-bithiophene-5-boronic Acid [1369328-62-0]</p>	<p>B3200 1g</p>  <p>5-Bromo-5'-(4,4,5,5-tetramethyl-1,3,2-dioxaborolan-2-yl)-2,2'-bithiophene [578715-23-8]</p>	
<p>B3363 5g</p>  <p>5,5'-Bis(4,4,5,5-tetramethyl-1,3,2-dioxaborolan-2-yl)-2,2'-bithiophene [239075-02-6]</p>	<p>D4972 1g 5g</p>  <p>Dithieno-[3,2-<i>c</i>:2',3'-<i>e</i>]oxepine-4,6-dione [1043023-52-4]</p>	<p>T1196 1g 5g</p>  <p>2,2':5',2''-Terthiophene (purified by sublimation) [1081-34-1]</p>	<p>T1805 1g</p>  <p>2,2':5',2''-Terthiophene-5-carboxaldehyde [7342-41-8]</p>	<p>T2639 1g</p>  <p>2,2':5',2''-Terthiophene-5,5''-dicarboxaldehyde [13130-50-2]</p>	
<p>B2058 100mg 1g 5g</p>  <p>5''-Bromo-2,2':5',2''-terthiophene-5-carboxaldehyde [161726-69-8]</p>	<p>D4050 1g 5g</p>  <p>5,5''-Dibromo-2,2':5',2''-terthiophene [98057-08-0]</p>	<p>B4050 1g 5g</p>  <p>3'-Bromo-2,2':5',2''-terthiophene [105125-00-6]</p>	<p>Q0078 100mg</p>  <p>α-Quarterthiophene [5632-29-1]</p>	<p>Q0079 100mg 500mg</p>  <p>α-Quinquethiophene [5660-45-7]</p>	
<p>S0504 100mg 1g</p>  <p>α-Sexithiophene (purified by sublimation) [88493-55-4]</p>	<p>S0505 100mg</p>  <p>α-Septithiophene [86100-63-2]</p>	<p>O0313 100mg</p>  <p>α-Octithiophene [113728-71-5]</p>			
<h2 style="margin: 0;">Pyrroles</h2>			<p>P0574 25mL 100mL 500mL</p>  <p>Pyrrole [109-97-7]</p>	<p>P1899 1g 5g</p>  <p>Pyrrole-3-carboxylic Acid [931-03-3]</p>	<p>N0502 1g 5g</p>  <p>3-Nitropyrrole [5930-94-9]</p>



































H0809 1g 5g 3-Heptylpyrrole [878-11-5]	O0241 1g 5g 3-n-Octylpyrrole [118799-18-1]	M0414 25mL 500mL 1-Methylpyrrole [96-54-8]	E0444 1g 5g 1-Ethylpyrrole [617-92-5]	O0242 5g 1-n-Octylpyrrole [50966-65-9]
O0219 1g 5g 1-Octadecylpyrrole [89601-24-1]	C1347 5g 1-(2-Chloroethyl)pyrrole [77200-24-9]	B1622 1g 5g 1-(2-Bromoethyl)pyrrole [78358-86-8]	H0826 1g 5g 1-(2-Hydroxyethyl)pyrrole [6719-02-4]	B1623 1g 5g 1-(3-Bromopropyl)pyrrole [100779-91-7]
H0827 1g 1-(3-Hydroxypropyl)pyrrole [50966-69-3]	B2335 5g 1-Benzylpyrrole [2051-97-0]	P1122 5g 25g 1-Phenylpyrrole [635-90-5]	N0584 1g 1-(4-Nitrophenyl)pyrrole [4533-42-0]	A1022 1g 5g 1-Aminopyrrole [765-39-9]
T1677 5mL 25mL 1-(Triisopropylsilyl)pyrrole [87630-35-1]	T1386 5g 1-(p-Toluenesulfonyl)pyrrole [17639-64-4]	D3876 1g 2,5-Di(2-thienyl)-1H-pyrrole [89814-62-0]		
<h2>Carbazoles</h2>				
D4548 5g 25g 3,6-Diiodocarbazole [57103-02-3]	D2116 1g 5g 3,6-Diaminocarbazole [86-71-5]	D3952 1g 5g 3,6-Di-tert-butylcarbazole [37500-95-1]	D4275 200mg 1g 3,6-Diethynylcarbazole [909342-65-0]	D4433 1g 5g 3,6-Diphenylcarbazole [56525-79-2]
D3932 1g 5g 2,7-Dibromocarbazole [136630-39-2]	T2932 200mg 1g 1,3,6,8-Tetrabromocarbazole [55119-09-0]	D4473 200mg 1g 7H-Dibenzo[c,g]carbazole [194-59-2]	M0960 5g 25g 9-Methylcarbazole [1484-12-4]	V0021 25g 9-Vinylcarbazole [1484-13-5]
A2051 1g 5g 9-Acetyl-3,6-diiodocarbazole [606129-89-9]	E0071 25g 500g 9-Ethylcarbazole >97.0%(GC) [86-28-2]	E0247 25g 9-Ethylcarbazole >99.0%(GC) [86-28-2]	D2982 1g 5g 3,6-Dibromo-9-ethylcarbazole [33255-13-9]	D5055 1g 5g 2,7-Dibromo-9-hexylcarbazole [654676-12-7]

























<p>D4482 200mg 1g</p>  <p>3,6-Dibromo-9-<i>n</i>-octylcarbazole [79554-93-1]</p>	<p>D3935 200mg</p>  <p>2,7-Dibromo-9-<i>n</i>-octylcarbazole [726169-75-1]</p>	<p>D4830 1g 5g</p>  <p>3,6-Dibromo-9-(2-ethylhexyl)carbazole [173063-52-0]</p>	<p>D4735 200mg</p>  <p>2,7-Dibromo-9-(2-ethylhexyl)carbazole [544436-46-6]</p>	<p>D4661 200mg</p>  <p>2,7-Dibromo-9-dodecylcarbazole [544436-47-7]</p>
<p>D4490 200mg 1g</p>  <p>2,7-Dibromo-9-(9-heptadecyl)carbazole [955964-73-5]</p>	<p>O0428 200mg 1g</p>  <p>9-(9-Heptadecanyl)-2,7-bis(4,4,5,5-tetramethyl-1,3,2-dioxaborolan-2-yl)carbazole [958261-51-3]</p>	<p>B2864 1g 5g</p>  <p>9-Benzoylcarbazole [19264-68-7]</p>	<p>B3558 5g 25g</p>  <p>9-Benzylcarbazole [19402-87-0]</p>	<p>B2805 100mg</p>  <p>9-Benzylcarbazole-3,6-dicarboxaldehyde [200698-05-1]</p>
<p>B4348 1g</p>  <p>9-Benzyl-3,6-dibromocarbazole [118599-27-2]</p>	<p>P1492 1g 5g 25g</p>  <p>9-Phenylcarbazole [1150-62-5]</p>	<p>D4832 200mg 1g</p>  <p>2,7-Dibromo-9-phenylcarbazole [444796-09-2]</p>	<p>D2981 1g 5g</p>  <p>3,6-Dibromo-9-phenylcarbazole [57103-20-5]</p>	<p>D4543 1g 5g</p>  <p>3,6-Diiodo-9-phenylcarbazole [57103-21-6]</p>
<p>T2227 1g 5g</p>  <p>9-(<i>p</i>-Tolyl)carbazole [19264-73-4]</p>	<p>D5001 200mg 1g</p>  <p>3,6-Dibromo-9-(<i>p</i>-tolyl)-9<i>H</i>-carbazole [357437-74-2]</p>	<p>D5000 200mg 1g</p>  <p>3,6-Dibromo-9-(4-<i>tert</i>-butylphenyl)-9<i>H</i>-carbazole [741293-42-5]</p>	<p>D4563 200mg 1g</p>  <p>3,6-Dibromo-9-(4-bromophenyl)carbazole [73087-83-9]</p>	<p>D4978 200mg 1g</p>  <p>9-(3,5-Dibromophenyl)carbazole [750573-26-3]</p>
<p>P2350 200mg 1g</p>  <p>9-Phenyl-3,6-bis(4,4,5,5-tetramethyl-1,3,2-dioxaborolan-2-yl)carbazole [618442-57-2]</p>				
<h2>Fluorenes</h2>				
<p>D3557 5g 25g</p>  <p>2,7-Dibromo-9-fluorenone [14348-75-5]</p>	<p>D3859 5g 25g</p>  <p>2,7-Dibromo-9,9-dimethylfluorene [28320-32-3]</p>	<p>D5032 1g 5g</p>  <p>2,7-Diiodo-9,9-dimethylfluorene [144981-86-2]</p>	<p>D4569 1g 5g</p>  <p>2,7-Dibromo-9,9-dipropylfluorene [157771-56-7]</p>	<p>D3933 5g 25g</p>  <p>2,7-Dibromo-9,9-dihexylfluorene [189367-54-2]</p>
<p>B4029 1g</p>  <p>2,7-Bis(4,4,5,5-tetramethyl-1,3,2-dioxaborolan-2-yl)-9,9-dihexylfluorene [254755-24-3]</p>	<p>D3934 5g 25g</p>  <p>2,7-Dibromo-9,9-di-<i>n</i>-octylfluorene [198964-46-4]</p>	<p>B3582 1g</p>  <p>2,7-Bis(4,4,5,5-tetramethyl-1,3,2-dioxaborolan-2-yl)-9,9-di-<i>n</i>-octylfluorene [196207-58-6]</p>	<p>D4671 1g</p>  <p>9,9-Di-<i>n</i>-octylfluorene-2,7-dicarboxaldehyde [380600-91-9]</p>	<p>D4719 200mg 1g</p>  <p>2,7-Diamino-9,9-di-<i>n</i>-octylfluorene [851042-10-9]</p>

<p>D5215 5g</p>  <p>2,7-Dibromo-9,9-bis(2-ethylhexyl)fluorene (This product is only available in Japan.) [188200-93-3]</p>	<p>B4892 200mg 1g</p>  <p>9,9-Bis(2-ethylhexyl)-2,7-bis(4,4,5,5-tetramethyl-1,3,2-dioxaborolan-2-yl)-fluorene [357219-41-1]</p>	<p>B4588 1g</p>  <p>2,7-Bis(4,4,5,5-tetramethyl-1,3,2-dioxaborolan-2-yl)-9,9-didecylfluorene [711026-06-1]</p>	<p>D4052 5g 25g</p>  <p>2,7-Dibromo-9,9-didodecylfluorene (This product is only available in Japan.) [286438-45-7]</p>	<p>B4624 1g</p>  <p>2,7-Bis(4,4,5,5-tetramethyl-1,3,2-dioxaborolan-2-yl)-9,9-didodecylfluorene [749900-93-4]</p>
<p>B3552 1g 5g</p>  <p>2,7-Dibromo-9,9-bis(6-bromohexyl)fluorene [570414-33-4]</p>	<p>D4550 200mg 1g</p>  <p>2,7-Dibromo-9,9-bis[3-(dimethylamino)propyl]fluorene [673474-73-2]</p>	<p>D4547 1g 5g</p>  <p>2,7-Dibromo-9,9-diphenylfluorene [186259-63-2]</p>	<p>B3645 1g 5g</p>  <p>9,9-Bis(4-bromophenyl)fluorene [128406-10-0]</p>	<p>B1549 25g 100g</p>  <p>9,9-Bis(4-aminophenyl)fluorene [15499-84-0]</p>
<p>B2654 1g</p>  <p>9,9-Bis(4-aminophenyl)fluorene (purified by sublimation) [15499-84-0]</p>	<p>D3781 1g</p>  <p>2,7-Dibromo-9,9'-spirobi[9H-fluorene] [171408-84-7]</p>	<p>D3872 1g</p>  <p>2,2'-Dibromo-9,9'-spirobi[9H-fluorene] [67665-47-8]</p>	<p>T2597 1g 5g</p>  <p>2,2',7,7'-Tetrabromo-9,9'-spirobi[9H-fluorene] [128055-74-3]</p>	<p>T2759 1g</p>  <p>2,2',7,7'-Tetrabromo-9,9'-bifluorenylidene [27192-91-2]</p>
<p>D4573 1g</p>  <p>2,8-Dibromoindeno[1,2-b]fluorene-6,12-dione [853234-57-8]</p>	<p>Aryls</p>			
		<p>B3863 200mg 1g</p>  <p>Benzocyclobutene [694-87-1]</p>	<p>D0170 25g 100g 500g</p>  <p>1,4-Dibromobenzene [106-37-6]</p>	<p>B4633 5g 25g</p>  <p>1,4-Benzenediboronic Acid Bis(pinacol) Ester [99770-93-1]</p>
<p>B5070 1g 5g</p>  <p>1,3-Benzenediboronic Acid Bis(pinacol) Ester [196212-27-8]</p>	<p>D0608 5g 25g</p>  <p>1,4-Diiodobenzene [624-38-4]</p>	<p>D1864 10g 25g</p>  <p>1,4-Dibromo-2,5-difluorobenzene [327-51-5]</p>	<p>D2359 5g 25g</p>  <p>1,4-Dibromonaphthalene [83-53-4]</p>	<p>D4154 1g</p>  <p>2,6-Dibromonaphthalene [13720-06-4]</p>
<p>D3624 5g 25g</p>  <p>2,7-Dibromonaphthalene [58556-75-5]</p>	<p>D4660 1g 5g</p>  <p>1,5-Dibromonaphthalene [7351-74-8]</p>	<p>D1695 25g 250g</p>  <p>4,4'-Dibromobiphenyl [92-86-4]</p>	<p>B2490 1g 5g</p>  <p>4,4'-Biphenyldiboronic Acid [4151-80-8]</p>	<p>B4166 1g 5g</p>  <p>4,4'-Bis(4,4,5,5-tetramethyl-1,3,2-dioxaborolan-2-yl)-biphenyl [207611-87-8]</p>
<p>B3151 1g</p>  <p>4,4'-Bis(5,5-dimethyl-1,3,2-dioxaborolan-2-yl)-biphenyl [5487-93-4]</p>	<p>D2315 25g 250g</p>  <p>4,4'-Diiodobiphenyl [3001-15-8]</p>	<p>B2860 1g</p>  <p>4'-Bromo-4-biphenylboronic Acid [480996-05-2]</p>	<p>D4560 200mg</p>  <p>4,4'-Dibromo-2,2'-diiodobiphenyl [852138-89-7]</p>	<p>D4608 1g 5g</p>  <p>4,4'-Diiodo-2,2'-dimethylbiphenyl [69571-02-4]</p>

Conducting Polymer Research Reagents

<p>D2759 5g 25g</p>  <p>4,4'-Diiodo-3,3'-dimethylbiphenyl [7583-27-9]</p>	<p>D1696 1g 5g</p>  <p>4,4'-Dibromooctafluorobiphenyl [10386-84-2]</p>	<p>D4363 1g 5g</p>  <p>2,7-Dibromophenanthrene-9,10-dione [84405-44-7]</p>	<p>D3171 1g 5g</p>  <p>2,6-Dibromoanthracene [186517-01-1]</p>	<p>D3182 1g 5g</p>  <p>2,6-Dibromoanthraquinone [633-70-5]</p>
<p>D0166 25g 100g 500g</p>  <p>9,10-Dibromoanthracene [523-27-3]</p>	<p>D5050 200mg 1g</p>  <p>2,7-Dibromophenanthrene [62325-30-8]</p>	<p>D3372 1g 5g</p>  <p>4,4''-Dibromo-p-terphenyl [17788-94-2]</p>	<p>D3534 1g 5g</p>  <p>4,4''-Diiodo-p-terphenyl [19053-14-6]</p>	<p>D4801 200mg 1g</p>  <p>2,7-Dibromotriphenylene [888041-37-0]</p>
<p>D4236 200mg 1g</p>  <p>6,12-Dibromochrysene [131222-99-6]</p>	<p>D3169 1g</p>  <p>2,7-Dibromopyrene [102587-98-4]</p>	<p>B3956 1g 5g</p>  <p>2,7-Bis(4,4,5,5-tetramethyl-1,3,2-dioxaborolan-2-yl)pyrene [688756-58-3]</p>	<p>D4187 1g 5g</p>  <p>1,3-Dibromo-7-tert-butylpyrene [1005771-04-9]</p>	<p>D3167 1g 5g</p>  <p>1,6-Dibromopyrene [27973-29-1]</p>
<p>T2716 5g</p>  <p>1,3,6,8-Tetrabromopyrene [128-63-2]</p>	<p>D4254 1g 5g</p>  <p>10,10'-Dibromo-9,9'-bianthracene [121848-75-7]</p>	<p>N0369 25g 250g</p>  <p>Naphthalene-1,4,5,8-tetracarboxylic Dianhydride [81-30-1]</p>	<p>N0755 1g 5g</p>  <p>Naphthalene-1,4,5,8-tetracarboxylic Dianhydride (purified by sublimation) [81-30-1]</p>	<p>D4339 1g 5g</p>  <p>2,6-Dibromonaphthalene-1,4,5,8-tetracarboxylic Dianhydride [83204-68-6]</p>
<p>N0536 1g 5g</p>  <p>1,4,5,8-Naphthalenetetracarboxydiimide [5690-24-4]</p>	<p>D4946 1g 5g</p>  <p>2,6-Dibromo-N,N'-di-n-octyl-1,8,4,5-naphthalenetetracarboxydiimide [926643-78-9]</p>	<p>D4945 1g 5g</p>  <p>2,6-Dibromo-N,N'-bis(2-ethylhexyl)-1,8,4,5-naphthalenetetracarboxydiimide [1088205-02-0]</p>	<p>P0972 25g 100g 500g</p>  <p>3,4,9,10-Perylenetetracarboxylic Dianhydride [128-69-8]</p>	<p>P2102 1g</p>  <p>3,4,9,10-Perylenetetracarboxylic Dianhydride (purified by sublimation) [128-69-8]</p>
<p>P0984 25g</p>  <p>3,4,9,10-Perylenetetracarboxylic Diimide [81-33-4]</p>	<p>D3947 200mg 1g</p>  <p>N,N'-Ditridecyl-3,4,9,10-perylenetetracarboxylic Diimide [95689-92-2]</p>	<p>D5211 1g 5g</p>  <p>4,4'-Dibromobenzophenone [3988-03-2]</p>	<p>B3996 1g 5g</p>  <p>Bis(4-bromophenyl)acetylene [2789-89-1]</p>	<p>B4989 1g 5g</p>  <p>Bis(2-bromophenyl)acetylene [38399-13-2]</p>
<p>D3094 1g 5g</p>  <p>4,4'-Diiodo-trans-stilbene [201861-91-8]</p>	<p>D4724 200mg 1g</p>  <p>3,6-Dibromophthalic Anhydride [25834-16-6]</p>	<p>D4591 1g</p>  <p>1,3-Diiodoazulene [36044-42-5]</p>		
<h2>Dichlorosilanes</h2>				
<p>D0358 25mL 100mL 500mL</p>  <p>Dichlorodimethylsilane [75-78-5]</p>				
<p>D1976 5g 25g</p>  <p>Dichlorodiethylsilane [1719-53-5]</p>				
<p>D2603 5g 25g</p>  <p>Dichlorodiisopropylsilane [7751-38-4]</p>				

D2704 5g 25g  Dibutylchlorosilane [3449-28-3]	D2469 5g  Di-tert-butylchlorosilane [18395-90-9]	D2710 5g  Dichlorodipentylsilane [18037-39-3]	D2712 5g 25g  Dichlorodihexylsilane [18204-93-8]	D4336 1g  Dichloro-di- <i>n</i> -octylsilane [18416-07-4]
H0855 25mL  Dichlorohexylmethylsilane [14799-94-1]	O0265 25mL  Dichloro(methyl)- <i>n</i> -octylsilane [14799-93-0]	D2348 25mL  Dichlorodecylmethylsilane [18051-88-2]	D2349 25mL  Dichlorododecylmethylsilane [18407-07-3]	O0264 25mL  Dichloro(methyl)octadecylsilane [5157-75-5]
C1399 5g 25g  Dichlorocyclohexylmethylsilane [5578-42-7]	P0790 25g 250g  Dichloro(methyl)phenylsilane [149-74-6]	D0362 25g 100g 500g  Dichlorodiphenylsilane [80-10-4]	D2144 5g 25g  1,2-Dichlorotetramethyl-disilane [4342-61-4]	T1590 5g 25g  1,1,2,2-Tetrachloro-1,2-dimethyldisilane [4518-98-3]
D4729 5g 25g  Dichloro(methyl)-(3,3,3-trifluoropropyl)silane [675-62-7]				
<h2 style="background-color: #e91e63; color: white; padding: 10px;">Others</h2>				
S0829 1g 5g  Selenophene [288-05-1]	F0782 1g 5g  2-Formylselenophene (stabilized with HQ) [25109-26-6]	D4028 1g 5g  2,5-Dibromoselenophene [1755-36-8]		
D3960 5g  2,5-Dibromothiazole [4175-78-4]	D4559 200mg  4,4'-Dinonyl-2,2'-bithiazole [180729-91-3]	D4558 200mg  5,5'-Dibromo-4,4'-dinonyl-2,2'-bithiazole [172100-44-6]	D4605 200mg  4,4'-Dibromo-2,2'-bis(triisopropylsilyl)-5,5'-bithiazole [1223559-98-5]	B0921 10g 25g  2,1,3-Benzothiadiazole [273-13-2]
D3842 1g 5g 25g  4,7-Dibromo-2,1,3-benzothiadiazole [15155-41-6]	B3573 1g 5g  4,7-Bis(4,4,5,5-tetramethyl-1,3,2-dioxaborolan-2-yl)-2,1,3-benzothiadiazole [934365-16-9]	D4750 200mg 1g  4,7-Dibromo-5-fluoro-2,1,3-benzothiadiazole [1347736-74-6]	D4734 200mg  4,7-Dibromo-5,6-difluoro-2,1,3-benzothiadiazole [1295502-53-2]	B4440 200mg  5,6-Bis(<i>n</i> -octyloxy)-2,1,3-benzothiadiazole [1254353-37-1]
D4529 200mg  4,7-Dibromo-5,6-di- <i>n</i> -octyloxy-2,1,3-benzothiadiazole [1192352-08-1]	D4622 200mg  4,7-Dibromo-5,6-bis(dodecyloxy)-2,1,3-benzothiadiazole [1313876-00-4]	D4442 200mg 1g  4,7-Dibromo[1,2,5]-thiadiazolo[3,4- <i>c</i>]pyridine [333432-27-2]	B4473 1g 5g  2,1,3-Benzoxadiazole [273-09-6]	B4313 200mg 1g  2,1,3-Benzoselenadiazole [273-15-4]

<p>D3944 1g</p>  <p>4,7-Dibromo-2,1,3-benzoselenadiazole [63224-42-0]</p>	<p>O0430 200mg 1g</p>  <p>2-<i>n</i>-Octylbenzotriazole [112642-69-0]</p>	<p>D4514 200mg</p>  <p>4,7-Dibromo-2-(6-bromohexyl)benzotriazole [890704-02-6]</p>	<p>D4461 200mg</p>  <p>5,8-Dibromo-2,3-diphenylquinoxaline [94544-77-1]</p>	<p>D4653 200mg</p>  <p>5,8-Dibromo-2,3-bis(4-methoxyphenyl)quinoxaline [162967-90-0]</p>
<p>D4821 1g 5g</p>  <p>2,8-Dibromodibenzofuran [10016-52-1]</p>	<p>D4835 200mg 1g</p>  <p>2,8-Diiododibenzofuran [5943-11-3]</p>	<p>D4308 1g 5g</p>  <p>9,9-Dimethyl-9H-9-silafluorene [13688-68-1]</p>	<p>D4631 1g 5g</p>  <p>2,5-Dibromo-1,1-dimethyl-3,4-diphenylsilole [686290-22-2]</p>	<p>D3209 200mg 1g</p>  <p>3,8-Dibromo-1,10-phenanthroline [100125-12-0]</p>
<p>B5087 200mg 1g</p>  <p>Bis(4-bromophenyl)diphenylsilane [18733-91-0]</p>	<p>B2685 1g 5g</p>  <p>Bis(4-iodophenyl)amine [20255-70-3]</p>	<p>D1910 1g 5g</p>  <p>4,4'-Dibromotriphenylamine [81090-53-1]</p>	<p>B2843 1g</p>  <p>Bis(4-formylphenyl)phenylamine (purified by sublimation) [53566-95-3]</p>	<p>D3303 1g</p>  <p>4,4'-Dibromo-4''-phenyltriphenylamine [884530-69-2]</p>
<p>D3302 1g</p>  <p>4,4'-Dibromo-4''-cyclohexyltriphenylamine [1397196-59-6]</p>	<p>B1336 1g</p>  <p><i>N,N'</i>-Bis(4-chlorophenyl)-<i>N,N'</i>-diphenyl-1,4-phenylenediamine [113703-66-5]</p>	<p>B3991 200mg</p>  <p>4-[<i>N,N'</i>-Bis(4-bromophenyl)amino]benzaldehyde [25069-38-9]</p>	<p>B4793 200mg 1g</p>  <p><i>N,N'</i>-Bis(4-formylphenyl)-<i>N,N'</i>-diphenylbenzidine [137911-28-5]</p>	<p>T1361 5g 25g</p>  <p>Tris(4-bromophenyl)amine [4316-58-9]</p>
<p>T1724 5g 25g</p>  <p>Tris(4-iodophenyl)amine [4181-20-8]</p>	<p>T2310 500mg</p>  <p>Tris(4-formylphenyl)amine [119001-43-3]</p>	<p>T2332 1g 5g</p>  <p>Tris(4-aminophenyl)amine [5981-09-9]</p>	<p>T3315 200mg</p>  <p>Tris[4-(5-bromothiophen-2-yl)phenyl]amine [339985-36-3]</p>	

Ordering and Customer Service

TCI AMERICA

Tel : 800-423-8616 / 503-283-1681
 Fax : 888-520-1075 / 503-283-1987
 E-mail : Sales-US@TCIchemicals.com

TCI EUROPE N.V.

Tel : +32 (0)3 735 07 00
 Fax : +32 (0)3 735 07 01
 E-mail : Sales-EU@TCIchemicals.com

TCI Deutschland GmbH

Tel : +49 (0)6196 64053-00
 Fax : +49 (0)6196 64053-01
 E-mail : Sales-DE@TCIchemicals.com

Tokyo Chemical Industry UK Ltd.

Tel : +44 (0)1865 784560
 Fax : +44 (0)1865 784561
 E-mail : Sales-UK@TCIchemicals.com

TCI Chemicals (India) Pvt. Ltd.

Tel : 1800 425 7889 / 044-2262 0909
 Fax : 044-2262 8902
 E-mail : Sales-IN@TCIchemicals.com

梯希爱(上海)化成工业发展有限公司

Tel : 800-988-0390 / 021-67121386
 Fax : 021-6712-1385
 E-mail : Sales-CN@TCIchemicals.com

TOKYO CHEMICAL INDUSTRY CO., LTD.

Tel : +81 (0)3-5640-8878
 Fax : +81 (0)3-5640-8902
 E-mail : globalbusiness@TCIchemicals.com

Availability, price or specification of the listed products are subject to change without prior notice. Reproduction forbidden without the prior written consent of Tokyo Chemical Industry Co., Ltd.