



磁性金属配合物


β-二酮配合物




Bis(2,4-pentanedionato)-cobalt(II)
25g [B2681]




Bis(2,4-pentanedionato)-cobalt(II) Dihydrate
25g / 500g [C0373]




Bis(trifluoro-2,4-pentanedionato)-cobalt(II) Hydrate
5g [T0746]




Bis(hexafluoroacetylacetonato)-cobalt(II) Hydrate
5g [H0553]




Tris(2,4-pentanedionato)iron(III)
25g / 100g / 500g [I0079]




Tris(trifluoro-2,4-pentanedionato)iron(III)
5g [T0750]




Tris(hexafluoroacetylacetonato)iron(III)
1g [H0555]




Tris(dibenzoylmethanato) Iron
5g / 25g [T1686]




Bis(2,4-pentanedionato)copper(II)
25g / 250g [C0384]




Bis(trifluoro-2,4-pentanedionato)copper(II)
5g [T0752]




Bis(hexafluoroacetylacetonato)-copper(II) Hydrate
1g / 5g [H0554]




Bis(2,4-pentanedionato)-manganese(II) Dihydrate
25g [M0042]




Bis(hexafluoroacetylacetonato)-manganese(II) Hydrate
1g / 5g [H0557]




Bis(trifluoro-2,4-pentanedionato)manganese(II)
5g [T0685]




Tris(2,4-pentanedionato)-manganese(III)
25g / 100g / 500g [M0043]



Bis(2,4-pentanedionato)-vanadium(IV) Oxide
25g [V0016]




Tris(trifluoro-2,4-pentanedionato)chromium(III)
1g [T0932]




Tris(2,4-pentanedionato)-chromium(III)
25g / 100g [C0336]


卟啉·酞菁配合物




Cobalt(II) Phthalocyanine
5g / 25g [P0887]




Cobalt(II) Tetraphenylporphyrin
1g / 5g [C2986]




[5,10,15,20-Tetrakis(4-methoxyphenyl)-porphyrinato]cobalt(II)
1g / 5g [T1861]




Copper(II) Phthalocyanine (purified by sublimation)
1g [P1628]



Copper(II) Phthalocyanine (α-form)
25g / 250g [P1005]




Copper(II) Phthalocyanine (β-form)
25g / 100g / 500g [P1006]




Copper(II) 2,9,16,23-Tetra-tert-butylphthalocyanine
100mg [T1256]


二硫烯配合物




Bis(tetrabutylammonium) Bis(maleonitriledithiolato)nickel(II)
1g [B1371]




Tetrabutylammonium Bis(maleonitriledithiolato)nickel(III)
1g [T1415]



Bis(dithiobenzil)nickel(II)
1g / 25g [B1350]




Bis[4,4'-dimethoxy(dithiobenzil)]nickel(II)
1g / 5g [B4360]




Bis(4-dimethyl-amino-dithiobenzil)nickel(II)
1g [B4361]




Tetrabutylammonium Bis(4-methyl-1,2-benzenedithiolato)nickel(III)
200mg / 1g [T3220]




Tetrabutylammonium Bis(3,6-dichloro-1,2-benzenedithiolato)nickel(III)
200mg [T3204]




Tetrabutylammonium Bis(3,4,6-trichloro-1,2-benzenedithiolato)nickel(III)
200mg [T3205]



Tetrabutylammonium Bis(1,3-dithiole-2-thione-4,5-dithiolato)nickel(III)
100mg [T1272]



Tetrabutylphosphonium Bis(1,3-dithiole-2-thione-4,5-dithiolato)nickel(III)
1g [T1416]




Bis(tetrabutylammonium) Bis(1,3-dithiole-2-thione-4,5-dithiolato)palladium(II)
100mg / 1g [B1437]




Bis(tetrabutylammonium) Bis(1,3-dithiole-2-thione-4,5-dithiolato)platinum(II)
100mg [B1438]


二硫代氨基甲酸配合物




Iron(III) Dimethyl-dithiocarbamate
25g / 500g [D1267]




Iron(III) Diethyl-dithiocarbamate
25g [D0489]



Copper(II) Dimethyl-dithiocarbamate
25g [D1022]



Copper(II) Diethyl-dithiocarbamate
5g / 25g [D0487]




Copper(II) Bis(2-hydroxyethyl)dithiocarbamate
5g / 25g [B0491]


其它配合物




Iron(III) Hexacyanoferrate(II)
25g [I0849]




Copper(II) Acetate Monohydrate
25g / 500g [C2346]




Chromium(III) Pyridine-2-carboxylate
5g / 25g [P1347]



Bis(8-quinolinolato)copper(II)
25g [B1677]



Bis(8-quinolinolato)copper(II) (purified by sublimation)
1g [B3834]



N,N'-Ethylenebis(salicylidene-iminato)cobalt(II)
25g / 100g / 500g [S0318]

更多信息, 请查看我们的主页: www.TCIchemicals.com

磁性金属

