

多糖類

Alginate Acid	25g / 500g [A0733]
Alginate Acid Sodium Salt 100-200	25g / 500g [A0205]
Alginate Acid Sodium Salt 300-400	25g / 500g [O0583]
Alginate Acid Sodium Salt 500-600	25g / 500g [O0584]
Alginate Acid Sodium Salt 800-900	25g / 500g [O0585]
Alginate Acid Calcium Salt	25g / 500g [A0738]
Amylopectin Hydrate (Amylose free), from Waxy Corn	25g / 500g [A0456]
Arabic Gum	100g [A3553]
(+)-Arabinogalactan from Larch Wood	25g / 100g [A1328]
Carboxymethyl Cellulose Sodium (n=approx. 500)	25g / 500g [C0045]
Carboxymethyl Cellulose Sodium (n=approx. 1,050)	25g / 500g [C0603]
Carboxymethyl dextran Sodium Salt (Mw.=ca. 10000)	1g / 5g [C3250]
Carboxymethyl dextran Sodium Salt (Mw.=ca. 40000)	1g / 5g [C3251]
κ-Carrageenan	25g / 500g [C1804]
ι-Carrageenan	25g / 500g [C1805]
λ-Carrageenan (Low-viscosity)	1g / 5g [C2871]
λ-Carrageenan (High-viscosity)	25g / 500g [C3313]
Chitin	25g / 250g [C0072]
Chitosan (5-20mPa·s, 0.5% in 0.5% Acetic Acid at 20°C)	25g / 100g / 500g [C2395]
Chitosan (20-100mPa·s, 0.5% in 0.5% Acetic Acid at 20°C)	25g / 100g [C2396]
Chitosan (200-600mPa·s, 0.5% in 0.5% Acetic Acid at 20°C)	25g / 500g [C0831]
Chondroitin Sulfate Sodium Salt	25g / 100g [C0335]
Curdlan	5g / 25g [C3859]
Dermatan Sulfate Sodium Salt	20mg / 100mg [D3672]
Dextran 40 (Mw.=ca. 40000)	25g / 100g / 500g [D1448]
Dextran 70 (Mw.=ca. 70000)	25g / 100g / 500g [D1449]
Dextran Sulfate Sodium from Dextran of Mw. Approx. 8000	25g [D5143]
Dextran Sulfate Sodium from Dextran of Mw. Approx. 40000	5g / 25g [D5144]
Dextrin	100g / 500g [D4657]
Dextrin (Soluble fiber)	100g / 500g [D5658]
Ethyl Cellulose [9-11mPa·s, 5% in Toluene + Ethanol (80:20) at 25°C]	25g / 500g [E0265]
Ethyl Cellulose [18-22mPa·s, 5% in Toluene + Ethanol (80:20) at 25°C]	25g / 500g [E0072]
Ethyl Cellulose [45-55mPa·s, 5% in Toluene + Ethanol (80:20) at 25°C]	25g / 500g [E0266]
Ethyl Cellulose [90-110mPa·s, 5% in Toluene + Ethanol (80:20) at 25°C]	25g / 500g [E0290]
Glucan from Black Yeast	1g / 5g [G0331]
Guar Gum	25g / 500g [G0478]

Heparin Sodium Salt from Hog intestine	100mg / 1g [H0393]
Hyaluronic Acid from Cockscomb	1g [H0595]
Hyaluronic Acid Sodium Salt from Cockscomb	100mg / 1g [H0603]
Hyaluronic Acid Potassium Salt from Cockscomb	1g [H0652]
Hyaluronic Acid from Bacteria	1g [H1807]
Hyaluronic Acid Sodium Salt from Bacteria	100mg / 1g [H1791]
Hyaluronic Acid Potassium Salt from Bacteria	1g [H1808]
Hyaluronic Acid, High-Sulfated	10mg [H1739]
Hyaluronic Acid, Low-Sulfated	10mg [H1740]
Hydroxyethyl Cellulose (200-300mPa-s, 2% in Water at 20°C)	25g / 500g [H0242]
Hydroxyethyl Cellulose (800-1500mPa-s, 2% in Water at 20°C)	25g / 500g [H0418]
Hydroxyethyl Cellulose (4500-6500mPa-s, 2% in Water at 25°C)	25g / 500g [H0392]
Hydroxypropyl Cellulose (3-6mPa-s, 2% in Water at 20°C)	25g / 500g [H0473]
Hydroxypropyl Cellulose (6-10mPa-s, 2% in Water at 20°C)	25g / 500g [H0474]
Hydroxypropyl Cellulose (150-400mPa-s, 2% in Water at 20°C)	25g / 500g [H0386]
Hydroxypropyl Cellulose (1000-4000mPa-s, 2% in Water at 20°C)	25g / 500g [H0475]
Inulin (by Enzymatic synthesis)	25g / 500g [I1067]
Konjac Glucomannan	25g / 100g [K0075]
Laminaran from <i>Eisenia Bicyclis</i>	1g / 25g [L0088]
Locust Bean Gum	25g / 100g [L0459]
Methyl Cellulose (13-18mPa-s, 2% in Water at 20°C)	25g / 500g [M0290]
Methyl Cellulose (20-30mPa-s, 2% in Water at 20°C)	25g / 500g [M0291]
Methyl Cellulose (80-120mPa-s, 2% in Water at 20°C)	25g / 500g [M0292]
Methyl Cellulose (350-550mPa-s, 2% in Water at 20°C)	25g / 500g [M0293]
Methyl Cellulose (1000-1800mPa-s, 2% in Water at 20°C)	25g / 500g [M0294]
Methyl Cellulose (3500-5600mPa-s, 2% in Water at 20°C)	25g / 500g [M0185]
Methyl Cellulose (7000-10000mPa-s, 2% in Water at 20°C)	25g / 500g [M0295]
Pectin from Citrus	25g / 500g [P0024]
Pullulan	25g / 100g / 500g [P0978]
Tamarind Gum from Tamarind seed, Polysaccharide	25g / 500g [T0909]
Xanthan Gum	25g / 100g / 500g [X0048]
Xylan from Corn Core	25g / 100g [X0078]

これらの製品はすべて“試薬”です。試験・研究用にご使用ください。

上記以外の化合物についても取り揃えています。各製品の詳細はTCIのウェブサイトでも ▶▶▶ TCI ポリサッカライド

東京化成工業株式会社

試薬製品について

■本社営業部 〒103-0001 東京都中央区日本橋小伝馬町 16-12 T-PLUS 日本橋小伝馬町8階
Tel: 03-3668-0489 Fax: 03-3668-0520 E-mail: Sales-JP@TCIchemicals.com

■大阪営業部 〒541-0041 大阪府大阪市中央区北浜 1-1-21 第2中井ビル1階
Tel: 06-6228-1155 Fax: 06-6228-1158 E-mail: osaka-s@TCIchemicals.com

スケールアップ、受託サービス(合成・開発・製造)について

□化成品営業部 〒103-0001 東京都中央区日本橋小伝馬町 16-12 T-PLUS 日本橋小伝馬町8階
Tel: 03-5651-5171 Fax: 03-5640-8021 E-mail: finechemicals@TCIchemicals.com

弊社製品取扱店

本誌掲載の化学品は試験・研究用にのみ使用するものです。化学知識のある専門家以外の方のご使用はお避けください。品目や製品情報等、掲載内容の変更を予告なく行う場合があります。内容の一部または全部の無断転載・複製はご遠慮ください。